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Installation Qualification of Sterilization

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According to ISO11135-1: 2007 and the standard of GB18279-2000 (*Medical devices –Validation and routine control of ethylene oxide sterilization*) requirements, did the effective IQ of sterilization chamber, to demonstrate that the sterilization chamber had been finished the effective installation process, and can run normally, meet the requirements of ISO11135-1: 2007.

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There are three qualification items of IQ of sterilization chamber:

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System installation check: check the correctness of the equipment and the design of equipment is right or wrong and can run or not.

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Information collect: make sure we have all documents and calibration certificates, collect and appropriate storage before doing the installation test.

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Installation test: test the equipment detail function can use normally or not, and whether fits the regulate requirement.

3.1 系统安装的检查 System installation check:

3.1.1 安装的正确性 The correctness of installation

序号 Item	鉴定项目 Qualification Item	检查记录 Check Record	是否符合要求 Accepted or Rejected
1	箱体安装位置 Chamber installation position		☐ Acc ☐ Rej
2	柜体尺寸 Chamber size		☐ Acc ☐ Rej
3	箱体安装水平 Level of chamber installation		☐ Acc ☐ Rej
4	箱体倾斜度 Gradient of chamber		☐ Acc ☐ Rej
5	箱体周围无障碍 Fraise of chamber round		☐ Acc ☐ Rej
6	门与箱体装配 Door and chamber installation		☐ Acc ☐ Rej
7	门的活动无障碍 Door can moving and no fraise		☐ Acc ☐ Rej
8	管道安装平直 Straight of pipeline installation		☐ Acc ☐ Rej

Performed by:

Verified by:

1.2 安装的完整、准确性 The integrity and veracity of installation

序号 Item	鉴定项目 Qualification Item	检查记录 Check Record	是否符合要求 Accepted or Rejected
1	控制系统完整性、准确性 Integrity and veracity control system		☐ Acc ☐ Rej
2	供电系统完整性、准确性 Integrity and veracity electric system		☐ Acc ☐ Rej
3	供水系统完整性、准确性 Integrity and veracity water system		☐ Acc ☐ Rej
4	供气系统完整性、准确性 Integrity and veracity t gas system		☐ Acc ☐ Rej
6	加湿系统完整性、准确性 Integrate and correct control system		☐ Acc ☐ Rej
8	管道标识符合性、准确性 Acceptance and veracity pipeline sign		☐ Acc ☐ Rej
9	仪表及传感器位置和数量 The placement and number of meters and sensors		☐ Acc ☐ Rej

Performed by:

Verified by:

1.3 电器控制系统安装 Electric controlling system installation

项目名称 Item name	规格\型号\技术指标的符合性 The acceptance of spec, model, technique requirement	安装准确性 Installation veracity	标识准确性 Sign veracity	检查记录 Check record
开关 switch	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
按键 pushbutton	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
链接 interlinkage	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
传感器 sensor	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	

仪表 meter	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
动力器具 motivity apparatus	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
加热器件 Heating device	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
加湿器件 Humidity device	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	
安全接地 safe earth	☐ Yes ☐ No	☐ Yes ☐ No	☐ Yes ☐ No	

Performed by:

Verified by:

1.4 计算机系安装 Computer system installation

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Verified by:

3. 安装测试 Installation test:

3.1 辅助设备的运行测试 Assistant equipment run test

辅助设备 Assistant equipments	运行时间 Run time	噪 音 Noise	旋转方向 Circumgyrate direction	异常记录 Unconventional record
真空泵 Vacuum pump	10min			
气泵 Pump	10min			
循环泵 Circulate pump	20min			
气化泵 Gasification pump	20min			
加热系统 (蒸气、水箱) Heating system (steam and water tank)	30min		N/A	
蒸汽发生器 Steam produce equipment	30min	N/A	N/A	

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Conclusion: After finishing the installation qualification, the sterilization chamber and each assistant part's (vacuum system, cycle system, heating system, compress air system, electric controlling system, steam system) documents are complete and install correctly, system position nicety, the blueprint of equipment, technique documents and interrelated documents also are collected, and it meets the requirements of 5.2 of ISO11135-1994, the calibration certificates as follows attachment 1.

Prepared By: _____ Date: _____
ETO Sterilization Process Engineer

Accepted By: _____ Date: _____
ETO Sterilization Assurance Engineer

Accepted By: _____ Date: _____
Sterilization Manager3

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Operational Qualification of Sterilization

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According to ISO11135-1:2007 and the standard of GB18279-2000 (Medical devices –Validation and routine control of ethylene oxide sterilization) requirements, did the effective OQ of sterilization chamber, to demonstrate that installation process of the sterilization chamber is effective, the equipment can run normally, and meet the requirements of ISO11135-1:2007.

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There are five qualification items of OQ of sterilization chamber:

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Vacuum speed test: to calculate the vacuum speed by setting vacuum parameter during the vacuum process.

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Pressure leak test: to observe and calculate the leak speed of chamber by keeping pressure during the subatmospheric and superatmospheric process.

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Temperature uniformity test: to determine the temperature in a number of locations distributed throughout the sterilization area to test the temperature uniformity of sterilization chamber.

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Humidity test: to determine the validity of humidity injection system and calculate the variety relation between humidity and pressure by injecting humidity into the sterilization chamber.

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Set the vacuum parameter -50KP, record the time that reach to the setting parameters, and calculate the vacuum speed.

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The data record as follows attachment 1.

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Accords with the requirements of ISO11135-1:2007, do the leak test of the chamber during the subatmospheric and superatmospheric process; subatmospheric: -50KP, keep pressure 60mins; superatmospheric: 50KP, keep pressure 60mins, observe and calculate the leak speed of the sterilization chamber of the keeping pressure process.

2. 2 物理参数运行结果 The run results of physical parameters

真空度 Vacuum	开始压力 Starting pressure (kpa)	结束压力 End pressure (kpa)	压力变化 Pressure changes (kpa)	所用时间 Using time (min)	泄漏速率 Leak speed (kpa/min)	泄漏速率要求 Leak speed Requirements (kpa/min)
50Kpa						≤0.1000
-50Kpa						≤0.1000

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The data record as follows attachment 7.

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Accords with the requirements of ISO11135-1: 2007, the volume of sterilization chamber is no more than 1 m³, should use 3ea temperature sensors to monitor the sterilization chamber temperature, add a minimum of one temperature sensor for adding 3m³ of sterilization chamber if the chamber volume is more than 1 m³ during the empty loading condition; and it need 3ea temperature sensors at least to

monitor the sterilization chamber temperature of 1 m³ sterilization chamber (in fact, we used 20ea temperature sensors), the temperature sensor placement diagram, please see the attachment 4.

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Data and temperature trend, sterilization data and trend as follows attachment 2

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The inner wall temperature uniformity of sterilization chamber is good by the operational qualification.

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
Accords with the requirements of ISO11135-1:2007, the volume of sterilization chamber is no more than 3 m^3 , should use 3 ea temperature sensors to monitor the sterilization chamber temperature, add a minimum of one temperature sensor for adding 20 m^3 of sterilization chamber if the chamber volume is more than 5 m^3 during the empty loading condition; and it need 20ea temperature sensors at least to monitor the sterilization chamber temperature of 1 m^3 sterilization chamber (in fact, we used 20ea temperature sensors), the temperature sensor placement diagram, please see the attachment 5.

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Data and temperature trend, sterilization data and trend as follows attachment 4.

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The space temperature uniformity of sterilization chamber is good by the operational qualification.

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
Inject humidity into the sterilization chamber when the temperature of chamber is ____℃ and pressure is ____KP, record the pressure variety when the humidity rises from 30%RH to 80%RH.. record and analyze data, the Humidity sensor placement diagram, please see the attachment 6.

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
PQ of Sterilization Chamber Of ETO Department

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According to ISO11135-1:2007 and the standard of GB18279-2000 (Medical devices –Validation and routine control of ethylene oxide sterilization) requirements, did the effective PQ of sterilization chamber, the equipment can run normally, and meet the requirements of ISO11135-1:2007.

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Project in allusion to new Sterilization as follows

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Temperature/Humidity uniformity test

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This method is that keeping all parameters in same condition except time; ensure that after the shortest EO sterilization time, no bacteria survive, and repeat once same test to certify this shortest sterilization time. 2 times test should show that there are no bacteria in BI. According to regulation the normal sterilization time should be at least 2 times of the shortest sterilization

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Follow above sterilization techniques to do the validation, BI are negative. Repeat twice.

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According to above validation parameters, reduce the sterilization time to be 60/120 mins for sub lethal cycle test, when some of BI showed negative, test is finished.

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Follow above validation parameters to do full cycle test.

验证和结论 Validation and result

2、负载温度/湿度均匀性确认 Temperature/Humidity validation when products is in chamber

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Accords with the requirements of ISO11135-2007, the volume of sterilization chamber is no more than 3 m³, should use 3ea temperature sensors to monitor the sterilization chamber temperature, add a minimum of one temperature sensor for adding 1m³ of sterilization chamber if the chamber volume is more than 3 m³ during the empty loading condition; and it need 3ea temperature sensors at least to monitor the sterilization chamber temperature of 20 m³ sterilization chamber (in fact, we used 20ea temperature sensors), the temperature sensor placement diagram, please see the attachment1.

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Set the sterilization temperature °C, record and analyze data.

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Drawing about distribution of temperature sensors for temperature uniformity validation in autoclave

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Validation of biological performance is the validation of validity of sterilization process with autoclave in sterilization cycle.

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According to ISO11135-1:2007, Atcc7372, which is strong against ETO gas, is used as our biological indicator during the validation of biological performance. The original biological organism content is _____cfu. Before preheated process, according to uniformity distribution rule, we put 40pc BI into the chamber, the distribution drawing is in attached document 2.

根据灭菌负载的产品构造、包装形式及其在灭菌器内的分布方式, 确定微生物性能确认时的相关灭菌工艺如下: According to products structure, package form and distribution in the chamber, we decide the relative sterilization techniques as following:

产品 Product

灭菌工序的灭菌 Sterilization process condition

微生物性能验证工艺参数	
Parameters Specifications/Tolerances/Ranges	
参数 Parameter	参数容限 Specifications/Tolerances/Ranges
预处理 Preconditioning	
1.1 时间 Time	
1.2 预热温度 Temperature in the cell	
1.3 相对湿度 Relative Humidity in the cell	
灭菌 Sterilization	
2.1 灭菌温度 Temperature	
2.2 保温时间 Temperature conditioning time	
2.3 初始真空度 Initial vacuum	
2.4 保压时间及最大允许泄漏量 Maximum allowed leak rate during a set time leak test	
2.5 相对湿度 Relative humidity in the chamber	
2.6 灭菌剂浓度 Theoretical calculated gas Concentration based on validation	
2.7 灭菌剂温度 Gas temperrature at the admission to the chamber	
2.8 灭菌时间 Stirilization Cycle time	
2.9 换气真空度 Pressure at rinsing	
3.0 换气次数 Number of air ringsing steps	
3.1 通风时间 Aeration time	

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sterilization techniques according to the result if BI has bacteria or not after cultivation.
For details please see table 1.

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From the result of table 1, in these sterilization parameters, if the product structure, package form and load distribution keep in the same conditions, the sterilization process is validity. The sterilization critical time for this product (shortest effective sterilization time)

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To summarized, after validating this autoclave (20m³) about installation, operation and performance validation, this autoclave validation process and result are accordance with ISO11135-1:2007.

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Additionally, this result is got from half cycle method; we recommend for product A, the sterilization time is_____mins for total cycle to guarantee the validity sterilization.

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To certify the validity of sterilization process, we use half cycle method; keep above sterilization techniques, change sterilization time, and culture BI which is used in different sterilization time in sterile environment, after that evaluating the validity of sterilization techniques according to the result if BI has bacteria or not after cultivation. For details please see table 2.

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