



U.S. Department of Health & Human Services

Food and Drug Administration

SAVE REQUEST

USER: (szs)
FOLDER: K103184 - 354 pages
COMPANY: HOME SKINOVATIONS LTD. (HOMESKIN)
PRODUCT: LIGHT BASED OVER-THE-COUNTER HAIR REMOVAL (OHT)
SUMMARY: Product: SILK'N FLASH N GO

DATE REQUESTED: Apr 19, 2012

DATE PRINTED: Apr 19, 2012

Note: Printed



K103184

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Home Skinovations Ltd.

Silk'n Flash N Go

NOV 10 2010

This summary of safety and effectiveness information is being submitted in accordance with the requirements of the SMDA 1990 and 21 CFR 807.92.

Submitter's information

Name: Home Skinovations Ltd.
Address: Apolo building, POB 533, Yokneam 20692, Israel
Contact: Dr. Amir Waldman VP Regulatory Affairs

Device information

Trade/Proprietary name: Silk'n Flash N Go
Common/Usual name: Light based hair removal device
Classification name: Laser surgical instrument for use in general and plastic surgery and in dermatology (21CFR §878.4810)
Product code: GEX

Predicate devices

- Flash N Go (K082298), by Home Skinovations Ltd.
- TRIA Laser Hair Removal System (K090820), by Tria Beauty, Inc.
- EpiLight and PhotoDerm HR (K991935), by ESC Medical system Inc.

Intended use:

The Flash N Go device is an over the counter device intended for the removal of unwanted hair. Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

Device Description:

The Flash N Go hair removal system is a pulsed light system composed of a base unit and hand held applicator.

Performance data:

The device complies with the following U.S. Food and Drug Administration performance standards: 21CFR § 1040.10 & 1040.11. Clinical data was collected in a prospective multisite clinical study.

Substantial Equivalence:

The Flash N Go system is substantial equivalent to its predicate devices. The data in this 510k submission demonstrate that the Flash N Go system is identical to the cleared Flash N' Go, and shares the same intended use as other predicate devices. Therefore is substantial equivalent to its predicate devices. Details are provided in Substantial equivalent section of this submission.

Based upon an analysis of the overall performance characteristic for the device, Home Skinovations Ltd. believes that no significant differences exist. Therefore the Silk'n Flash N Go should raise no new issues of safety or effectiveness.

October 21, 2010



Date

Dr. Amir Waldman,
VP Regulatory Affairs
Home Skinovations Ltd.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Home Skinovations Ltd.
% Dr. Amir Waldman
Vice President, Regulatory Affairs
Apolo Building, P.O. Box 533
Yokneam 20692, Israel

NOV 10 2010

Re: K103184

Trade/Device Name: Silk'n Flash N Go
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology
Regulatory Class: Class II
Product Code: OHT
Dated: October 21, 2010
Received: October 29, 2010

Dear Dr. Waldman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson
Director
Division of Surgical, Orthopedic
And Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known) K103184.

Device Name Flash N Go

Indications For Use:

NOV 10 2010

Flash N Go is an over the counter device intended for the removal of unwanted hair.

Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over The Counter Use X

Neil R. O'Brien, MD
(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

(Optional Format 1-2-96)

510(k) Number K103184



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Home Skinovations Ltd.
% Dr. Amir Waldman
Vice President, Regulatory Affairs
Apolo Building, P.O. Box 533
Yokneam 20692, Israel

NOV 10 2010

Re: K103184

Trade/Device Name: Silk'n Flash N Go
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology
Regulatory Class: Class II
Product Code: OHT
Dated: October 21, 2010
Received: October 29, 2010

Dear Dr. Waldman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson
Director
Division of Surgical, Orthopedic
And Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known) K103184

Device Name Flash N Go

Indications For Use:

NOV 10 2010

Flash N Go is an over the counter device intended for the removal of unwanted hair.

Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over The Counter Use X

Neil R. Ogden for MxM
(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

(Optional Format 1-2-96)

510(k) Number K103184



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

U.S. Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center 6 WO66-G609
10903 New Hampshire Avenue
Silver Spring, MD 20993-0002

November 01, 2010

HOME SKINOVATIONS LTD.
APOLO BUILDING POB 533
YOKNEAM ILLIT
ISRAEL 20692
ATTN: AMIR WALDMAN

510k Number: K103184

Received: 10/29/2010

Product: SILK'N FLASH N GO

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification, (510(k)), you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product and for the above referenced 510(k) submitter. Please note, if the 510(k) submitter is incorrect, please notify the 510(k) Staff immediately. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. We will notify you when the processing of your 510(k) has been completed or if any additional information is required. **YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.**

Please remember that all correspondence concerning your submission **MUST** be sent to the Document Mail Center (DMC) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official 510(k) submission.

On September 27, 2007, the President signed an act reauthorizing medical device user fees for fiscal years 2008 - 2012. The legislation - the Medical Device User Fee Amendments of 2007 is part of a larger bill, the Food and Drug Amendments Act of 2007. Please visit our website at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Overview/MedicalDeviceUserFeeandModernizationActMDUFMA/default.htm>

for more information regarding fees and FDA review goals. In addition, effective January 2, 2008, any firm that chooses to use a standard in the review of ANY new 510(k) needs to fill out the new standards form (Form 3654) and submit it with their 510(k). The form may be found at <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm>.

We remind you that Title VIII of the Food and Drug Administration Amendments Act of 2007 (FDAAA) amended the PHS Act by adding new section 402(j) (42 U.S.C. § 282(j)), which expanded the current database known as ClinicalTrials.gov to include mandatory registration and reporting of results for applicable clinical trials of human drugs (including biological products) and devices. Section 402(j) requires that a certification form <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm> accompany 510(k)/HDE/PMA submissions. The agency has issued a draft guidance titled: "Certifications To Accompany Drug, Biological

Product, and Device Applications/Submissions: Compliance with Section 402(j) of The Public Health Service Act, Added By Title VIII of The Food and Drug Administration Amendments Act of 2007"

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm134034.htm>. According to the draft guidance, 510(k) submissions that do not contain clinical data do not need the certification form.

Please note the following documents as they relate to 510(k) review: 1) Guidance for Industry and FDA Staff entitled, "Interactive Review for Medical Device Submissions: 510(k)s, Original PMAs, PMA Supplements, Original BLAs and BLA Supplements". This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>. Please refer to this guidance for information on a formalized interactive review process. 2) Guidance for Industry and FDA Staff entitled, "Format for Traditional and Abbreviated 510(k)s". This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm084365.htm>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRH's e-Copy Program, you may replace one paper copy of any premarket submission (e.g., 510(k), IDE, PMA, HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/ucm134508.html>. In addition, the 510(k) Program Video is now available for viewing on line at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm070201.htm>.

Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.

Lastly, you should be familiar with the regulatory requirements for medical devices available at Device Advice <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have questions on the status of your submission, please contact DSMICA at (301)796-7100 or the toll-free number (800)638-2041, or at their internet address <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have procedural questions, please contact the 510(k) Staff at (301)796-5640.

Sincerely,

510(k) Staff



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

U.S. Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center – WO66-0609
10903 New Hampshire Avenue
Silver Spring, MD 20993-0002

510(k) Submission for Silk'n Flash N Go

Home Skinovations Ltd.

October 21, 2010

K-S

Silk'n Flash N'Go
510(k) Notification

Table of contents

| | | |
|-----|--|-------------|
| 1. | Medical device user fee cover sheet (form FDA 3601) | |
| 2. | CDRH Premarketing Review Submission Cover Sheet | |
| 3. | 510(k) Cover Letter | |
| 4. | Indication For Use Statement | |
| 5. | 510(k) Summary Of Safety and Effectiveness | |
| 6. | Premarket Notification Truthful and Accuracy Statement | |
| 7. | Class III Summary and Certification | NA |
| 8. | Financial Certification or Disclosure Statement | NA |
| 9. | Declaration of Conformity and Summary Reports | NA |
| 10. | Executive Summary | NA |
| 11. | Device Description | Appendix 1* |
| 12. | Substantial Equivalent Discussion | |
| 13. | Proposed labeling | |
| 14. | Sterilization and Shelf Life | NA |
| 15. | Biocompatibility | Appendix 1* |
| 16. | Risk Analysis and Software Documents | Appendix 1* |
| 17. | Electromagnetic Compatibility and Electrical Safety | Appendix 1* |
| 18. | Performance Testing - Bench | Appendix 1* |
| 19. | Performance Testing - Animal | NA |
| 20. | Performance Testing - Clinical | |

*Appendix 1 contains the Flash N Go (cleared) submission and all communication with the FDA.

Chapter 1

Medical device user fee cover sheet (form FDA 3601)

Form Approved: OMB No. 0910-511 Expiration Date: January 31, 2010. See Instructions for OMB Statement.

| | | | | | | | |
|--|--|--|--|--|--|---|---|
| DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET | | PAYMENT IDENTIFICATION NUMBER: MD6052074-956733 Write the Payment Identification number on your check. | | | | | |
| A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: http://www.fda.gov/oc/mdufma/coversheet.html | | | | | | | |
| 1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code) Home Skinovations Ltd. Apolo Buld POB 533 Yokneam 20692 IL 1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) | | 2. CONTACT NAME Yoram Sadeh 2.1 E-MAIL ADDRESS yorams@silkn.com 2.2 TELEPHONE NUMBER (include Area code) 972-49097476 2.3 FACSIMILE (FAX) NUMBER (include Area code) | | | | | |
| 3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: http://www.fda.gov/oc/mdufma) <table border="0"> <tr> <td> <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice </td> <td> 3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP) </td> </tr> </table> | | | | <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice | 3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP) | | |
| <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice | 3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP) | | | | | | |
| 4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business 4.1 If Yes, please enter your Small Business Decision Number: | | | | | | | |
| 5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see http://www.fda.gov/cdrh/mdufma for additional information) | | | | | | | |
| 6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <table border="0"> <tr> <td><input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates</td> <td><input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population</td> </tr> <tr> <td><input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only</td> <td><input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially</td> </tr> </table> | | | | <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates | <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population | <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only | <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially |
| <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates | <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population | | | | | | |
| <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only | <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially | | | | | | |
| 7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA)). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO | | | | | | | |
| 8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION \$4,348.00 | | | | | | | |

Form FDA 3601 (01/2007)

["Close Window"](#) [Print Cover sheet](#)

Chapter 2

CDRH Premarketing Review Submission Cover Sheet

| | | | |
|--|--|---|--|
| DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION CDRH PREMARKET REVIEW SUBMISSION COVER SHEET | | Form Approval OMB No. 0910-0120 Expiration Date: August 31, 2010. See OMB Statement on page 5. | |
| Date of Submission October 21, 2010 | | User Fee Payment ID Number MD6052074-956733 | |
| FDA Submission Document Number (if known) | | | |

| SECTION A TYPE OF SUBMISSION | | | | |
|--|--|---|--|--|
| PMA <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement | PMA & HDE Supplement <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other | PDP <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP | 510(k) <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party | Meeting <input type="checkbox"/> Pre-510(K) Meeting <input type="checkbox"/> Pre-IDE Meeting <input type="checkbox"/> Pre-PMA Meeting <input type="checkbox"/> Pre-PDP Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Other (specify): |
| IDE <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement | Humanitarian Device Exemption (HDE) <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment | Class II Exemption Petition <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information | Evaluation of Automatic Class III Designation (De Novo) <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information | Other Submission <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission): |

Have you used or cited Standards in your submission? ☐ Yes ☐ No (If Yes, please complete Section I, Page 5)

| SECTION B SUBMITTER, APPLICANT OR SPONSOR | | | |
|--|------------------|--|-------------------|
| Company / Institution Name Home Skinovations Ltd. | | Establishment Registration Number (if known) 3006948334 | |
| Division Name (if applicable) | | Phone Number (including area code) +972-547404616 | |
| Street Address Apolo building POB 533 | | FAX Number (including area code) +972-49097471 | |
| City Yokneam Illit | State / Province | ZIP/Postal Code 20692 | Country ISRAEL |
| Contact Name Amir Waldman Ph.D. | | | |
| Contact Title VP regulatory affairs | | Contact E-mail Address amirw@silkn.com | |

| SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above) | | | |
|---|------------------|------------------------------------|---------|
| Company / Institution Name | | | |
| Division Name (if applicable) | | Phone Number (including area code) | |
| Street Address | | FAX Number (including area code) | |
| City | State / Province | ZIP Code | Country |
| Contact Name | | | |
| Contact Title | | Contact E-mail Address | |

| SECTION D1 | | | REASON FOR APPLICATION - PMA, PDP, OR HDE | | |
|--|---|---|---|--|--|
| <input type="checkbox"/> New Device <input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site | <input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software / Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other (specify below) | <input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager <input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment | | | |
| <input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Packaging <input type="checkbox"/> Sterilization <input type="checkbox"/> Other (specify below) | <input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance Characteristics <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other (specify below) | <input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address | | | |
| <input type="checkbox"/> Response to FDA correspondence: | | | | | |
| <input type="checkbox"/> Other Reason (specify): | | | | | |
| SECTION D2 | | | REASON FOR APPLICATION - IDE | | |
| <input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access | <input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent / Applicant <input type="checkbox"/> Design / Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor <input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final | <input type="checkbox"/> Response to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing | | | |
| <input type="checkbox"/> Other Reason (specify): | | | | | |
| SECTION D3 | | | REASON FOR SUBMISSION - 510(k) | | |
| <input checked="" type="checkbox"/> New Device | <input type="checkbox"/> Additional or Expanded Indications | <input type="checkbox"/> Change in Technology | | | |
| <input type="checkbox"/> Other Reason (specify): | | | | | |

SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS

| | | | | | |
|--|-----|---|--|--|--|
| Product codes of devices to which substantial equivalence is claimed | | | | Summary of, or statement concerning, safety and effectiveness information | |
| 1 | GEX | 2 | | 3 | |
| 5 | | 6 | | 7 | |
| | | | | <input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement | |

Information on devices to which substantial equivalence is claimed (if known)

| | 510(k) Number | Trade or Proprietary or Model Name | Manufacturer |
|---|---------------|------------------------------------|-------------------------|
| 1 | K082298 | Flash N Go | Home Skinovations Ltd. |
| 2 | K090820 | TRIA Laser Hair Removal System | Tria Beauty Inc. |
| 3 | K991935 | EpiLight and PhotoDerm HR | ESC Medical system Inc. |
| 4 | | | |
| 5 | | | |
| 6 | | | |

SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS

Common or usual name or classification name

Laser Instruments, Surgical powered- general and Plastic Surgery - Class II, 79 GEX

| | Trade or Proprietary or Model Name for This Device | Model Number |
|---|--|--------------|
| 1 | Silk'n Flash N Go | 1 |
| 2 | | 2 |
| 3 | | 3 |
| 4 | | 4 |
| 5 | | 5 |

FDA document numbers of all prior related submissions (regardless of outcome)

| | | | | | |
|---|---|---|----|----|----|
| 1 | 2 | 3 | 4 | 5 | 6 |
| 7 | 8 | 9 | 10 | 11 | 12 |

Data Included in Submission

☐ Laboratory Testing☐ Animal Trials☐ Human Trials**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

| | | |
|----------------------|--------------------------------|---|
| Product Code | C.F.R. Section (if applicable) | Device Class |
| GEX | | <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified |
| Classification Panel | | |

Indications (from labeling)

The Flash N' Go device is an over the counter device intended for the removal of unwanted hair.
Flash N' Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

Note: Submission of this information does not affect the need to submit a 2891 or 2891a Device Establishment Registration form.

FDA Document Number (if known)

SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION

| | | | |
|--|--|--|---|
| <input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete | Facility Establishment Identifier (FEI) Number | <input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer | <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler |
| Company / Institution Name Home Skinovations Ltd. | | Establishment Registration Number 3006948334 | |
| Division Name (if applicable) | | Phone Number (including area code) +972-547404616 | |
| Street Address Apolo Buld. POB 533 | | FAX Number (including area code) +972-49097471 | |
| City Yokneam Illit | | State / Province | ZIP Code 20692 |
| | | Country Israel | |
| Contact Name Amir Waldman Ph.D. | Contact Title VP Regulatory affairs | Contact E-mail Address amirw@silkn.com | |

| | | | |
|---|--|---|---|
| <input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete | Facility Establishment Identifier (FEI) Number | <input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer | <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler |
| Company / Institution Name | | Establishment Registration Number | |
| Division Name (if applicable) | | Phone Number (including area code) | |
| Street Address | | FAX Number (including area code) | |
| City | | State / Province | ZIP Code |
| | | Country | |
| Contact Name | Contact Title | Contact E-mail Address | |

| | | | |
|---|--|---|---|
| <input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete | Facility Establishment Identifier (FEI) Number | <input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer | <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler |
| Company / Institution Name | | Establishment Registration Number | |
| Division Name (if applicable) | | Phone Number (including area code) | |
| Street Address | | FAX Number (including area code) | |
| City | | State / Province | ZIP Code |
| | | Country | |
| Contact Name | Contact Title | Contact E-mail Address | |

SECTION I**UTILIZATION OF STANDARDS**

Note: Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

| | | | | | |
|---|--------------------------------|---|---|------------------------|--------------------|
| 1 | Standards No. IEC 60601-1 | Standards Organization International Electrotechnical Commission | Standards Title Medical Electrical Equipment" Part 1: "General requirements for safety" | Version 2 edition | Date 01/01/1988 |
| 2 | Standards No. IEC 60601-1-2 | Standards Organization International Electrotechnical Commission | Standards Title Medical Electrical Equipment part 1-2 collateral standard electromagnetic compatibility - requirement and tests | Version 2.1 edition | Date 01/01/2004 |
| 3 | Standards No. | Standards Organization | Standards Title | Version | Date |
| 4 | Standards No. | Standards Organization | Standards Title | Version | Date |
| 5 | Standards No. | Standards Organization | Standards Title | Version | Date |
| 6 | Standards No. | Standards Organization | Standards Title | Version | Date |
| 7 | Standards No. | Standards Organization | Standards Title | Version | Date |

Please include any additional standards to be cited on a separate page.

Public reporting burden for this collection of information is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Department of Health and Human Services
Food and Drug Administration
Office of the Chief Information Officer (HFA-710)
5600 Fishers Lane
Rockville, Maryland 20857

An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.

Chapter 3
510(k) Cover Letter



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

October 21, 2010

U.S. Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center – WO66-0609
10903 New Hampshire Avenue
Silver Spring, MD 20993-0002

FDA CDRH DMC

OCT 29 2010

Received

V-S

Subject: 510(k) Submission for Flash N Go

Dear Sir/Madam,

Home Skinovations Ltd. is submitting this premarket notification for the purpose of obtaining FDA marketing clearance for its product. A table of contents for this 510(k) application is presented immediately after this cover letter. The following is information required under 21 CFR 807.87:

This is to notify you of the intention of Home Skinovations Ltd. to manufacture and market the following device:

Please send all correspondence by fax (+972-4-9097471), e-mail amirw@silkn.com in addition to mail

| | |
|--------------------------------|--|
| Trade/Proprietary name: | Flash N Go |
| Common name: | Light based hair removal system |
| Classification name: | The subject of this application is a Laser surgical instrument for use in general and plastic surgery and in dermatology (21 CFR 878.4810) |
| Product codes: | GEX |
| Class of device: | The device should be considered a class II device, as are the predicate devices. |
| Reason for Application: | New device. |



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

| | |
|---|---|
| Name and address of | Home Skinovations Ltd. |
| Manufacturer: | Home Skinovations Ltd. Apolo building POB 533 Yokneam 20692 ISRAEL Tel: +972(4)9097440 Fax: +972(4)9097471 |
| Contact person: | Dr. Amir Waldman VP regulatory affairs Tel. +972-547-404616, +972-4-9097440 Fax +972-4-9097471 E-mail - amirw@silkn.com |
| Applicable standards: | Action taken by the person required to register to comply with the requirements of the Act under section 514 for performance standards or section 513 for special controls: There are no mandatory performance standards applicable to this device. There are voluntary performance standards that the Flash N' Go complies with as discussed in chapter 4. |
| Labeling/Promotional material: | The labeling and promotional materials for the Flash N' Go, and for the substantially equivalent predicate devices, are presented in later sections of this application. |
| Substantial equivalence predicate devices: | Flash N' Go is substantially equivalent to the following devices: <ul style="list-style-type: none">• Flash N' Go: The Flash N' Go is an over the counter device intended for the removal of unwanted hair. Manufactured by Home Skinovations Ltd. and subject of K082298.• TRIA Laser Hair Removal System: TRIA is an over-the-counter device intended for adjunctive |

use with shaving for hair removal sustained with periodic treatments. TRIA is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime. Manufactured by Tria Beauty, Inc. and subject of K090820.

- EpiLight and PhotoDerm HR: EpiLight and Photoderm HR are used for removal of unwanted hair. EpiLight and PhotoDerm HR are also intended to effect stable long-term, or permanent hair reduction in skin types I-V through selective targeting of melanin in hair follicle. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs regrowing after a treatment regime. Manufactured by ESC Medical system Inc. and subject of K991935.

Confidentiality: The material in this application is considered confidential and is not meant for public release.

Truth and Accuracy Statement: A truthful and accurate statement, according to 21CFR 807.87(k), is presented after a separate page with the indications for use.

SMDA Safety and Effectiveness Summary: Summary of Safety and is presented following the table of contents.

The Flash N' Go is a pulsed light device (b) (4) cleared Flash N' Go device K082298, and has technological characteristics that are substantially equivalent to other predicate devices. Therefore all safety features are identical to the cleared Flash N' Go device. The indication for use are identical to the other 2 cleared devices, and (b) (4) (b) (4).



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

We trust that the information included in this application will be adequate to allow for its prompt review. Please contact me if there are any questions.

Sincerely yours,

Dr. Amir Waldman
VP regulatory affairs
Home Skinovations Ltd.

A handwritten signature in cursive script, appearing to read "Amir Waldman", is written over the typed name.

Chapter 4
Indication For Use Statement

Chapter 5

510(k) Summary of Safety and Effectiveness

510(k) Number (if known)_____.

Device Name Flash N Go

Indications For Use:

Flash N Go is an over the counter device intended for the removal of unwanted hair.

Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use_____

(Per 21 CFR 801.109)

OR

Over The Counter Use X

(Optional Format 1-2-96)

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Home Skinovations Ltd.

Silk'n Flash N Go

This summary of safety and effectiveness information is being submitted in accordance with the requirements of the SMDA 1990 and 21 CFR 807.92.

Submitter's information

Name: Home Skinovations Ltd.
Address: Apolo building, POB 533, Yokneam 20692, Israel
Contact: Dr. Amir Waldman VP Regulatory Affairs

Device information

Trade/Proprietary name: Silk'n Flash N Go
Common/Usual name: Light based hair removal device
Classification name: Laser surgical instrument for use in general and plastic surgery and in dermatology (21CFR §878.4810)
Product code: GEX

Predicate devices

- Flash N Go (K082298), by Home Skinovations Ltd.
- TRIA Laser Hair Removal System (K090820), by Tria Beauty, Inc.
- EpiLight and PhotoDerm HR (K991935), by ESC Medical system Inc.

Intended use:

The Flash N Go device is an over the counter device intended for the removal of unwanted hair. Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

Device Description:

The Flash N Go hair removal system is a pulsed light system composed of a base unit and hand held applicator.

Performance data:

The device complies with the following U.S. Food and Drug Administration performance standards: 21CFR § 1040.10 & 1040.11. Clinical data was collected in a prospective multisite clinical study.

Substantial Equivalence:

The Flash N Go system is substantial equivalent to its predicate devices. The data in this 510k submission demonstrate that the Flash N Go system is identical to the cleared Flash N' Go, and shares the same intended use as other predicate devices. Therefore is substantial equivalent to its predicate devices. Details are provided in Substantial equivalent section of this submission.

Based upon an analysis of the overall performance characteristic for the device, Home Skinovations Ltd. believes that no significant differences exist. Therefore the Silk'n Flash N Go should raise no new issues of safety or effectiveness.

October 21, 2010

Date



Dr. Amir Waldman,
VP Regulatory Affairs
Home Skinovations Ltd.

Chapter 6

Premarket Notification Truthful and Accuracy Statement

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS
PREMARKET NOTIFICATION TRUTHFUL AND
ACCURACY STATEMENT**

(As required by 21 CFR 807.87(j))

I certify that in my capacity as VP of Regulatory Affairs of Home Skinovations Ltd., I believe, to the best of my knowledge that all data and information submitted in this premarket notification are truthful and accurate and that no material fact has been omitted.



Dr. Amir Waldman

October 21, 2010

Chapter 12
Substantial Equivalent Discussion

Substantial Equivalence Discussion

Flash N' Go indication for use: The Flash N Go device is an over the counter device intended for the removal of unwanted hair. Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

The Flash N Go is substantially equivalent (SE) to the following devices:

- Flash N Go: The Flash N' Go is an over the counter device intended for the removal of unwanted hair. Manufactured by Home Skinovations Ltd. and subject of K082298.
- TRIA Laser Hair Removal System: TRIA is an over-the-counter device intended for adjunctive use with shaving for hair removal sustained with periodic treatments. TRIA is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime. Manufactured by Tria Beauty, Inc. and subject of K090820.
- EpiLight and PhotoDerm HR: EpiLight and Photoderm HR are used for removal of unwanted hair. EpiLight and PhotoDerm HR are also intended to effect stable long-term, or permanent hair reduction in skin types I-V through selective targeting of melanin in hair follicle. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs regrowing after a treatment regime. Manufactured by ESC Medical system Inc. and subject of K991935.

Labeling and promotional material of the substantially equivalent devices are presented in Appendices 5-1, 5-2, 5-3 of this chapter.

Substantial equivalence is demonstrated according to the criteria outlined in the 510(k) "Substantial Equivalence" Decision Making Process, as published in the FDA 510(k) Guidance (FDA 90-4158).

A comparison table of the subject and predicate devices is included in this discussion.

Technological equivalence

The Flash N Go is a pulsed light device (b) (4)

(b) (4) and has technological characteristics that are substantially equivalent to other predicate devices. (b) (4)

EpiLight and PhotoDerm HR are intense pulsed light devices. Tria hair removal laser is a laser with wavelength of 810 nm.

A comparison is provided in the table below:

| Technological Characteristic | Flash N Go(new) | Flash N Go(cleared) | Tria | EpiLight / PhotoDerm HR | Comments |
|---|--------------------------------------|--------------------------------------|--------------------|--------------------------------------|--|
| Nature and Source of Energy Imparted to Patient | Pulsed light via external applicator | Pulsed light via external applicator | Pulsed Diode laser | Pulsed light via external applicator | All devices use optical energy as the energy source for hair removal treatment |
| Wavelength range (nm) | 475-1200 | 475-1200 | 810 | 645-1200 | All devices use Visible / Near IR Wavelengths range to achieve reduced hair growth |
| Maximum Light Energy (J/cm ²) | 5 | 5 | 20 | 45 | Same output parameters as the cleared Flash N' Go |

No new types of safety or effectiveness questions are raised and differences in the technological characteristics between the Flash N Go and the predicate devices can be assessed for their effects on safety and effectiveness.

Intended use equivalence

The Flash N Go device is an over the counter device intended for the removal of unwanted hair. Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

The Flash N Go employs pulses of low energy light to achieve permanent hair reduction. (b) (4)

(b) (4)

The Flash N Go (new) intended use is identical to the Tria, EpiLight and PhotoDerm HR. Both the Flash N Go and the Tria devices are over the counter device: "TRIA is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime". While the EpiLight and PhotoDerm HR are prescription use, but have the same intended use.

Safety and Effectiveness equivalence


The Flash N Go does not present any new safety and effectiveness issues as compared with the predicate devices.

(b) (4)

As far as effectiveness, since the energy pulse duration as well as spectrum is within similar range to the predicate devices efficacy is expected to be similar. (b) (4)

(b) (4)

In light of the results, we conclude that the Flash N Go device is safe and effective as an over the counter device indicated for the removal of unwanted hair and also for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.



enter email for exclusive offers [sign up](#)

NOW SHIPPING online activation my account shopping basket (0)

acne & breakouts see the results news & events what is TRIA quick shop

hair removal

now save \$100—new lower prices



TRIA Laser Hair Removal System

stop shaving. start showing off.

Imagine waking up to silky smooth skin every day and being free forever from endless shaving and painful waxing! The TRIA Laser is the only way to get permanent hair removal at home—guaranteed.

- See results in 3 months; in as little as 6 months, you're done! [How To Use](#)
- Get baby soft skin—like the hair was never there
- Created by the same scientists that invented diode laser hair removal for dermatologists
- Get the same results as professional laser hair removal and save thousands of dollars
- FDA-cleared safe, effective, permanent

See what [TRIA users](#) say and [read the buzz](#) on why *Allure*, *InStyle*, *Women's Wear Daily* and *The Doctor's* have named the TRIA Laser a beauty "must-have".

Find out if the TRIA Laser is [right for you](#).

NOW SHIPPING

TRIA Laser Hair Removal System \$495 (was \$595)
new lower price!

[add to basket](#)

☐ half price today with [easypay](#) | \$249 down and six payments of \$49. [details](#)

share: [| more](#)

you may also like:



TRIA Skin Clarifying System, \$295



Stay Clear 3-Piece Refill Kit, \$80



Super Clear 2-Piece Refill Kit, \$60

how to use

what to expect

is it right for you?


included



3 simple steps. amazing results.

step 1 : get ready

Cleanse, shave, and towel-dry the skin.



step 2 : get comfortable

Select the highest setting that suits your personal comfort level.



step 3 : treat

Apply the tip of the laser to the skin and listen for the beep. Repeat until you have covered the entire area to be treated. Treat twice a month for the first 3 months, then once per month for the next several months until the hair stops coming back. Bikini line and underarms take 5-10 minutes to treat. Legs may take 20-30 minutes each.

*For detailed instructions check out our [step-by-step How To Use Video](#).

Questions? [Contact Us](#)

stay connected: [f](#) [t](#) [v](#)


1-877-321-TRIA (8742) | customer care

[about us](#) | [store locator](#) | [contact us](#) | [affiliate program](#) | [japan](#) | [korea](#)

© 2010 TRIA Beauty, Inc.

Flash N' GoTM

Personal Hair Removal at Home
with Home Pulsed LightTM




Flash N' Go™ is a physician recommended and innovative light-based device for hair removal in the privacy of your home. Safe, easy to use and cost effective, **Flash N' Go™** achieves marked clinical results while offering personal convenience.

THE HPL™ INNOVATION


The Flash N' Go™ device utilizes the highly sophisticated and new **Home Pulsed Light™ (HPL™)** technology, developed by Home Skinovations Ltd. especially for safe use at home.

HPL™ Epilation is based on the theory of selective photothermolysis in which **optical** energy is used to disable hair growth. Epilation with Home Skinovations' proprietary HPL™ technology is further benefited by a unique **acoustic** effect that enhances the normal process of epilation by photothermolysis.


To safely achieve the desired clinical results the hair shaft needs to selectively absorb HPL™ light energy and transform it into heat. Melanin is the pigment in the hair shaft that is responsible for the absorption of light, which generates the heat at the follicle. Selectivity is achieved when the HPL™ optical energy delivered to the tissue is mostly absorbed by pigment in the hair shaft, while the epidermis and the surrounding tissue remain cool. Additionally HPL™ technology creates an acoustic effect that enhances the overall hair removal result. Together the optical and acoustic effects of HPL™ technology are what enable the Flash N' Go™ device to be safely used by patients for home based hair removal.



● Before using Flash N' Go™



● During Flash N' Go™ pulse



● After Flash N' Go™ application

Flash N' Go™ HAIR REMOVAL CLINICAL STUDY

Flash N' Go™ has been tested in a multi-center clinical study that took place in North America and Israel and involved over 150 female patients.

Treatment Protocol:

| | |
|-----------------|--|
| Method: | Self treatment at the physician's office |
| Treatment plan: | 3 treatments performed at 2 week intervals |
| Hair count: | Initial hair count and follow up hair counts 3 and 6 months post last treatment |
| Investigators: | Dr. Tina Alster, Washington DC. Dr. Stephen Mulholland, Los Angeles Dr. James Shaoul, Tel Aviv |

Result

| Body site | No. of patients | Average reduction after 3 months | Average reduction after 6 months |
|-----------|-----------------|----------------------------------|----------------------------------|
| Axilla | 65 | 54% | 41% |
| Legs | 34 | 65% | 54% |
| Bikini | 56 | 56% | 43% |
| Arms | 20 | 58% | 52% |

WHAT EXPERTS ARE SAYING ABOUT Flash N' Go™

"As a leading investigator in the development of state-of-the-art lasers and IPL systems for hair removal, I was surprised by the marked clinical results that were achieved with Flash N' Go™. Patients who were previously treated with professional hair removal systems used in my clinic had a chance to use Flash N' Go™ in clinical trials. They reported that self-treatment was easy, relatively painless, and resulted in noticeable hair reduction after each treatment."

Tina S. Alster, M.D.

Director, Washington Institute of Dermatologic Laser Surgery
Clinical Professor of Dermatology, Georgetown University Medical Center

"After supporting the development of leading professional IPL and laser systems for hair removal, and using them in my clinic for many years, I introduced Flash N' Go™ to my patients. After 4 to 6 treatments with the Flash N' Go™ device they were extremely satisfied to see dramatic hair reduction."


R. Stephen Mulholland, MD FRCS(C)

Cosmetic Plastic Surgeon

TYPICAL RESULTS OF Flash N' Go™



● Before



● After

Flash N' Go™ ADVANTAGES

● Physician Recommended

Tested over 12 months in a multi-center clinical study involving over 150 female patients.

● Small and Portable

Flash N' Go™ is about the size of a hairdryer and it runs off of regular electric power, so treatment can be done anywhere, anytime and non-stop for as long as necessary.

● Safe

When used according to instructions, hair removal with Flash N' Go™ is safe and with minimal discomfort compared to waxing or laser.

● Easy to Use

Flash N' Go™ is easy to set up, and treatments are simple, clean and very quick.

● Cost Effective

Using Flash N' Go™ will save you money when compared to extended use of professional laser treatments, waxing or shaving.

● Marked clinical Results

Although small and portable, the HPL™ technology in Flash N' Go™ captures all the clinical benefits of light-based hair removal.

TECHNICAL SPECIFICATIONS

| | |
|---------------------------------|--|
| Spot size | 2cm x 3cm [6cm ²] |
| Speed | 1 pulse every 3.5 second: 1.7 cm ² /sec |
| Technology | Home Pulsed Light™ |
| Max Energy Level | 5 J |
| Wavelength | 475-1200 nm |
| Charging time / Power source | Continuous operation |
| Time needed to treat lower legs | 30 minutes |
| Operation and safety | Safety sensor tip enables maximum safety & control. Easy to use. |

CONTACT INFORMATION

Home Skinovations Ltd.

Apolo Building, Shaar Yokneam, POB 533 Yokneam 20692 Israel

Tel. +972-4-9097470, Fax +972-4-9097471


info@silkn.com | www.silkn.com

ISO 13485



IPL[®] Quantum[™]

Intense Pulsed Light[™] and Laser Technology



Expandable...To Grow With Your Practice


One System: A multitude of applications

Build the system to suit your practice needs— today and into the future


Whether it's a full-time system for skin treatments using photorejuvenation or a single system that can treat a multitude of applications, the IPL Quantum allows you the flexibility to build the system that best meets the needs of your practice today, and will expand to meet your emerging requirements.

| Clinical Indication | SR | DL | HR | QS |
|---|-----|--------------|-----|-------------------------|
| | IPL | Nd:YAG laser | IPL | Q-Switched Nd:YAG laser |
| Rosacea | ✓ | | | |
| Telangiectasia | ✓ | ✓ | | |
| Vascular lesions | ✓ | ✓ | | |
| Hemangiomas | ✓ | ✓ | | |
| Port wine stains | ✓ | ✓ | | |
| Leg veins | ✓ | ✓ | | |
| Solar lentigines | ✓ | | | ✓ |
| Epidermal pigmented lesions | ✓ | | | ✓ |
| Dermal pigmented lesions | ✓ | | | ✓ |
| Poikiloderma of Civatte | ✓ | | | |
| Skin treatments using photorejuvenation | ✓ | | | |
| Hair removal and permanent hair reduction | | | ✓ | |
| Tattoo removal | | | | ✓ |

A winning combination: Multi-Synchronized Pulsing™ using a broad spectrum of light



High peak-power, Multi-Synchronized Pulsing™ provides optimum treatment of vascular and pigmented lesions while sparing the skin.




Patented advanced technology provides a unique high-powered, broad spectrum of light. The result is the flexibility of IPL with the intensity and optimized results only possible with high peak powers—usually afforded only by lasers.

From the Pioneers...

IPL technology proven in over 10 years of clinical use

Clinically proven in over 10 years of use, Lumenis pioneered and perfected IPL. While some competitive products are similar to our first generation, Lumenis IPL has evolved well beyond these technologies. Lumenis provides the most advanced technology and the confidence of over 10 years of clinical efficacy documented in numerous peer-reviewed papers.

True Versatility



"IPL Photorejuvenation using the Lumenis product has provided my patients the lasting results that they are looking for—with a series of treatments, there is minimal to no 'down-time'. The IPL devices simultaneously and consistently reduce unwanted redness and brown spots, and can even improve skin wrinkling."


—Michael Kulick, M.D.

IPL Quantum SR IPL Skin Treatments Using Photorejuvenation



Our unique IPL photorejuvenation head provides the highest peak-power, most advanced Multi-Syncronized Polses. Engineered specifically for the highest level of safety and efficacy, it offers your patients a remarkable treatment for telangiectasias, broken capillaries and benign brown pigmentation from sun damage and photaging—treatments to restore a youthful texture and complexion to the entire face.


- Fully utilizes advantages of IPL for treating age spots such as sun induced freckles, solar lentigines and telangiectasia—dramatically improving the appearance of photaged skin
- Provides a smoother appearance in 4-6 treatments over four months
- No downtime—patients can resume normal activities immediately
- Treats the entire face for an even and satisfying effect
- Contact cooling for maximum safety and patient comfort



Three months post 3 IPL treatments
Robert A. Weiss, MD


Three IPL treatments over 10 weeks
Robert A. Weiss, MD

IPL Quantum HR Upgrade IPL Photoepilation



The IPL Photoepilation head uses lamp configurations for outputs optimized in the red-to-near infrared part of the spectrum. It offers your patients a proven method of achieving permanent hair reduction by treating all skin types including the dark skin type VI, and nearly every hair color.


- Adds the versatility of permanent hair reduction to the IPL Quantum system
- Multiwavelength and broad pulsewidths allow safe and effective hair removal on all skin types I-VI and pseudofolliculitis barbae (PFB)
- Treats all body areas, hair depths and a broad spectrum of colors
- High speed, large spot-size system supports high patient throughput



One year post 5 IPL treatments
Franklin Johnson, M.D.

Nine months post 6 IPL treatments
Franklin Johnson, M.D.

IPL Quantum DL 1064 nm Nd:YAG Laser to Treat Leg Veins and Benign Cutaneous Vascular Lesions



The 1064 nm laser light has been clinically proven to provide effective treatment of deep and larger telangiectasias, as well as benign cutaneous vascular lesions such as hemangiomas, and melasma.


- Proven, highly effective Nd:YAG technology for the treatment of leg veins up to 4.0 mm and contact cooling for maximum safety and comfort
- Multi-Syncronized Pulsing for optimal efficacy while protecting the skin



Three weeks post 1 Laser treatment
Neil Sauter, M.D., F.A.C.P.


Three weeks post 3 IPL and 1 Laser treatment
Adam Steiner, M.D.

IPL Quantum QS Upgrade Q-Switched 1064 nm Nd:YAG Laser to Treat Tattoos and Benign Pigmented Lesions



Ultra-short, high peak-power nanosecond pulses fragment pigment particles allowing safe removal by the body's natural healing processes. IPL Quantum QS effectively treats dark tattoos and a wide range of benign pigmented epidermal lesions.

- State-of-the-art laser technology for the safe and effective removal of dark tattoos and benign pigmented lesions
- Added versatility for the growing laser tattoo-removal market



Four weeks post 6 Q-Switched Nd:YAG treatments
M. Capoliccio, M.D.

Six weeks post 1 Q-Switched Nd:YAG treatment
M. Capoliccio, M.D.

IPL Quantum™

Intense Pulsed Light™ and Laser Technology

| Technical Specifications | SR | DL | HR | QS |
|--------------------------|---|--------------|----------------------------|----------------------|
| Light Source | IPL | Nd:YAG laser | IPL | Nd:YAG laser |
| Standard Spectrum | 515–1200 nm | 1064 nm | 695–1200 nm | 1064 nm |
| Optional Spectra | 590–1200 nm 640–1200 nm | | 645–1200 nm 755–1200 nm | — |
| Fluence | 15–45 J/cm² | 90–150 J/cm² | 20–45 J/cm² | 2–12.7 J/cm² |
| Pulse Duration | 6–26 ms | 5–38 ms | 6–18 ms | 6–8 ns |
| Pulse Delays | 5–60 ms | 5–100 ms | 10–150 ms | — |
| Spot Size | 34 x 8 mm | 6 mm | 34 x 8 mm | 2, 2.5, 3.5 and 5 mm |
| Repetition Rate | 0.5 Hz | 0.5 Hz | 0.5 Hz | Up to 5 Hz |
| Integrated Skin Cooler | Yes | Yes | Yes | No |
| Electrical Requirements | 200–240 VAC +/- 10%, 16A, 50–60 Hz. Dedicated line. | | | |
| Physical Dimensions | 40 x 40 x 100 cm / 16 in x 16 in x 39 in (WDH) | | | |
| Weight | 75 kg / 165 lbs | | | |

LUMENIS CUSTOMER CARE

Worry-Free Service and Support

At Lumenis we understand the importance of delivering service and technical support of the highest quality. When you buy a laser or light-based system from Lumenis, you make an investment in the world's most innovative medical technology in its category backed by the most comprehensive service and support organization in the industry. Customer satisfaction is paramount and the ultimate goal of our business.

Lumenis Global Sales Offices

Lumenis Inc.
2400 Condensa Street
Santa Clara, CA 95051 USA
Tel + 408 764 3000
877 586 3647
Fax + 408 764 3999
800 505 1133

Lumenis (GERMANY) GmbH
Dreieich-Dreieichenhain,
Germany
Tel + 49 6103 8335 0
Fax + 49 6103 8335 300

Lumenis (FRANCE) SARL
Orsay, France
Tel + 33 1 69 33 14 20
Fax + 33 1 60 19 57 23

Lumenis (ITALY) Srl.
Formello (RM), Italy
Tel + 39 06 90 75 230
Fax + 39 06 90 75 269

Lumenis (UK) Ltd.
London, UK
Tel + 44 20 8324 4200
Fax + 44 20 8324 4222

Lumenis (JAPAN) Co. Ltd.
Tokyo, Japan
Tel + 81 3 5789 8300
Fax + 81 3 5789 8310

Lumenis (CHINA) Ltd.
Beijing, China
Tel + 86 10 6510 2620
Fax + 86 10 6510 2621

Lumenis (HK) Ltd.
Hong Kong
Tel + 852 2174 2800
Fax + 852 2722 5151

Internet: www.lumenis.com

Email: information@lumenis.com



Caution: U.S. law restricts this device to sale by or on the order of a physician.

©2005, the Lumenis group of companies. Lumenis, its logo, Quantum, IPL and Multi-Synchronized Pulsing are trademarks or registered trademarks of the Lumenis group of companies. Product specifications are subject to change without notice. All rights reserved. Printed in the U.S.A. PB000137 Rev. C

Chapter 13

Proposed labeling

The only change in this labeling compared with the previously cleared Flash N Go is the new indication for use on page 4.



Flash N' Go User manual

Copyright © 2008 Home Skinovations Ltd. All rights reserved.

Print date: February 2008

Home Skinovations Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd. is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd.

Data is subject to change without notification.

Home Skinovations Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd.

Specifications are subject to change without notice.

Home Skinovations Ltd.
Apolo Building, Shaar Yokneam
POB 533
Yokneam 20692
ISRAEL

Website: www.silkn.com
e-mail: info@silkn.com

Table of Contents

| | |
|--|-----------|
| 1. What is Flash N' Go™? | 3 |
| 1.1. Flash N' Go™ device description | 3 |
| 1.2. Package contents | 4 |
| 1.3. The Disposable Lamp Cartridge | 4 |
| .2 Flash N' Go™ Intended Use | 4 |
| .3 Safety with Flash N' Go™ | 5 |
| 3.1. With Flash N' Go™ Safety Comes First | 5 |
| 4. Contraindications | 5 |
| 5. When to avoid using Flash N' Go™? | 7 |
| 6. Precautions – How to use Flash N' Go™ Safely | 8 |
| 7. Reducing the risk of injury | 9 |
| 8. Risks of using Flash N' Go Flash N' Go™ | 10 |
| 9. Long Term Hair Removal the Flash N' Go™ Way | 13 |
| 9.1. How does light remove hair? | 13 |
| 9.2. How does the hair growth cycle impact light-based hair removal? | 13 |
| 9.3. Plan your Flash N' Go™ hair removal for best results | 14 |
| 10. How To Use Flash N' Go™ | 16 |
| 10.1. Flash N' Go™ Setup | 16 |
| 10.2. Treating with Flash N' Go™ for the first time | 16 |
| 10.3. What to Expect when treating with Flash N' Go™? | 18 |
| 10.4. After treating with Flash N' Go™ | 19 |
| 11. Maintenance of Flash N' Go™ | 20 |
| 11.1. Cleaning Flash N' Go™ | 20 |
| 11.2. Replacing the LAMP CARTRIDGE | 20 |
| 12. Troubleshooting | 22 |
| 12.1. "My Flash N' Go™ does not start." | 22 |
| 12.2. "A light pulse is not emitted when I press the PULSE SWITCH." | 22 |
| 13. Customer Service | 23 |
| 14. Frequently Asked Questions | 24 |
| 15. Specifications | 26 |
| 16. Labeling | 27 |
| EC Declaration of Conformity | 28 |

Before using Flash N' Go™ for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures.

We recommend you re-familiarize yourself with this User Manual before each use of Flash N' Go™.



Flash N' Go™ is a powerful electrical device. As such, it should be used with special attention to safety.



Please read all warnings and safety precautions before use, and strictly follow them when using Flash N' Go™

1. What is Flash N' Go™?

Flash N' Go™ is a light-based device for long-term hair removal designed for home-use.

1.1. Flash N' Go™ device description

Flash N' Go consists of a Base Unit, APPLICATOR and Disposable Lamp Cartridges.



1. BASE UNIT
2. APPLICATOR
3. PULSE BUTTON
4. DISPOSABLE LAMP CARTRIDGE
5. SKIN COLOR SENSOR
6. CONTROL PANEL
7. POWER ON/OFF SWITCH
8. POWER INDICATOR LIGHT
9. ENERGY LEVEL SETTING BUTTONS
10. ENERGY LEVEL INDICATOR LIGHTS
11. READY INDICATOR LIGHT
12. AUDIO ON/OFF SWITCH
13. CARTRIDGE 90% USAGE INDICATOR LIGHT
14. CARTRIDGE 100% USAGE INDICATOR LIGHT
15. SYSTEM WARNING INDICATOR LIGHT

Your Flash N' Go™ consists of a **BASE UNIT** and an **APPLICATOR**.

On the **BASE UNIT** you can find the Flash N' Go™ **CONTROL PANEL** including **POWER ON/OFF SWITCH** and **POWER INDICATOR LIGHT**, the **ENERGY LEVEL SETTING BUTTONS** and **ENERGY LEVEL INDICATOR LIGHTS**, the **AUDIO ON/OFF SWITCH**, **READY INDICATOR**, **SYSTEM WARNING INDICATOR LIGHT** and **CARTRIDGE 90% and 100% USAGE INDICATOR LIGHTS**.

On the Flash N' Go™ **APPLICATOR** you can find the **PULSE BUTTON**. The **DISPOSABLE LAMP CARTRIDGE** is located at the **APPLICATOR TIP**.

1.2. Package contents

Upon opening the Flash N' Go™ package, you will find the following parts:

- Flash N' Go™ **BASE UNIT** and **APPLICATOR**
- An AC cord
- A second **DISPOSABLE LAMP CARTRIDGE**
- This User's Manual and a Quick Start Guide leaflet
- An Instructional DVD

1.3. The Disposable Lamp Cartridge

The Flash N' Go™ **DISPOSABLE LAMP CARTRIDGE** can fire 750 light pulses (regardless of the energy level of these pulses). For average body size 750 light pulses could cover for example:

One treatment sessions of 2 Legs, 2 Arms, 2 Underarms, and Bikini line

Or

Two treatments sessions of 2 legs and 2 Underarms

Or

Ten treatments sessions of 2 Underarms and Bikini line

Depending on the body areas you wish to treat, you will typically need to use from one to seven **DISPOSABLE LAMP CARTRIDGES** during the first year of using Flash N' Go.

2. Flash N' Go™ Intended Use

Flash N' Go™ is intended for removal of unwanted hair. Flash N' Go™ is also intended to effect long term, or permanent hair reduction. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs re-growing after a treatment regime.

Flash N' Go™ may be used to remove unwanted body hair. Ideal body areas for Flash N' Go™ use includes the underarms, bikini line, arms and legs.

3. Safety with Flash N' Go™

3.1. With Flash N' Go™ Safety Comes First

- **HPL™ technology in Flash N' Go™ - Superior safety with lower energy level**

Home Pulsed Light™ technology is able to achieve long-term hair removal results at a fraction of the energy level used in other light-based hair removal equipment. The low energy used in Flash N' Go™ reduces its potential to cause harm or complications, and contributes to your overall safety.

- **Flash N' Go™ protects your eyes**

The Flash N' Go™ *APPLICATOR* has a **built-in safety feature** for eye protection. It has been designed so that a light pulse can not be emitted when the *APPLICATOR* is facing open air. The safety switch is activated only when the *APPLICATOR TIP* is in full contact with the tissue and pressed.

- **Flash N' Go™ protects your skin**

Flash N' Go™ comes with a *SKIN COLOR SENSOR* enabling use only on lighter skin complexions to ensure skin safety.

Furthermore, the 2X3cm² *OPTICAL LENS* through which pulses of light are delivered is recessed inside the *LAMP CARTRIDGE* at the *APPLICATOR TIP*. This enables Flash N' Go™ to protect your skin by avoiding direct contact between the *OPTICAL LENS* and the skin.

4. Contraindications

Important Safety Information – Read Before Use!

Flash N' Go™ is not designed for everyone. Please read and consider the information in the following section before use. For further information and personalized advice you may also visit www.silkn.com or your local Flash N' Go™ domain.


If you are unsure if Flash N' Go™ is safe for you to use, please ask your doctor or dermatologist!

- **DO NOT USE Flash N' Go™ on naturally dark skin complexion!**

Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™.

Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

Not Safe



DO NOT USE Flash N' Go™ on naturally dark skin!

A unique *SKIN COLOR SENSOR* is embedded in Flash N' Go™ to measure the treated skin complexion at the beginning of each session and occasionally during the session. *SKIN COLOR SENSOR* ensures that pulse will be emitted on suitable skin tones.

- **DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure!**

Tanned skin particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk when using Flash N' Go™.

Using Flash N' Go™ on skin that has been recently exposed to sunlight can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).



DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure! Such use can cause serious burns or skin injury. Avoid exposure to the sun for 4 weeks before your Flash N' Go™ treatment!

The unique Flash N' Go™ *SKIN COLOR SENSOR* will also help you avoid treating sun exposed skin.

- **DO NOT USE Flash N' Go™ on the face or neck**

Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair.

Using Flash N' Go™ to remove facial hair may cause serious eye injury and may stimulate facial hair growth.



Flash N' Go™ is not recommended for use on the face or neck.

- **NOTE!** Flash N' Go™ is not effective on light hair

The Flash N' Go™ device is not effective on naturally white, grey, blond and red body hair. If your body hair is of these colors, Flash N' Go™ will not work on you.

5. When to avoid using Flash N' Go™?

Certain conditions may prevent the use of Flash N' Go™ temporarily. **DO NOT USE Flash N' Go™** if any of the following currently apply to you:

- If you are pregnant or nursing (lactating).
- If you were exposed to strong sunlight or an artificial tanning machine during the past 28 days.
- If you have a tattoo or permanent makeup on the area to be treated.
- If you have dark brown or black spots, such as large freckles, birth marks, moles or warts on the area to be treated.
- If you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sore in the area to be treated. Wait for the effected area to heal before using Flash N' Go™.
- If you have a history of keloidal scar formation, a known sensitivity to light (photosensitivity) or are taking medication that makes the skin more sensitive to light, including non-steroidal anti-inflammatory agents, (e.g., aspirins, ibuprofens, acetaminophen), tetracyclines, phenothiazines, thiazide, diuretics, sulfonluraes, sulfonamides, DTIC, fluorouracil, vinblastine, griseofulvin, Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, Accutane® and/or topical retinoids.
- If you have abnormal skin conditions caused by diabetes, for example, or other systemic or metabolic diseases
- If you are currently or have recently been treated with Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, topical retinoids or azelaic acid.

- If you have been treated with Accutane® (isotretinoin) within the past 6 months.
- If you have been on a steroid regimen within the past 3 months.
- If you have a history of herpes outbreaks in the area of treatment, unless you have consulted your physician and received preventative treatment before using Flash N' Go™.
- If you suffer from epilepsy.
- If you have an active implant, such as a pacemaker, incontinence device, insulin pump, etc.
- If you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.
- If you have a history of skin cancer or areas of potential skin malignancies.
- If you have received radiation therapy or chemotherapy treatments within the past 3 months.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.

If you are unsure if Flash N' Go™ is safe for you to use, please consult with your doctor or dermatologist.

6. Precautions – How to use Flash N' Go™ Safely

- **Choose your energy levels CAREFULLY!**

Energy level refers to the intensity of the light pulse that is projected on your skin during use, from the lowest level (-) to the highest level (+). *INDICATOR LIGHTS* on the *CONTROL PANEL* illustrate the energy level at which the machine is set. As the energy level increases, so do the results of Flash N' Go as well as the risk of side effects (see "Possible Side Effects" below).

Always begin your first use of Flash N' Go™ at the lowest energy setting (one light at "-")!

Only if you experience little or no discomfort during and after use of Flash N' Go™ at the lowest energy setting, raise the energy level by one *INDICATOR LIGHT* the next time you use Flash N' Go™, and so on for each subsequent hair removal session.

For detailed instructions on energy level setting see "Energy Level" box in section 4.4 - "Treating with Flash N' Go™ for the first time".

- **Avoid adverse effects!**

Do not treat the same area of skin more than once per hair removal session!
Try to avoid overlapping pulses!
If your skin blisters or burns, **STOP USE IMMEDIATELY!**

- **Avoid complications after use of Flash N' Go™!**

Do not expose treated areas of skin to the sun!

Sun exposure includes unprotected exposure to direct sunlight of over 15 minutes constantly, or unprotected exposure to diffuse sunlight of over 1 hour constantly.

To protect recently treated skin when exposed to sunlight, be sure to thoroughly apply sunscreen SPF 30 or higher, for 2 weeks after each hair removal session.

- Always shave the area to be treated and make sure that the skin is clean and dry before using the Flash N' Go™.
- Cover birthmarks and tattoos before Flash N' Go™ application.
- Cover dark brown or black spots, such as large freckles, birth marks, moles or warts before Flash N' Go™ application.
- Never look directly at the light coming from the Flash N' Go™ *APPLICATOR* and *LAMP CARTRIDGE*.
- Do not use Flash N' Go™ on nipples and genitals (male or female).
- Do not use Flash N' Go™ on any body site where you might later want hair.
- Do not use Flash N' Go™ for any purpose other than hair removal.
- Never point the Flash N' Go™ *APPLICATOR* in an attempt to emit a light pulse into open space. Always make sure that the *APPLICATOR* is pointed at, and in full contact with the skin during application.
- Remove the Flash N' Go™ *APPLICATOR* from the skin if either the skin or the *APPLICATOR* is too hot.
- Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone to clean the skin before using Flash N' Go™.
- Use of Flash N' Go™ may cause temporary pigmentation changes (See "Possible Side Effects" below).
- Keep this device out of the reach of children. Do not use Flash N' Go™ on children or allow children to use it.

7. Reducing the risk of injury

As with any electrical device, certain precautions must be taken in order to ensure your safety when using Flash N' Go™.

- **Keep Flash N' Go™ away from water!**
Flash N' Go™ is an electrical device. As such it should always be kept away from water.
Do not place or store Flash N' Go™ where it can fall or be pushed into a tub, sink or any other vessel containing water. Do not place in, or drop into water or any other liquid.
This may cause severe electrocution.
Do not use Flash N' Go™ while bathing.
Do not use Flash N' Go™ if it becomes damp or wet.
Do not reach for Flash N' Go™ if it has fallen into water.
Unplug Flash N' Go™ immediately if it has fallen into water.



Keep Flash N' Go™ away from water!

- **Never open Flash N' Go™!**
Do not attempt to open or repair your Flash N' Go™ device. Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

- Trying to open Flash N' Go™ may also damage the device and will void your warranty.
Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.
- Use Flash N' Go™ only for its intended use and as described in its manual.
- Flash N' Go™ should never be left unattended when plugged into an outlet.
- Do not operate Flash N' Go™ if it has a damaged cord or plug and keep the power cord away from heated surfaces.
- Do not use Flash N' Go™ if you see or smell smoke when it is in use.
- Do not use Flash N' Go™ if it is not working properly or if it appears damaged.
- Do not use Flash N' Go™ if the fan vent in its *APPLICATOR* is cracked, coming off or missing.
- Do not use Flash N' Go™ if the *SKIN COLOR SENSOR* in its *APPLICATOR* is cracked, or broken.
- Do not use Flash N' Go™ if the outer shell is cracked or is coming apart.
- Do not use Flash N' Go™ with a damaged *DISPOSABLE LAMP CARTRIDGE*, or if its *OPTICAL LENS* is cracked, chipped or missing
- Always unplug Flash N' Go™ from the electrical outlet immediately after use.
- Unplug Flash N' Go™ before cleaning.
- Do not use Flash N' Go™ with any attachments or accessories not recommended by Home Skinovations Ltd.

8. Risks of using Flash N' Go™

When used according to the instructions, side effects and complications associated with use of Flash N' Go™ are uncommon. However every cosmetic procedure,

including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that can occur with pulsed light hair removal systems designed for home use.

| Adverse Event | Likelihood of the Adverse Event | Adverse Effect | Likelihood of the Adverse Effect |
|--|---------------------------------|-----------------------------------|----------------------------------|
| Stacking or Overlapping of multiple pulses on the same skin spot | Minor | Minor Skin Discomfort | Minor |
| | | Skin Redness | Minor |
| | | Increased Sensitivity of the Skin | Minor |
| | | Skin Wounds and Burns | Rare |
| | | Scarring | Rare |
| | | Pigment Changes | Rare |
| | | Excessive Redness and Swelling | Rare |
| | | Infection | Negligible |
| | | Bruising | Negligible |

- **Minor Skin Discomfort**

Although home pulsed light hair removal is generally very well-tolerated, most users do feel some mild discomfort during use, usually described as being a mild stinging sensation on the treated skin areas. The stinging sensation usually lasts during the time of the application itself or for a few minutes thereafter. Anything beyond this minor discomfort is abnormal and means that either you should not continue to use Flash N' Go™ because you are unable to tolerate the hair removal application, or that the energy level setting is too high.

- **Skin Redness**

Your skin may become red right after using Flash N' Go™ or within 24 hours of using Flash N' Go™. Redness generally clears up within 24 hours. See your doctor if redness does not go away within 2 to 3 days.

- **Increased Sensitivity of the Skin**

The skin of the treated area is more sensitive so you may encounter dryness or flaking of the skin.

- **Skin Wounds and Burns**

Very rarely, burns or wounds to the skin can occur following the application. The burn or wound can require a few weeks to heal and, extremely rarely, may leave a noticeable permanent scar.

- **Scarring**

Although very rare, permanent scarring may occur. Usually when scarring occurs it is in the form of a flat and white lesion on the skin (hypotrophic). However, it can be large and red (hypertrophic) or large and extend beyond the margins of the injury itself (keloid). Subsequent aesthetic treatments may be required to improve the appearance of the scar.

- **Pigment Changes**

Flash N' Go™ targets the hair shaft, in particular the pigmented cells in the hair follicle and the hair follicle itself. Nevertheless there is risk of temporary hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. This risk of changes in skin pigmentation is higher for people with darker skin tones.

Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.

- **Excessive Redness and Swelling**

In rare cases treated skin may become very red and swollen. This is more common in sensitive areas of the body. The redness and swelling should subside within 2 to 7 days and should be treated with frequent applications of ice. Gentle cleansing is OK, but one should avoid exposure to sun.

- **Infection**

Infection of the skin is exceedingly rare but is still a possible risk following a skin burn or wound caused by Flash N' Go™.

- **Bruising**


Very rarely, use of Flash N' Go™ may cause blue-purple bruising that can last 5 to 10 days. As the bruising fades, there may be a rust-brown discoloration of the skin (hyperpigmentation) that can be permanent.

9. Long Term Hair Removal the Flash N' Go™ Way

Flash N' Go™ is a personal light-based device for long-term hair removal. The process of laser and light-based hair removal is well known and established. It has been proven in clinical use around the world for over 15 years as a safe and effective way to achieve long-term hair reduction.

9.1. How does light remove hair?

Light-based hair removal is based on the theory of selective photothermolysis in which optical energy is used to disable hair growth. In order to achieve such thermal effect the hair shaft needs to selectively absorb light energy and transform it into heat. This selectivity is achieved when optical energy that is delivered to the tissue is mostly absorbed by hair shaft pigment, while the skin and the surrounding tissue stays cool. Melanin is the pigment in the hair shaft that is responsible for the absorption of the light, which generates the heat that eventually disables hair growth. Therefore the more melanin present in the hair (i.e. the darker the hair) the more light that can be absorbed and the more effective light can be at removing hair.




Before Flash N' Go™ Application
Go™ Application

During Flash N' Go™ Pulse After Flash N'

9.2. How does the hair growth cycle impact light-based hair removal?


Every hair in our body goes through the three phases of the hair growth cycle: Anagen, Catagen and Telogen. These phases have an important impact on how the process of light-based hair removal works.

Anagen is the hair growth phase while Catagen and Telogen, both are resting phases.



Anagen—Growth phase Catagen—Resting phase Telogen—Resting phase

The time it takes to complete a full hair growth cycle varies from person to person and the location of the hair on the body, but is typically 18-24 months. At any given moment the majority of the hair follicles in any skin area are in the resting phases. These resting hairs cannot be affected by Flash N' Go™. However, hairs in the growing Anagen phase will respond to Flash N' Go™ applications. It is important to understand that it may take a full hair growth cycle to realize complete hair removal results with Flash N' Go™.



9.3. Plan your Flash N' Go™ hair removal for best results

A typical full hair growth cycle may take 18-24 months. During this time multiple Flash N' Go™ sessions may be required in order to achieve long term hair removal.

The efficiency of hair removal varies from person to person according to body area, hair color, and how Flash N' Go™ is used.

Typical Flash N' Go™ hair removal plan during a full hair growth cycle:

- The first 3-4 hair removal sessions with Flash N' Go™ will be approximately two weeks apart.
- Hair removal sessions 5-7 with Flash N' Go™ will be approximately four weeks apart.
- After that you will typically use Flash N' Go™ again from time to time if and when needed, until long-term results are achieved.

This is the recommended treatment schedule that has proven to produce the best results, but you may plan your personal treatment schedule differently and still get satisfactory results.

Note that treating the same body part more often than once in two weeks, will not improve hair removal results.

Typical maintenance with Flash N' Go™ after achieving long-term hair removal:

Due to hormonal or other physiological changes dormant hair follicles may become active. Maintenance hair removal sessions with Flash N' Go™ may be required from time to time.

10. How To Use Flash N' Go™

10.1. Flash N' Go™ Setup

1. Remove Flash N' Go™ *BASE UNIT*, *APPLICATOR* and other components from box.
 2. Insert the *APPLICATOR* into its cradle in the Flash N' Go™ *BASE UNIT*.
 3. Verify that the *DISPOSABLE LAMP CARTRIDGE* is inserted correctly into the *APPLICATOR*.
 4. Plug the power cord into the Flash N' Go™ *BASE UNIT* socket.
 5. Plug the other end of the power cord into an electrical outlet.
- Your Flash N' Go™ is now ready to start.

10.2. Treating with Flash N' Go™ for the first time

The skin should be **shaved, clean, dry** and **free** of any powders, antiperspirants or deodorants.

1. Press the *POWER ON/OFF SWITCH*. The *POWER INDICATOR LIGHT* will turn on and a fan sound (similar to the sound of a hairdryer) will start.
2. Approximately 3.5 seconds after pressing the *POWER ON/OFF SWITCH*, the *READY INDICATOR LIGHT* will turn on and the device beeps. The device sets itself to deliver the lowest energy level pulses. The device is then ready for you to trigger the first pulse.
3. If the energy level should be higher than the lowest, press the *ENERGY LEVEL SETTING BUTTONS* using the "-" or "+" to respectively decrease or increase the energy level, until the desired energy level is set. *ENERGY LEVEL INDICATOR LIGHTS* will coincide with the energy level setting. (See "Energy Level" frame).
4. Apply the *APPLICATOR* firmly to the skin, making sure the skin is spread evenly and smoothly. Make sure the *APPLICATOR TIP* is in full contact with the skin.
5. Press the *PULSE BUTTON*. The device will first determine the color of your skin.
If the color of the skin is light enough for safe application, the device will flash a pulse of light onto your skin. You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You may feel a mild sensation of warmth and tingling.
Flash N' Go™ will immediately recharge for the next pulse. After 3.5 seconds the *READY INDICATOR LIGHT* will turn on again and the device will beep.
If you see no light pulse, and the 5 green energy level LEDs blink for 3 seconds, your skin tone is too dark for safe application. Try using the device on a different body part or contact the Flash N' Go support.
6. Remove the *APPLICATOR TIP* from the treated area of skin.
7. Move the *APPLICATOR TIP* to another area of skin.
Use the pressure marks the device just made in your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses (See "Covering skin areas" frame).
Warning: Do not treat the same area of skin more than once per hair removal

session! Treating the same area of skin more than once per session increases the likelihood of adverse effects.

8. Repeat the process starting again with Step 5.

Energy Level

Energy level determines the intensity of the Flash N' Go™ light pulse delivered to your skin, from the lowest level (-) to the highest level (+). Corresponding *INDICATOR LIGHTS* on the *CONTROL PANEL* represent the increases in energy. As energy level increases, so does hair removal results as well as the risk of possible side effects and complications.

Always start your first Flash N' Go™ hair removal session at the lowest energy setting.

If you experience little or no discomfort during and after the hair removal session using the lowest energy level, raise the energy level by one *INDICATOR LIGHT* at the next session, and so on for each subsequent session.



Do not raise the energy level if you experience abnormal discomfort during or after treating with Flash N' Go™ (See section 2.6 - "Possible Side Effects"). Do not raise the energy level during hair removal session even if you experience no discomfort. Discomfort may also appear some time after the session.

Whenever Flash N' Go™ is turned on its energy level will automatically be reset to the lowest energy level. Only one *ENERGY LEVEL INDICATOR LIGHT* will be on.


To set the energy level, press the *ENERGY LEVEL SETTING BUTTONS* using "-" or "+" to respectively decrease or increase the energy level. The number of *ENERGY LEVEL INDICATOR LIGHTS* will coincide with the change in energy level.

For your safety, when using Flash N' Go™ for the very first time, the system is automatically set to deliver the first 50 pulses at the lowest energy level, and the next 200 pulses at up to level 3.

To cancel these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 2 consecutive "beep" sounds. To **restore** these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 3 consecutive "beep" sounds.

Covering skin areas

Flash N' Go™ pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas



When applied to the skin, the Flash N' Go™ *APPLICATOR TIP* is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.



Try to avoid overlapping pulses!

Do not treat the same area of skin more than once per hair removal session!



If your skin blisters or burns, STOP USE IMMEDIATELY!

10.3. What to Expect when treating with Flash N' Go™?

For many people, using Flash N' Go™ may be their first experience with a light-based device designed for home use. Flash N' Go™ is simple to use, and hair removal sessions go by quickly. During a Flash N' Go™ session it is **normal** to experience and feel:

- **A Flash of Light** – The bright light of Flash N' Go™ will not harm the eyes when applied to non-facial sites, and special eye-protection is not needed when using Flash N' Go™.
- **A Fan Noise** – The cooling fan in Flash N' Go™ makes noise similar to a hairdryer. This is normal.
- **A Pop Sound with Each Pulse** – When a pulse of HPL™ light is activated, it is normal to hear a subtle pop sound simultaneously with the flash of light.
- **Moderate Pressure of the APPLICATOR** – This is necessary and helpful for placement of adjacent pulses of light, and is part of the unique safety feature of Flash N' Go™.
- **A Sensation of Warmth and Tingling** – During each pulse of light it is normal to feel a mild sensation of warmth and tingling from the light energy. Remember it is important to always use low energy settings for initial hair removal sessions. You may feel some warmth for up to an hour after your Flash N' Go™ session.
- **Some Mild Red or Pink Color** – During and just after your Flash N' Go™ session it is not uncommon to see some very mild, pink-like color of the skin. This is usually most noticeable around the hairs themselves. However if you see full redness of the skin, blistering or burns stop use of Flash N' Go™ immediately.

10.4. After treating with Flash N' Go™

- When Flash N' Go™ session has been completed turn Flash N' Go™ off by pressing the *POWER ON/OFF SWITCH*. (Be sure to remember the last energy level setting you used, as it will not be restored when turning Flash N' Go™ on again.)
- Unplug the power cord from the electrical outlet.
- After each hair removal session it is recommended that you clean your Flash N' Go™ device, especially the *APPLICATOR TIP* (See section 5.1: "Cleaning Flash N' Go™").
- After cleaning, it is recommended to store your Flash N' Go™ device in its original box, and keep it away from water.



Skin care following hair removal session

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last Flash N' Go™ session.

Side effects and complications

Some patients may experience pigmentation changes resulting from treating with Flash N' Go™. These effects, if they occur, are generally mild and transient.



In case you experience any complication (See section 2.6 Possible Side Effects using Flash N' Go™) ***please contact your physician immediately.***

11. Maintenance of Flash N' Go™

11.1. Cleaning Flash N' Go™

After each hair removal session, it is recommended to clean your Flash N' Go™ device, and especially the *APPLICATOR TIP*.

- Unplug Flash N' Go™ before cleaning.
- Use a dry, clean cloth and a specially formulated cleaner for electronic equipment to gently wipe Flash N' Go™ surface, and especially the *APPLICATOR TIP*.
Never immerse Flash N' Go™ or any of its parts in water!

11.2. Replacing the LAMP CARTRIDGE


The Flash N' Go™ *DISPOSABLE LAMP CARTRIDGE* can fire 750 light pulses that would typically cover 2 Legs, Arms, Under arms, and Bikini line. Pulse intensity is determined only according to the energy level setting of the device. There is no decrease of energy during the usable lifetime of the *LAMP CARTRIDGE*.

When a *DISPOSABLE LAMP CARTRIDGE* has reached 90% of its possible lifetime the *CARTRIDGE 90% USAGE INDICATOR LIGHT* on the *CONTROL PANEL* will be activated, indicating that the *DISPOSABLE LAMP CARTRIDGE* should soon be replaced.

Once all 750 pulses in a *DISPOSABLE LAMP CARTRIDGE* have been used the *CARTRIDGE 100% USAGE INDICATOR LIGHT* will be activated, indicating that further pulses can not be emitted.

In order to continue the hair removal session, it will now be necessary to replace the *DISPOSABLE LAMP CARTRIDGE*.

Note: A *DISPOSABLE LAMP CARTRIDGE* should also be replaced if large spots appear inside it or if its *LENS* is broken.



To replace a *DISPOSABLE LAMP CARTRIDGE*:

1. Press the *POWER ON/OFF SWITCH* to turn Flash N' Go™ off.
2. Unplug the power cord from the electrical outlet.
3. Grasp the used *DISPOSABLE LAMP CARTRIDGE* on both sides, pull it out of the

socket and discard with normal trash.

4. Unwrap a new *DISPOSABLE LAMP CARTRIDGE*.

5. Line up the holes in the *APPLICATOR* socket with the metal and plastic hinges in the new *DISPOSABLE LAMP CARTRIDGE*, and push it gently into place.

If correctly installed, the *DISPOSABLE LAMP CARTRIDGE* will spring back when pushed in and released.

12. Troubleshooting

12.1. "My Flash N' Go™ does not start."

- Make sure the power cord is properly connected to the Flash N' Go™ device.
- Make sure the power cord is plugged into an electrical outlet on the wall.

12.2. "A light pulse is not emitted when I press the PULSE SWITCH."

- Make sure that you have good contact with the skin and that the *APPLICATOR TIP* is evenly and firmly pressed to the skin. For your safety, the *PULSE BUTTON* will activate a pulse only if the *APPLICATOR TIP* is firmly pressed against the skin.
- Skin color sensor stopped the light pulse, your skin may be too dark for using the Flash N' Go. Try to apply Flash N' Go™ again.
- Check the *CARTRIDGE 100% USAGE INDICATOR*. If it is on, disconnect Flash N' Go™ and replace the *DISPOSABLE LAMP CARTRIDGE*.
- Make sure that the *READY INDICATOR* is on.
 - a. If within 10 seconds the *READY INDICATOR* remains off turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - b. If the problem persists, contact your local Flash N' Go™ Customer Service Center.
- Check the *SYSTEM WARNING INDICATOR LIGHT*.
 - c. If the light is on turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - d. If the problem persists, contact your local Flash N' Go™ Customer Service Center.

If these problems persist, contact your local Flash N' Go™ Customer Service Center.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.

Trying to open Flash N' Go™ may also damage the device and will void your warranty.

Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.

13. Customer Service

For more information about Flash N' Go™ please enter www.silkn.com.

If your Flash N' Go™ is broken, damaged; in need of repair, or for any other Flash N' Go™ user assistance, please contact Flash N' Go™ Customer Service:

1-877-DO-SILKN / 1-877-367-4556

contact@silkn.com

14. Frequently Asked Questions

1. ***Does Flash N' Go™ really work?***
Yes. In clinical trials held by physicians, Flash N' Go™ was proven to safely achieve long-term hair removal results.
2. ***Where on my body can I use Flash N' Go™?***
The Flash N' Go™ device has been designed for body hair removal anywhere below the neck. The most common areas treated with Flash N' Go™ are: legs, underarms, arms and bikini line. It is not recommended to use Flash N' Go™ on the face. Flash N' Go™ is not suitable for everyone. We recommend that you read all the Flash N' Go™ warnings and contraindications in this User's Manual.
3. ***How long does a Flash N' Go™ treatment session take?***
The time can vary depending on the area of the body treated. A full leg can take up to 30 minutes, or two underarms could take 10 minutes. Because Flash N' Go™ runs on regular electric power it can be used for as long as needed to complete a full hair removal session of the desired body part(s).
4. ***Is Flash N' Go™ safe?***
Flash N' Go™ has been designed with your safety in mind, and tested and approved by top dermatologists and plastic surgeons to meet their safety standards for a home-use device. But like any skin product or electronic device, one must use according to the operating instructions and user precautions.
5. ***Will Flash N' Go™ hurt?***
When used properly most users of Flash N' Go™ report feeling a slight sensation of heat at the time of the pulse of light. Users with thicker and darker hairs may feel slightly more discomfort, but this discomfort subsides once the hair removal session is completed. For your convenience Flash N' Go™ has five setting levels that can be used according to your sensitivity.
6. ***How often should I use Flash N' Go™?***
Hair removal sessions with Flash N' Go™ should be spaced every two weeks for the first three to four sessions. After that hair removal sessions should be done if hairs have grown back, until the desired results are achieved.
7. ***Is Flash N' Go™ effective on white, grey or blonde hairs?***
Flash N' Go™ works best on darker hair types, or hair that contains more melanin. Melanin is the pigment that gives hair and skin their color, and will absorb light energy. Black and dark brown respond the best, although brown and light brown hairs will also respond but typically require more hair removal sessions. Red may show some response. White, grey or blonde hairs usually don't respond to Flash N' Go™ though some users have noted results after multiple hair removal sessions.
8. ***Can I use Flash N' Go™ on brown or black skin***
Do not use Flash N' Go™ on naturally dark skin complexion! Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™. Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).
9. ***When will I see results from Flash N' Go™?***
As with any light-based or laser hair removal device, results are not immediate, and in fact you may not think anything happened at all. Hair may sometimes appear to be growing back after a hair removal session, but typically after two weeks many of these hairs will simply fall out. Additionally hairs grow in three different stages and only hairs in an active growth stage will be affected by Flash N' Go™. This is one of the main reasons that multiple sessions are required to achieve the desired result.

10. **Can a man use Flash N' Go™?**

Though designed exclusively for women, Flash N' Go™ may be suitable for use by men. However, hairs on men, typically those on the chest, will require more hair removal sessions than that of women to get the desired results. As with women, using Flash N' Go™ on men's facial hair such as beard and mustache, is not recommended.

11. **Why is my hair growing, even though I treated it a week ago?**

It is quite common for hair to appear as if it is still growing up to two weeks after a hair removal session with Flash N' Go™. This process is known as "ejection" and at around two-weeks you'll see that these hairs simply fall out or slide out with a slight tug. (We don't however recommend pulling on the hairs – just let them come out naturally.) It is also possible that some hairs, due to missed application or different stages of growth, were not affected by Flash N' Go™. These hairs will be treated in follow-up sessions, and hence the reason multiple hair removal sessions are needed to get the best result with Flash N' Go™.

12. **I've heard that some hairs grow back lighter and finer after light based hair removal?**

This phenomenon is well documented amongst aestheticians and doctors using light and laser devices for hair removal. It is possible that some hairs will grow back lighter and finer after hair removal with Flash N' Go™. Usually these hairs are a fraction of what was originally there, and continued use may have a desirable effect on them.

13. **Why can't I treat myself if I have an "active" suntan?**

Do not use Flash N' Go™ on tanned skin or after sun exposure! Tanned skin and particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk of adverse effects when using Flash N' Go™ including burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

14. **Is long-term use of Flash N' Go™ dangerous for my skin?**

The use of light and laser energy in aesthetic medicine has been well documented for over 15 years in professional peer-reviewed journals, and by well respected institutions like the Mayo Clinic. These journals and institutions have not reported any side-effects or damage from long-term use of light and laser device.

15. **Can I use Flash N' Go™ to remove my chin hair or elsewhere on my face?**

It is not recommended to use Flash N' Go™ on the face or neck. Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair. In addition, using Flash N' Go™ to remove facial hair may cause serious eye injury.

16. **How long should I wait to treat with Flash N' Go™ after unprotected exposure to the sun?**

One should wait 4 weeks before using Flash N' Go™ after unprotected exposure to the sun. However, if there is ever any uncertainty about sun exposure please contact your physician or Home Skinovations customer support.

17. **Should I do anything before using Flash N' Go™?**

Before any Flash N' Go™ session it is important to avoid sun exposure on the treated area for at least four weeks. A high level UV Sun Screen (SPF 50+) will help, as will clothing covering the treated area. The area to be treated should also be cleaned with mild soap and water, and the hairs shaved down to skin level.

18. **How should I care for the treatment area after using Flash N' Go™?**

The area treated with Flash N' Go™ can be cleaned and maintained with standard skin care products. Special care must be taken to avoid unprotected sun exposure. Strong sunscreens (50+ SPF) and covering clothing are suitable for protection from the sun.


19. **Should I pull the hairs out after treatment?**

No, let the hairs gradually fall out on their own. This may take up to 2 weeks.







15. Specifications

| | |
|---------------------------------|--|
| Spot size | 2cm x 3cm [6cm ²] |
| Speed | 1 pulse every 3.5 second: 1.7 cm ² /sec |
| Technology | New Home Pulsed Light™ |
| Max Energy Level | Max 5J/cm ² |
| Wavelength | 475-1200nm |
| Charging time / Power source | Continuous operation |
| Electrical requirements | 100-240VAC, 2A |
| Time needed to treat lower legs | 30 minutes |
| Operation and Safety | SKIN COLOR SENSOR seamlessly ensures use only on appropriate skin types. |
| package size | Height 9 inch, Width 9 inch, Depth 5.2 inch |
| System weight | 4 pounds |
| | |


16. Labeling



| | |
|------------|------------------------|
| MODEL : | FlashN'Go |
| PART NO. : | AS100756A |
| POWER: | ~100Vac ;2A 50-60Hz |



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation



000001001

008024

Made in ISRAEL

LB100756A

Class II equipment

Degree of protection against electric shock: type BF applied part

WEEE - Waste Electrical and Electronic Equipment

CE Mark

Follow operating instructions

This device comply with part 15 of the (FCC) Federal Communications Commission.

CSA Mark for USA and Canada

Degree of protection against ingress water: ordinary

This device is not suitable for use in the presence of flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.





Brand Name: **Home Skinovations**

Model Name: **Flash N' Go, Disposable Lamp Cartridge**

Models Number **AS100001E,AS100674B,AS100756A**

This declaration is based on the full compliance of the products with the following European standards:

| | |
|--|---|
| EMC: | EN 61000-3-2:2000; (A2:2005), EN 61000-3-3:1995; (A1:2001) EN 55014-1:2000; (A2:2002), EN55014-2:1997; (A1:2001) |
| Electrical safety: | EN 60335-1:2002; (A1:2004; A12:2006), EN 60335-2-23:2003, EN 60335-2-27:2003 |
| Human Exposure to electromagnetic field radiated | EN 50364:2001 |

| | |
|--------------------------------|-----------------------------------|
| Manufacturer: | Authorized Representative: |
| Home Skinovations Ltd. | Home Skinovations GMBH |
| Apolo Building, Shaar Yokneam, | Dr. Kurt Huber Str. 6 |
| POB 533, Yokneam 20692, ISRAEL | D-82031 Grünwald, GERMANY |
| Tel: +972-4-9097470 | Tel. +49-89-64919530 |
| Fax: +972-4-9097471 | Fax +49-89-64919531 |

Signed: _____ Signed: _____

Gabi Lavi
General Manager
Yokneam, February 4, 2008

Dr. Amir Waldman
VP Regulatory & Clinical Affairs
Yokneam, February 4, 2008

For Information: On the basis of this declaration, these products and packaging will bear the following mark:



87

Flash N' Go™ is a **physician recommended** and innovative light-based device for hair removal in the **privacy of your home**. **Safe, easy to use** and **cost effective**, Flash N' Go™ achieves **marked clinical results** while offering personal convenience.

Flash N' Go™ utilizes the highly sophisticated and new **Home Pulsed Light™ (HPL™)** technology developed by **Home Skinovations Ltd.**

Chapter 20

Performance testing – Clinical

Clinical Information

Introduction

Clinical trials were carried out in order to verify the safety and effectiveness of the Silk'n Flash N Go.

A prospective, multi-center study was carried out on patients with different skin color, age. (b) (4), enrolled and completed protocol. Investigator determined patient eligibility to participate in the study according to inclusion / exclusion criteria.

After completing the informed consent form, the investigator performed a base line hair count. Patients treated according to the treatment protocol. (b) (4)

(b) (4)

The patients had various skin pigmentations which varied from skin type (b) (4) different hair color. Effectiveness (b) (4) post treatment results to those of the pre-treatment. Percentage of hair reduction was calculated, and the resulting average hair reduction for all patients was (b) (4) months, (b) (4). Safety was evaluated by monitoring immediate reaction and adverse effects. (b) (4)

(b) (4)

were no other adverse effects associated with the treatment.

Clinical study

Method

The information, which is presented in this report, was obtained from a multi center clinical study that worked under IRB approved protocol specified in the relevant section below. The study was carried out at 3 centers in the U.S., Canada, and Israel.

Study Sponsor:
Home Skinovations Ltd.

Study monitor:
Dr. Amir Waldman

Apolo Buld. P.O.B. 533
Yokneam Illit 20692 Israel

Home Skinovations Ltd.
Apolo Buld. P.O.B. 533
Yokneam Illit 20692 Israel

The centers that conducted the study:

(b) (4)



Treatment protocol

The medical device that was used for this study was Silk'n Flash N Go, used for the clinical study under IRB approval. The detailed protocol, informed consent form, Post treatment instruction can be found in Appendix 1. Highlights of the study protocol are presented in this section. The device was used with the following parameters: Light optical output power (b) (4)

This is a prospective clinical study of (b) (4) designed to determine the safety and efficacy of the Silk'n Flash N Go device is an over the counter device intended for the removal of unwanted hair. Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

Purpose and objective

Purpose: The purpose of this clinical investigation is to study and document the use of the Silk'n Flash N Go, as an over the counter device, for permanent hair reduction.

Objectives: The objectives of the study are to determine the safety and effectiveness of low energy light pulses for permanent hair reduction.

Patient selection

Each patient was evaluated by the investigator to ensure that they met the inclusion / exclusion criteria of the study.

Inclusion Criteria

- 1- Presence of unwanted hairs on the body, (b) (4)
- 2- Skin Type (b) (4) (Fitzpatrick)
- 3- Females older than 21 years of age, must be either post-menopausal or surgically sterilized, or using a medically acceptable form of birth control (i.e., oral contraceptives, IUD, contraceptive implant, barrier methods with spermicide or abstinence).
- 4- Informed consent agreement by the subject.
- 5- Willingness to follow the treatment schedule and post treatment care

Exclusion Criteria

1. Malignant or pre-malignant pigmented lesions in the area to be treated
2. Scarring or infection of the area to be treated
3. Known photosensitivity
4. Pregnancy
5. Subjects with Diabetes (Type I or II)
6. Presence of a suntan in the area to be treated
7. Use of medication known to induce photosensitivity

8. Known anticoagulative or thromboembolic condition
9. Subjects with a pacemaker or internal defibrillator
10. Use of NSAIDS two weeks prior to, and two weeks following the treatment
11. Subjects that used waxing or other methods of photo-epilation 3 months prior to treatment

Clinical Procedure

Pre-Procedure Evaluation

- 1- The Investigator discussed with the subject other available treatment options including the Silk'n Flash N Go hair removal system.
- 2- During the first visit, the investigator obtained an informed consent from the subject, clearly indicating her understanding of the requirements and risks involved with study participation, including but not limited to:

(b) (4)

- a.
- b.
- c.
- d.
- e.
- f.
- g.
- h.
- i.

- 3- A full detailed subject chart recorded initial hair count in the area to be treated. Initial hair count was performed to assess treatment efficacy as the treatment progressed.

- 4- (b) (4)

If the subject met with all inclusion and exclusion criteria, the investigator obtained signed Informed Consent, enrolled the subject into the study, and recorded required data on the appropriate study data form. The pre-treatment data was stored in the computer under the subject study number ID.

A copy of the Informed Consent Form (ICF) and a copy of a Case Report Form (CRF) are presented as part of appendix 1.

Treatment procedure

The treatment performed in the investigator clinic by a trained Nurse or technician.

Treatment Parameters

Instructions:

All patients received treatments according with the user manual. (b) (4)

(b) (4)

Treatment test spot

To perform the Treatment test spot, follow sections 1-6 and wait 15 minutes to check if any side effects are noted:

1. In the area to be treated, hairs should be trimmed to 3/32 inch length, (1-2 mm), or shaved 3 days prior to the treatment.
2. Clean the skin from the trimmed hairs.
3. Hold the applicator and ensure that the light output window is clean.
4. Switch on the system and set the energy level.
5. Place the Applicator on the treatment area and apply slight pressure.
Press the trigger switch to emit light pulse.
6. Move the applicator to the next spot to ensure full coverage of the area to be treated.
7. Repeat the test spot procedure for 3 pulses.
8. If side effects are noted after 15 minutes, reduce the energy level by one level to lower setting.

9. Transient Erythema and follicular redness are normal reactions as well as smell of burned hairs.

10. If no side effects are noted after 15 minutes, complete the treatment to cover entire area to be treated.

Post Treatment Care.

Apply moisturizer after the treatment.

Use high factor sunscreen and protect the treated area from sunlight for at least 2 days after the treatment. Tanning after a treatment may cause hyper-pigmentation.

Repetitive treatments

Treatments will be performed (b) (4)

Photographs using the digital camera system should be taken before treatment and at

(b) (4)

Complications

If an unanticipated adverse event occurs at any time during or after the use of the system, the Investigator must report it to the sponsor. A serious adverse event or complication directly attributable to the use of the System will be reported to the Sponsor as soon as possible within twenty-four (24) hours or sooner.

Data Analysis

Hair loss at each follow up session was calculated as the ratio of the number of hairs after the treatment compared with the hair count prior to treatment. Hair count will be preformed by the investigator.

Results


A total of (b) (4) volunteered to participate in the study. All of them were treated under the study protocol. Treated area for all patients were (b) (4)

(b) (4) A fellow physician, that didn't enroll the patients, performed evaluation of the treatment results by hair count at the (b) (4) follow up session. Hair

clearance was calculated as the percentage of (b) (4) follow up count compared with the initial hair count.

Patient information summary

The population that was studied covered a range of ages, and skin pigmentations. These distributions are shown in the graphs in Figures 1, 2 and 3. Raw data tables showing patient information can be found in table 1.



(b) (4)

Treatment results summary

Hair reduction was calculated by dividing the (b) (4) follow up count by the initial baseline hair count. The overall hair clearance for all patients, and all sites, (b) (4)

Figure 2 summarizes distribution of participants by hair-reduction/clearance category.

(b) (4)


months and is in average (b) (4) while the majority of the participant experienced (b) (4) -

(b) (4) The (b) (4)

groups, which proves that the hair count is stable and therefore (b) (4)

of the hairs in the treated area achieved. (b) (4)

(b) (4)



(b) (4)

Complications

The immediate skin reaction (b) (4)

(b) (4) other side effects were noted.

Conclusion

The results presented above lead to the conclusion that the Silk'n Flash N Go is safe and effective as an over the counter device intended for the removal of unwanted hair, and for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

**Clinical Study to Determine the Safety and Efficacy of the Silk'n
Home Pulsed Light Device**

Clinical Investigation Plan

Date: (b) (4)

Protocol number: (b) (4)

Version (b) (4)

Sponsor:

Home Skinovations Ltd.
P.O.B. 533 Yokneam Illit
20692 Israel
Tel: 972-4-9097471
Fax: 972-4-9097470

Table of Contents

| | |
|---|----|
| Table of Contents..... | 2 |
| Protocol Summary | 3 |
| Name and Intended Use of Device | 4 |
| Purpose and Objectives | 4 |
| Scope and Duration of the Study | 4 |
| Subject Selection | 5 |
| Clinical Procedure..... | 6 |
| Data Analysis | 9 |
| Materials | 10 |
| Training Requirements..... | 10 |
| Risk/ Benefit Analysis | 10 |
| Risks | 10 |
| Benefits | 10 |
| Alternatives | 10 |
| Records and Reports | 11 |
| Confidentiality | 11 |
| Monitoring Procedures..... | 11 |
| <u>Appendix 1 Study schematics</u> | |
| <u>Appendix 2 Inform consent form</u> | |
| <u>Appendix 3 Case follow up form</u> | |
| <u>Appendix 4 Post treatment instructions</u> | |

Protocol Summary

This is a clinical study to determine the efficacy and safety of low energy intense pulsed light, (IPL) device for hair removal. Treatment sites will include (b) (4) (b) (4).

The system uses two filtered Xenon lamps in a hand held applicator that generates a pulse of light. The principle of selective photo-thermolysis has been studied extensively for hair removal applications due to the wide-spread use of this technology. The light penetrates the epidermis to the dermis and is absorbed by the Melanin in the hair shaft. The heat generated by the absorbed light is then dissipated to the follicle and generates local thermal damage, causing reduction in hair growth. This technology has been used extensively by many devices to generate hair loss using lasers or filtered Xenon flash lamps, mostly in doctors' offices.

Present methods of hair removal:

1. Depilators
2. Shaving
3. Waxing
4. Laser or light based treatments.

The proposed method is based on the principal of selective photo-thermolysis using pulse light in (b) (4) optical fluences of up to 5 J/cm². Light spectrum is set by the optical filter to 475-1200 nm, which is highly absorbed by the hair shaft Melanin. The use of (b) (4) energy is expected to be safe even on darker skin types, and the ability to control the energy output will allow the user safe and effective treatment.

Name and Intended Use of Device

Proprietary name: Silk'n

Common name: Intense pulse light system

Intended Use: The device is indicated for body hair removal (b) (4)

Purpose and Objectives

Purpose: The purpose of this clinical investigation is to check the efficiency and safety of the Silk'n hair removal device.

Objectives: The objective of the study is to determine the ability of inducing hair loss with a low power photo-epilation system.

Scope and Duration of the Study

Scope: The clinical study will include (b) (4) subjects at (b) (4)

Duration: Eligible subjects will receive (b) (4) treatments according to the study protocol and will be followed for a period (b) (4) document safety and efficacy.

Substantial equivalence and risk level

The overall specifications, principle of operation, performance characteristics of the Silk'n device are substantially equivalent to those of the predicate device- Silk'n (Home Skinovations Ltd., K072906). The Silk'n, similar to its predicate device, is a non significant risk device. None the less, in order to establish the margins of the treatment parameters and determine the most effective and case specific parameters, we further wish to examine and optimize the safety and efficacy of the Silk'n device per clinical indication.

Subject Selection

Investigator will determine patients eligibility to participate in the study according to inclusion / exclusion criteria.

Inclusion Criteria

- 1- Presence of unwanted hairs on the body, (b) (4)
- 2- Skin Type (b) (4)
- 3- Males and females older than 21 years of age but not more than 60 years of age.
- 4- If female, must be either post-menopausal or surgically sterilized, or using a medically acceptable form of birth control (i.e., oral contraceptives, IUD, contraceptive implant, barrier methods with spermicide or abstinence).
- 5- Informed consent agreement by the subject.
- 6- Willingness to follow the treatment schedule and post treatment care

Exclusion Criteria

1. Malignant or pre-malignant pigmented lesions in the area to be treated
2. Scarring or infection of the area to be treated
3. Known photosensitivity

4. Pregnancy
5. Subjects with Diabetes (Type I or II)
6. Presence of a suntan in the area to be treated
7. Use of medication known to induce photosensitivity
8. Subject is on anticoagulative medication or thromboembolic condition
9. Subjects with a pacemaker or internal defibrillator
10. Use of NSAIDS two weeks prior to, and two weeks following the treatment
11. Subjects that use waxing or other methods of photo-epilation 3 months prior to treatment

Note: the device is not to be used on the face

Discontinuation

Subjects enrolled in the study will be free to discontinue their participation in the study at any time. A decision to discontinue participation will not prejudice their medical care. In such instances, Investigators will attempt to obtain clinical results and side effect evaluation concerning the subject prior to his/ her withdrawal.

Clinical Procedure

Pre-Procedure Evaluation

- 1- The Investigator will discuss with the subject many available treatment options including the Silk'n Hair Removal System and other treatment modalities.
- 2- During the first visit, the investigator will obtain an informed consent from the subject, clearly indicating his/her understanding of the requirements and risks involved with study participation, including but not limited to:
 - a. (b) (4)

- (b) (4)
- b.
 - c.
 - d.
 - e.
 - f.
 - g.
 - h.

3- A full detailed subject chart will record initial hair count in the area to be treated. Initial hair count will be used to assess treatment efficacy as the treatment progresses.

4- (b) (4)

If the subject meets all inclusion and exclusion criteria, the investigator will obtain signed Informed Consent, enroll the subject into the study, and record required data on the appropriate study data form. The pre-treatment data will be stored in the computer under the subject study number ID.

A copy of the Informed Consent Form (ICF) and a copy of a Case Report Form (CRF) are attached to this protocol.

Treatment

The treatment will be performed in the investigator clinic by a trained Nurse or technician.

Treatment Parameters

Energy level should match the subject skin type. Test pulse energy should be set according with the attached table:

| | Energy level for test pulse |
|-------------------|-----------------------------|
| Skin type (b) (4) | (b) (4) |
| Skin type | |
| Skin type | |

Treatment test spot

Treatment test pulse should be performed on the darker part of the area that is going to be treated by the Silk'n. If for example the arm is to be treated, the test pulse should be performed on the outer part of the arm where it is usually darker than the inner part of the arm.

If no skin reaction was noted 15 minutes after the test pulse, (see test pulse instruction in the next section), a second test pulse should be performed at one level higher than the first pulse. If no skin reaction was noted perform the treatment.

Treatment

To perform the Treatment, follow steps 1-6 and wait 15 minutes to check if any side effects are noted:

1. In the area to be treated, hairs should be trimmed to 3/32 inch length,(1-2 mm), or shaved 3 days prior to the treatment.
2. Clean the skin from the trimmed hairs.
3. Hold the applicator and ensure that the light output window is clean.
4. Switch on the system and set the energy level.
5. Place the Applicator on the treatment area and apply slight pressure. The light window should be in good contact with the skin, press the trigger switch to emit light pulse.

6. Move the applicator to the next spot to ensure full coverage of the area to be treated.
7. Repeat the test spot procedure for 3 pulses.
8. If side effects are noted after 15 minutes reduce the energy level by one level to lower setting.
9. Transient erythema and follicular redness are normal reactions as well as smell of burned hairs.

If no skin reaction noted after 15 minutes complete the treatment to cover entire area to be treated.

After the session is finished the applicator should be clean using isopropyl alcohol.

Post Treatment Care.

The investigator should supply moisturizer for subjects use after the treatment.

Use high factor sunscreen (SPF of 15 or greater), and protect the treated area from sunlight for at least 2 days after the treatment. Tanning after a treatment may cause hyper-pigmentation.

Repetitive treatments

(b) (4)



Data Analysis

Hair loss at each follow up session will be calculated as the ratio of the number of hairs after the treatment compared with the hair count prior to treatment. Hair count will be performed by the investigator or clinic staff.

Materials

Silk'n Hair Removal Device will be placed at the participating institution, maintained and serviced as needed. Home Skioventions Ltd. will provide the system to the institutions for the duration of the clinical investigation.

Training Requirements

Both the Investigator and the Sponsor, prior to any independent use of the device, will agree as to the Investigator's training requirements. Prior to the study, the sponsor will ensure that the investigators have received training regarding use of the device.

Risk/ Benefit Analysis

Risks

The potential risks for adverse effects of the treatment procedure include, but are not limited to, (b) (4)

(b) (4)

Benefits

Potential benefits include reduction of unwanted hair from the skin.

Alternatives

Present methods of skin treatments are by laser or light based devices as well as epilation, waxing, shaving, and depilators. Published studies confirm that these systems provide a reliable, usually complete treatment. Temporary dyschromias, and in rare instances, scarring have occurred as a result of this treatment.

Informed Consent

The Investigator will obtain written Informed Consent prior to the anticipated scheduled procedure. The investigator will inform the subject as to the experimental procedure to be applied.

The investigator will assure the subject that his or her decision regarding participation in the study will have no bearing on the quality of medical care received and that the decision whether to participate in or withdraw from the study is strictly voluntary.

If the system is used on a subject without obtaining an Informed Consent, the Investigator must report the use to the Sponsor and the reviewing IRB.

Records and Reports

Data from the pre-treatment evaluation, the treatment procedure (complication, if any), the post-procedural evaluations and all follow-up evaluations will be recorded directly on the clinical data forms prepared for the study.

All specified records and reports concerning the clinical study will be forwarded to the Sponsor within a reasonable amount of time. The Investigators shall retain copies of all documentation for a period of two years following the date on which the entire clinical investigation is terminated or discontinued.

Confidentiality

This study is confidential. Information made public regarding subjects enrolled in this study will not disclose the name of any individual.

Monitoring Procedures

The study will be monitored by representatives of Home Skinovations Ltd. by telephone, in writing and during on-site visits. At the very least, site visits will be scheduled prior to the initiation of the study, on the occasion of the initial use of investigational devices, and at the end of the study. The purpose of site visits will be to ensure compliance with the investigational plan, to ensure appropriate use of investigational devices, and to inspect and retrieve study data.

Investigator Certification

Home Skinovations Ltd. certifies that each Investigator who takes part in this study will read and sign the Investigator's Agreement. Original signed Agreements will be retained

by Home Skinovations Ltd. No Investigators will be added to the investigation until they have signed the Agreement.

APPENDIX - 1: STUDY SCHEMATICS

| | | | |
|---------------------------|---------|---|---|
| Visit# | (b) (4) | | |
| | | | |
| Week# | (b) (4) | | |
| | | | |
| Informed Consent Form | X | | |
| Inclusion/ Exclusion | X | | |
| Medical History | X | | |
| Test spot | X | | |
| Treatment | X | X | |
| Hair Count | X | | X |
| Photography | X | | X |
| Adverse Events Evaluation | X | X | X |

Appendix 3 - Case report form

Visit 1 – Study Entry Data

Patient Id: _____

Date of this assessment: _____/_____/_____

Subject Demographics:

Date of Birth

: _____/_____/_____

Male/Female

Gender (Enter appropriate response): _____

Subject History (Enter appropriate response):

Signed Informed Consent Form? _____

Skin Type: _____

Deep suntan present? _____

Known photosensitivity? _____

Pregnancy? _____

Use of acceptable birth control method? _____

Date of last menstrual period? _____

With Diabetes (Type I or II)? _____

Use of medication known to induce photosensitivity? _____

Scarring/ Infection present? _____

Known anticoagulative or thromboembolic condition? _____

Use of Accutane within past 6 months? _____

Use of aspirin or other NSAIDs two weeks prior to, and for two weeks following the treatment? _____

With a pacemaker or internal defibrillator? _____

Older than 21 years? _____

All questions above should be answered for subject to be enrolled in this study.

Other Relevant Medical Conditions or Further Comments on the Above:

Investigator Initials: _____ Date: _____

Visit 1 – Subject Medical History

Date of this assessment: ____/____/____

Additional relevant medical history, if any

Check here if none _____

| Relevant Medical History | Check if condition is ongoing |
|--------------------------|-------------------------------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

Comments:

Investigator Initials: _____ Date: _____

Treatment follow up form

Date of this assessment: ____/____/____

Investigator Assessment:

Clinical Observation: 0 = absent, 1 = mild, 2 = moderate, 3 = severe

Target Area : _____

Skin type

| | Erythema | Edema | comments |
|---------|----------|-------|----------|
| (b) (4) | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

INFORMED CONSENT FORM

Study Title: Clinical Study to Determine the Safety and Efficacy of the Silk'n Home Pulsed Light Device

Study Doctor

Name

Address

Telephone Number

INTRODUCTION:

This Informed Consent form is being presented to you, because after a discussion about this research study and being verbally informed, you have indicated that you are interested in participating in this research study for a device (Silk'n Home Use Hair Removal Device) for hair removal. The purpose of this document is to inform you in writing of the purpose of the research, a description of the procedures to be followed, possible risks and discomforts, possible benefits, and of basic ground rules which will govern this research project. After reading this from and having all your question answered, if you decide to participate, you should return this consent form to the study doctor's office, sign this form on the last page, initial and date each prior page in the presence of the study staff. Please read this form carefully before you make your decision. You may refuse to participate in this study and that decision will not be held against you, nor will it change any matters between you and this office, and you will continue to receive the same level of health care.

Subjects who may not participate are those who have:

- Scarring in the proposed treatment area.
- Have had a prior treatment of the area with another laser or light-based device.
- Are pregnant, or planning a pregnancy.
- Have diabetes.

(b) (4)

Protocol #: (b) (4)

Page 1 of 7

Subject's initials _____

Date _____

- Have a suntan in the area.
- Are using a medication that sensitizes skin to light or are taking blood-thinners or use pain relievers on a regular basis such as aspirin, Motrin®, Advil®, Naprosyn®, Aleve® two weeks prior to, and two weeks following the treatment.

PURPOSE:

The purpose of this research study is to determine the safety and effectiveness of the use of Silk'n device, for hair removal. The system is based on the principle of selective photothermolysis (breakdown of the hair follicle due to heat from intense pulse light), that has been studied and used extensively for hair removal uses.

The light penetrates the top layer of the skin to the deeper layer and is absorbed by the melanin (dark pigment) in the hair shaft. The heat generated by the absorbed light is then spread to the follicle and generates local heat damage, causing reduction in hair growth.

(b) (4)

PROCEDURES:

If you agree to participate (b) (4) study, you will be required to make as many as (b) (4)

(b) (4) receiving up to (b) (4), according to the study doctor's decision, and the achieved results. The first visit will determine if you are eligible to be in this study. (b) (4) prior to the first treatment, all other treatments of hair removal are prohibited in the area to be treated. Females of child-bearing potential must use of adequate birth control (i.e. oral contraceptives, IUD, contraceptive implant, barrier methods with spermicide or abstinence). You will be asked to provide a detailed medical history (particularly about any medications that affect the skin), have your skin and hairs examined to determine what skin type you have and to select the area for the study procedure. The actual treatment session should last approximately 20 minutes. All treatments will be administered by nurse or trained technician. Test spot will be performed on the area where the actual treatment is going to be (b) (4)

(b) (4). The test spot involves 3 applications and follow up of 15 minutes to assure that no adverse reaction is noted.

Version Date: (b) (4)

Protocol #: (b) (4)

Page 2 of 7

Subject's initials _____

Date _____

Before the treatment, photographs will be taken of the designated treatment area. These pictures may be used by the Sponsor (Home Skinovations Ltd.) in the future for educational and marketing purposes. Your identity will be masked in the picture. At the first treatment session, a hair count will be performed, and the best settings for the system output will be determined before the full treatment followed by a rest period to observe those effects. The area to be treated should be shaved 3 days prior to treatment. Afterwards, you will be given instructions on the care of the treated area and things to avoid.

In the event were any side effect are noted, a follow-up appointment will be made for the first or the second day after the treatment to re-examine your skin. You will be asked about any discomfort and reactions you may be experiencing. Another photograph of the treated area may be taken. The visit should last approximately 20 minutes.

(b) (4)

. An examination of the treated area will be done to evaluate your skin's condition. Each visit should last approximately 20 minutes.

(b) (4)

Photographs of the treated skin and a hair count will be taken. Each of the follow up sessions should last approximately 15 minutes.

(b) (4)

, unless you are experiencing a complication. If so, you will be asked to return for further examinations until it resolves. A photograph of the treated area may be taken and the results will be assessed. Your participation in the study will then end.

RISKS AND DISCOMFORTS:

The device may cause the following skin reactions: redness, pain, blistering, crusting, swelling, decrease or increase of skin color and scarring. There is a chance the area could become infected (especially if it is not kept clean or you pick at any crusts or scabs). The study doctor will advise you in a set of instructions how to avoid such complications. An antibiotic cream or ointment may be needed. If so, the study doctor will provide it to you at no charge. You should report any reactions or problems with the treated site to the study doctor promptly. Some risks are

Version Date (b) (4)

Protocol #: (b) (4)

Page 3 of 7

Subject's initials _____

Date _____

unforeseeable. Any new findings that develop during the course of this investigation, which may relate to your willingness to continue your participation, will be provided to you.

BENEFITS:

You may experience hair loss in the area that is treated. There is no guarantee that you may experience any improvement. Your participation may help in the development of a device for the treatment of unwanted hair.

ALTERNATIVES:

Current methods available for hair removal treatments include:

1. Depilators
2. Shaving
3. Waxing
4. Laser or light-based treatments
5. Electrolysis

You do not have to participate in this research study to receive hair removal treatment.

CONTACTS:

If you experience a study-related problem, or if you have any questions at any time about the study, contact Dr. [name] at [telephone number].

If you have any questions about your rights as a research participant or related concerns or complains, you may contact the (b) (4)

(b) (4). The (b) (4) is a committee that has reviewed this research study to help ensure that your rights and welfare as a research participant are protected and that the study will be carried out in an ethical manner. Review and approval of this study by the (b) (4) is not an endorsement of the study or its outcome.

RESEARCH PARTICIPATION INFORMATION

You can obtain information about participating in research studies from a number of sources. A few are:

Version Date: (b) (4)

Protocol #: (b) (4)

Page 4 of 7

Subject's initials _____

Date _____

- Center for Information and Study on Clinical Research Participation (CISCRP):
www.ciscrp.org
- Food and Drug Administration (FDA): www.fda.gov
- Office for Human Research Protections (OHRP): www.hhs.gov/ohrp
- National Institute of Health: clinical trials.gov
- National Cancer Institute: www.nci.nih.gov
- CenterWatch: www.centerwatch.com
- Various large university websites
- Various associations and societies concerned with specific diseases websites.”

FINANCIAL:

To help defray the costs of your participation, you will be paid (b) (4) for each visit. If you complete the study, you will be paid an additional bonus of (b) (4). The maximum compensation for subjects participation in the study is (b) (4).

COSTS:

There will be no costs to you for office visits, examinations, antibiotic cream, and procedures as part of this study.

COMPENSATION FOR INJURY:

If a research-related injury occurs, the sponsor, Home Skinovations Ltd. will pay for any treatment of research-related injury, which will be provided by the study doctor. However, the sponsor will not pay for pre-existing conditions or for any conditions arising after the study. Also, you will not receive compensation for wages associated with lost-time at your workplace. However, by signing this form you have not given up any of your legal rights.

CONFIDENTIALITY AND AUTHORIZATION TO USE AND DISCLOSE PERSONAL HEALTH INFORMATION:

A federal regulation called the “Health Insurance Portability and Accountability Act” (HIPAA), describes how your personal health information may be used, disclosed and made accessible to you. This privacy rule is designed to protect the confidentiality of your personal health information.

Version Date (b) (4)

Protocol #: (b) (4)

Page 5 of 7

Subject's initials _____

Date _____

This study can be performed only by collecting and using your personal health information. Your study records will be kept as confidential as possible under local, state and federal laws. Personnel from the following organizations may examine your study records: the sponsor, [Home Skinovations Ltd.], personnel associated with this study, regulatory agencies, such as the

(b) (4)

Because of the number of individuals who may see your records, absolute confidentiality cannot be guaranteed.

Personal health information that may be used and disclosed includes that which is obtained to determine your eligibility to participate and that which is collected from the procedures that are carried out. It may identify you by name, address, telephone number, Social Security Number, study number, date of birth or other identifiers. Once the information is disclosed, it is possible that it may be re-disclosed, at which time it may no longer be protected by federal regulations, but may be by state laws. If the final study data are prepared for publication and other reports, your identity will not be revealed. Under these federal privacy regulations, you have the right to see and copy any of the information gathered about you, until your study records are no longer kept by the study doctor. However, it may not be available until the study has been completed.

You may, by written notice to the study doctor, cancel your authorization to use or disclose your personal information at any time. If you withdraw your authorization, the information collected up to that time may still be used to preserve the scientific integrity of the study. By signing this consent form, you authorize these uses and disclosures of your personal information. If you do not authorize these uses and disclosures, you will not be able to participate in the study. This authorization does not have an expiration date.

VOLUNTARY PARTICIPATION/ WITHDRAWAL:

Your participation in this study is voluntary and, if you decide to withdraw at any time, you may do so without penalty or giving up any benefits to which you are otherwise entitled. You may be discontinued from the study by the study doctor for reasons of, but not limited to, (1) a severe adverse reaction, (2) the study is stopped, (3) your failure to follow instructions, (4) the study doctor determines it is your best interest.

Version Da (b) (4)

Protocol #:

Page 6 of 7

Subject's initials _____

Date _____

PHOTOGRAPH CONSENT

Photographs and digital images will be taken of your treated areas as a requirement for participating in this study. These photographs and digital images will not reveal your identity, but they may not be used without your written consent. By checking the box below, you give your permission to use photographs and digital images of your treated areas, which were taken during this study, for the purposes of professional publications, training, education or sales. You will automatically be withdrawn from the study, should you choose not to give your consent at this time.

I may cancel or withdraw my photography consent at any time during or after my participation in this study. I must do so in writing and submit that withdrawal of consent to Dr. _____ . I agree to have photographs taken as a part of this study as described in this paragraph.”

☐ YES _____ (initial)

☐ NO _____ (initial)

CONSENT:

I have read this consent form, have had any questions answered and understand the purpose, procedures, risks and other aspects of the study. I voluntarily agree to participate in this research study. I will be given a signed and dated copy of this consent form.

Subject's Name (printed)

Date: _____

Subject's Signature

Person conducting consent process

Date: _____

Study Doctor's Signature

Version Date: (b) (4)

Page 7 of 7

Protocol #: (b) (4)

Subject's initials _____

Date _____

POST-TREATMENT INSTRUCTIONS

Following Treatment:

There may be a stinging or sunburn sensation accompanying your treatment. This will subside over the next several hours and can be relieved with the application of a cool compress or an ice pack to the treated area.

Immediately following treatment with the device, there will be generalized redness, often accompanied by smell of burned hairs. These effects are generally temporary and should resolve over the next several hours.

Precautions:

Avoid direct sun exposure or tanning beds between treatments. A sunscreen with an SPF 15 or higher should be used when sun exposure cannot be avoided.

Avoid injury to the area. Do not rub, scratch, or pick at the treated skin. Treating the skin gently will decrease the chance of any problems.

Make-up application over the treatment area should be avoided if any scaling or crusting appears.

Shower or bathe as usual with a mild soap, but remember the treated area may be temperature sensitive. Avoid the use of exfoliants, loofa sponges, and aggressive scrubbing of the treatment area. While the area is still red or healing, use a soft cloth to pat the area dry.

Contact this office immediately; Dr. [Name], [Telephone number] with any questions.

**Clinical Study to Determine the Long Term Efficacy of the
Silk'n Home Pulsed Light Device**

Clinical Investigation Plan follow-up Study

Date: (b) (4)

Protocol number: (b) (4)

Version (b) (4)

Sponsor:

Home Skinovations Ltd.
P.O.B. 533 Yokneam Illit
20692 Israel
Tel: 972-4-9097471
Fax: 972-4-9097470

Home Skinovations Ltd. Proprietary and confidential Information

Page number: 1

(b) (4)

(b) (4)

(b) (4)

Protocol Summary

This is a clinical study to determine the long term efficacy of a low energy intense pulsed light, (IPL) device for hair removal. The study involves one follow up session on subjects that participated in the HR-2 study.

The system uses two filtered Xenon lamps in a hand held applicator that generates a pulse of light. The principle of selective photo-thermolysis has been studied extensively for hair removal applications due to the wide-spread use of this technology. The light penetrates the epidermis to the dermis and is absorbed by the Melanin in the hair shaft. The heat generated by the absorbed light is then dissipated to the follicle and generates local thermal damage, causing reduction in hair growth. This technology has been used extensively by many devices to generate hair loss using lasers or filtered Xenon flash lamps, mostly in doctors' offices.

Name of Device Used

| | |
|-------------------|--|
| Proprietary name: | Silk'n |
| Common name: | Intense pulse light system |
| Intended Use: | The device is indicated for permanent hair reduction |

Purpose and Objectives

| | |
|----------|---|
| Purpose: | The purpose of this clinical investigation is to check the long term efficiency of the Silk'n for permanent |
|----------|---|

hair reduction.

Objectives:

The objective of the study is to determine the ability of inducing permanent hair reduction with a low power photo-epilation system.

Scope and Duration of the Study

Scope:

The clinical study will include (b) (4)

(b) (4)

Duration:

Risk level

There is no risk involved in this study as no treatment will be performed. (b) (4)

(b) (4)

Subject Selection

Investigator will ask all participants in the previous (b) (4) to be part of this study.

Inclusion Criteria

- 1- Participation in the (b) (4)

Exclusion Criteria

- 2- (b) (4)

(b)

Discontinuation

Subjects enrolled in the study will be free to discontinue their participation in the study at any time. A decision to discontinue participation will not prejudice their medical care. In such instances, Investigators will attempt to obtain clinical results and side effect evaluation concerning the subject prior to his/ her withdrawal.

Procedure

Evaluation

If the subject meets all inclusion and exclusion criteria, the investigator will obtain signed Informed Consent, enroll the subject into the study, and record required data on the appropriate study data form. The evaluation data will be stored in the computer under the subject study number ID.

A copy of the Informed Consent Form (ICF) and a copy of a Case Report Form (CRF) are attached to this protocol.

Procedure

(b) (4)

Data Analysis

(b) (4)

Risks

(b) (4)

(b) (4)

Informed Consent

The Investigator will obtain written Informed Consent prior to the anticipated scheduled procedure. The investigator will inform the subject as to the evaluation procedure to be applied.

The investigator will assure the subject that his or her decision regarding participation in the study will have no bearing on the quality of medical care received and that the decision whether to participate in or withdraw from the study is strictly voluntary.

If the evaluation is performed on a subject without obtaining Informed Consent, the Investigator must report the violation to the Sponsor and the reviewing IRB.

Records and Reports

Data from the evaluation, the procedure (complication, if any), the post-procedural evaluations and all follow-up evaluations will be recorded directly on the clinical data forms prepared for the study.

All specified records and reports concerning the clinical study will be forwarded to the Sponsor within a reasonable amount of time. The Investigators shall retain copies of all documentation for a period of two years following the date on which the entire clinical investigation is terminated or discontinued.

Confidentiality

This study is confidential. Information made public regarding subjects enrolled in this study will not disclose the name of any individual.

Monitoring Procedures

The study will be monitored by representatives of Home Skinovations Ltd. by telephone, in writing and during on-site visits. At the very least, site visits will be scheduled prior to the initiation of the study, on the occasion of the initial use of investigational devices, and at the end of the study. The purpose of site visits will be to ensure compliance with the

investigational plan, to ensure appropriate use of investigational devices, and to inspect and retrieve study data.

Investigator Certification

Home Skinovations Ltd. certifies that each Investigator who takes part in this study will read and sign the Investigator's Agreement. Original signed Agreements will be retained by Home Skinovations Ltd. No Investigators will be added to the investigation until they have signed the Agreement.

APPENDIX - 1: STUDY SCHEMATICS

| | |
|---------|---|
| (b) (4) | |
| | |
| (b) (4) | |
| | |
| (b) (4) | X |
| | X |
| | X |
| | X |
| | |

Appendix 3 - Case report form

Visit 1 – Study Entry Data

Patient Id: _____

Date of this assessment: _____/_____/_____

Subject Demographics:

Date of Birth

: ____/____/____

Gender (Enter appropriate response):

Male/Female

Subject History (Enter appropriate response):

Signed Informed Consent Form? _____

Skin Type: _____

Deep suntan present? _____

Known photosensitivity? _____

Pregnancy? _____

Use of acceptable birth control method? _____

Date of last menstrual period? _____

With Diabetes (Type I or II)? _____

Use of medication known to induce photosensitivity? _____

Scarring/ Infection present? _____

Known anticoagulative or thromboembolic condition? _____

Use of Accutane within past 6 months? _____

Use of aspirin or other NSAIDs two weeks prior to, and for two weeks following the treatment? _____

With a pacemaker or internal defibrillator? _____

Older than 21 years? _____

All questions above should be answered for subject to be enrolled in this study.

Other Relevant Medical Conditions or Further Comments on the Above:

Investigator Initials: _____ Date: _____

Hair count:

Home Skinovations Ltd. Proprietary and confidential Information

Page number: 8

Date: (b) (4)

Protocol number: (b) (4)

(b) (4)

INFORMED CONSENT FORM

Study Title: Clinical Study to Determine the Long Term Efficacy of the Silk'n Home Pulsed Light Device

Study Doctor

Name

Address

Telephone Number

INTRODUCTION:

This Informed Consent form is being presented to you, because after a discussion about this research study and being verbally informed, you have verified that you previously participated in the (b) (4) research study for a device (Silk'n Home Use Hair Removal Device) for permanent hair removal. The purpose of this document is to inform you in writing of the purpose of the research, a description of the procedures to be followed, possible risks and discomforts, possible benefits, and of basic ground rules which will govern this research project. After reading this form and having all your question answered, if you decide to participate, you should return this consent form to the study doctor's office, sign this form on the last page, initial and date each prior page in the presence of the study staff. Please read this form carefully before you make your decision. You may refuse to participate in this study and that decision will not be held against you, nor will it change any matters between you and this office, and you will continue to receive the same level of health care.

Subjects who may not participate are those who have:

(b) (4)

Version Da (b) (4)

Protocol #:

Page 1 of 5

Subject's initials _____

Date _____

PURPOSE:

The purpose of this research study is to determine the long term effectiveness of the use of Silk'n device, for permanent hair reduction.

This study will include (b) (4)

(b) (4)

PROCEDURES:

If you agree to participate in this one follow-up visit, you will be required to come to the study doctor's office. (b) (4)

(b) (4)

Photographs will be taken of the designated treatment area. These pictures may be used by the Sponsor (Home Skinovations Ltd.) in the future for educational and marketing purposes. Your identity will be masked in the picture.

The examination should last approximately 15 minutes. After the visit your participation in the study will end.

ALTERNATIVE TO PARTICIPATION

(b) (4)

RISKS AND DISCOMFORTS:

(b) (4)

BENEFITS:

Your participation may help in the development of a device for the treatment of unwanted hair.

CONTACTS:

If you experience a study-related problem, or if you have any questions at any time about the study, contact Dr. [name] at [telephone number]. If you have any questions about your rights as a research participant or related concerns or complains, you may contact the (b) (4)

(b) (4)

You may also

(b) (4)

Version Date: (b) (4)

Protocol #: (b) (4)

Subject's initials _____

Date _____

study to help ensure that your rights and welfare as a research participant are protected and that the study will be carried out in an ethical manner. Review and approval of this study by the Essex IRB is not an endorsement of the study or its outcome.

RESEARCH PARTICIPATION INFORMATION

You can obtain information about participating in research studies from a number of sources. A few are:

- Center for Information and Study on Clinical Research Participation (CISCRP): www.ciscrp.org
- Food and Drug Administration (FDA): www.fda.gov
- Office for Human Research Protections (OHRP): www.hhs.gov/ohrp
- National Institutes of Health: www.nih.gov
- www.womancando.org
- National Cancer Institute: www.nci.nih.gov
- CenterWatch: www.centerwatch.com
- Various large university websites
- Various associations and societies concerned with specific diseases websites."

FINANCIAL:

To help defray the costs of your participation, you will be paid \$50 for your visit.

COSTS:

There will be no costs to you for the office visit, examination, or photographs as part of this study.

CONFIDENTIALITY AND AUTHORIZATION TO USE AND DISCLOSE PERSONAL HEALTH INFORMATION:

A federal regulation called the "Health Insurance Portability and Accountability Act" (HIPAA), describes how your personal health information may be used, disclosed and made accessible to you. This privacy rule is designed to protect the confidentiality of your personal health information.

This study can be performed only by collecting and using your personal health information. Your study records will be kept as confidential as possible under local, state and federal laws. Personnel from the following organizations may examine your study records: the sponsor,

Version Date (b) (4)

Protocol #:

Page 3 of 5

Subject's initials _____

Date _____

134

[Home Skinovations Ltd.], personnel associated with this study, regulatory agencies, such as the

(b) (4)

of the number of individuals who may see your records, absolute confidentiality cannot be guaranteed.

Personal health information that may be used and disclosed includes that which is obtained to determine your eligibility to participate and that which is collected from the procedures that are carried out. It may identify you by name, address, telephone number, Social Security Number, study number, date of birth or other identifiers. Once the information is disclosed, it is possible that it may be re-disclosed, at which time it may no longer be protected by federal regulations, but may be by state laws. If the final study data are prepared for publication and other reports, your identity will not be revealed. Under these federal privacy regulations, you have the right to see and copy any of the information gathered about you, until your study records are no longer kept by the study doctor. However, it may not be available until the study has been completed.

You may, by written notice to the study doctor, cancel your authorization to use or disclose your personal information at any time. If you withdraw your authorization, the information collected up to that time may still be used to preserve the scientific integrity of the study. By signing this consent form, you authorize these uses and disclosures of your personal information. If you do not authorize these uses and disclosures, you will not be able to participate in the study. This authorization does not have an expiration date.

VOLUNTARY PARTICIPATION/ WITHDRAWAL:

Your participation in this study is voluntary and, if you decide to withdraw at any time, you may do so without penalty or giving up any benefits to which you are otherwise entitled. You may be discontinued from the study by the study doctor for reasons of, but not limited to, (1) a severe adverse reaction, (2) the study is stopped, (3) your failure to follow instructions, (4) the study doctor determines it is your best interest.

This Space Intentionally Left Blank

Version Date (b) (4)

Protocol #:

Subject's initials _____

Date _____

Page 4 of 5

PHOTOGRAPH CONSENT

Photographs and digital images will be taken of your treated areas as a requirement for participating in this study. These photographs and digital images will not reveal your identity, and they may not be used without your written consent. By checking the box below, you give your permission to use photographs and digital images of your treated areas, which were taken during this study, for the purposes of professional publications, training, education or sales. You will automatically be withdrawn from the study, should you choose not to give your consent at this time.

I may cancel or withdraw my photography consent at any time during or after my participation in this study. I must do so in writing and submit that withdrawal of consent to Dr. _____. I agree to have photographs taken as a part of this study as described in this paragraph.”

☐ YES _____ (initial)

☐ NO _____ (initial)

CONSENT:

I have read this consent form, have had any questions answered and understand the purpose, procedures, risks and other aspects of the study. I voluntarily agree to participate in this research study. I will be given a signed and dated copy of this consent form.

Subject's Name (printed)

Date: _____

Subject's Signature

Person conducting consent process

Date: _____

Study Doctor's Signature

Version Date (b) (4)

Protocol #: _____

Subject's initials _____

Date _____

Appendix 1

Flash N Go (cleared) submission and all communication
with FDA

HomeSkinovations

Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
www.homeskinovations.com

(b) (4)

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, MD 20850

Re: Letter dated September, (b) (4)

(b) (4)

Thank you for your letter dated (b) (4) in response to our 510K submission (b) (4)

Please find the requested information according to the order asked in your letter:

1. (b) (4)

2. Software level of concern and all applicable software document are presented in appendix 1.

3. (b) (4)

4. (b) (4)

Therefore no new questions regarding biocompatibility should rise.

5. The Flash N Go should be used on clean skin area therefore the hairs should be shaved or trimmed prior to the procedure and the skin surface should be clean. As part of the clinical study with the Silk'n patient

HomeSkinovations


Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
www.homeskinovations.com

performed all part of the protocol by themselves including shaving/
trimming the area to be treated by the Silk'n.
6-15 The user manual was changed according to your recommendation and
presented in appendix 3.

We trust that the information included in this application will be adequate to allow for
its prompt review. Please contact me if there are any questions.

Sincerely yours,

Dr. Amir Waldman
VP regulatory affairs
Home Skinovations Ltd.





Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
www.homeskinovations.com

Appendix 1

Software level of concern and all required documents

(b) (4)

Flash N' Go Software Level of Concern

- Light pulse generation is enabling by hardware, no software fault can activate light pulse.
- Software fault may cause only higher light energy output pulse than requested by user. The max light energy output is limits by hardware.
- Color Sensor malfunction, software disable the light pulse.

Max light energy pulse energy is limits by hardware to the maximum charging voltage in case of software fault. (b) (4)

(b) (4)

In the most vulnerable user (dark skin type) were the energy setting was lowest (b) (4), a software fault may deliver double the energy. This may results in superficial burn with pigmentation changes and possible minor scarring.

Determining Level of Concern

According to the "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices"

| Major Level of Concern (If the answer to any <u>one</u> question below is Yes, the Level of Concern for the Software Device is likely to be Major) | | |
|---|---|--------|
| No. | Question | Answer |
| 1 | Does the Software Device qualify as Blood Establishment Computer Software? | No |
| 2 | Is the Software Device intended to be used in combination with a drug or biologic? | No |
| 3 | Is the Software Device an accessory to a medical device that has a Major Level of Concern? | No |
| 4.a | Does the Software Device control a life supporting or life sustaining function? | No |
| 4.b | Does the Software Device control the delivery of potentially harmful energy that could result in death or serious injury, such as radiation treatment systems, defibrillators, and ablation generators? | No |
| 4.c | Does the Software Device control the delivery of treatment or therapy such that an error or malfunction could result in death or serious injury? | No |
| 4.d | Does the Software Device provide diagnostic information that directly drives a decision regarding treatment or therapy, such that if misapplied it could result in serious injury or death? | No |
| 4.e | Does the Software Device provide vital signs monitoring and alarms for potentially life threatening situations in which medical intervention is necessary? | No |

| Moderate Level of Concern (If the Software Device is not Major Level of Concern and the answer to any one question below is Yes, the Level of Concern is likely to be Moderate.) | | |
|---|---|--------|
| No. | Question | Answer |
| 1 | Is the Software Device an accessory to a medical device that has a Moderate Level of Concern? | No |
| 2 | Prior to mitigation of hazards, could a failure of the Software Device result in Minor Injury, either to a patient or to a user of the device? | Yes |
| 3 | Could a malfunction of, or a latent design flaw in, the Software Device lead to an erroneous diagnosis or a delay in delivery of appropriate medical care that would likely lead to Minor Injury? | No |

Conclusion:

The Flash N' Go software level of concern is: **Moderate!**



Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go

Software Verification and Validation Plan

DO1009131

| | Role | Name | Signature | Date |
|---------|------|------|-----------|------|
| (b) (4) | | | | |

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Table of Contents

| | | |
|-----|---|----|
| 1 | Introduction and overview..... | 4 |
| 1.1 | Purpose of this document | 4 |
| 1.2 | Scope of the development Project | 4 |
| 1.3 | Definitions, Acronyms and Abbreviations | 4 |
| 1.4 | References | 4 |
| 1.5 | Overview of Document..... | 5 |
| 2 | Reviews; walkthroughs; inspections; and audits | 5 |
| 2.1 | Schedule of reviews, etc. | 5 |
| 2.2 | Procedure for reviews, etc. | 5 |
| 3 | Component Test Plans and Procedures | 5 |
| 3.1 | Component and Integration Test strategy overview | 5 |
| 3.2 | Component Test description | 6 |
| 3.3 | Integration test description | 8 |
| 4 | System Test Plans and Procedures | 14 |
| 4.1 | System test Strategy Overview | 14 |
| 4.2 | System Test Description | 14 |
| 5 | Defect tracking plans..... | 15 |
| 6 | Traceability from SRS to SDD | 15 |
| 7 | Test requirements cross-reference Matrix | 15 |
| 8 | Acceptance Test and Preparation for delivery | 15 |
| 8.1 | Procedure by which the software product will be acceptance tested | 15 |
| 8.2 | Specific acceptance criteria | 15 |
| 8.3 | Scenario by which the software product will be installed..... | 16 |

1 Introduction and overview

1.1 Purpose of this document

The purpose of this document is to provide an overview of the design and development of the Flash N' Go firmware. It will serve as a guide to the developers satisfying the client requirements and the specifications. This document will also be used in the verification and validation of the project.

1.2 Scope of the development Project

This project, named Flash N' Go, will consist of Base unit that contain the MCU firmware and the Applicator which contain the Disposable Lamp Cartridge and the Pulse button.

1.3 Definitions, Acronyms and Abbreviations

| Term | Definition |
|--------------|---|
| Validation | Assure that the software system (function and performance) conform to user needs and intended use(s). Confirmation by examination and provision of objective evidence that the particular requirements for a specific intended use can be consistently fulfilled |
| Verification | A systematic procedure of review, analysis, and testing, employed throughout the software development cycle to assure that the software design outputs meet the design inputs (all design requirements). Confirmation by examination and provision of objective evidence that specified requirements have been fulfilled. |
| Walkthrough | A type of design review in which several participants review a work product. The objective of the meeting is to discover faults in the product |
| SRS | System requirements specification |
| SwRS | Software requirements specification |
| VVP | Verification and validation plan |
| DHF | Design History File |

1.4 References

| ID | Reference |
|---------------|--|
| IEEE-Std-1012 | Std 1012-1998 IEEE Standard for Software Verification and Validation. IEEE Standards Collection: Software Engineering. The Institute of Electrical and Electronics Engineers, Inc., 1998 |

| Document | Location |
|---|----------|
| Software Requirements Specification documents | DHF |

1.5 Overview of Document

This document describes the testing procedures and requirements that will be implemented for our project. The first section goes over review, walkthroughs, inspections, and audits and gives a schedule for them the next section talks about how each component of the system will be tested. Then is a section of how the system will be tested as a whole. The following sections is about how to tract defects. After that we have a section that maps the requirements in the SRS to their description in the SwDD. We next have a section that gives the test requirements for each of the components. And finally we have the test requirements for the release of the project.

2 Reviews; walkthroughs; inspections; and audits

2.1 Schedule of reviews, etc.

Refer to project development plan document.

2.2 Procedure for reviews, etc.

The following procedure has been used to review the documents that we have completed so far. First ideas are contributes as to what the document should contain, then one person from the group puts all idea into the document (refer to design verification SW Life Cycle).

3 Component Test Plans and Procedures

3.1 Component and Integration Test strategy overview

Each software component should have the capability of being tested "individually" within the development environment. As a minimum, test cases should be developed to cover "black box" or functional requirements and component interfaces. "White box" testing to measure percentage code path coverage is also desirable.

Each component based on safety requirements shall have test case(s), and a test procedure(s).

Integration testing will be performed in 2 cycles:

The first cycle focus on the interface requirements for the various design entity that make up Flash N' Go Firmware. Special tools and concepts may be developed in order to help facilitate integration testing. Emulator was used to test the functionality of each design entity and the following test attributes were applied:

- Critical functions and Tasks (Normal Task; Interrupt Tasks)
- Boundary conditions.
- (b) (4)
- Incorrect and unexpected sequences and timing.
- Reaction of the firmware to system faults and failures.
- (b) (4)
- The Man Machine Interface.

The second cycle is integration with the real hardware and formal test procedures was followed and cover all the functional behavior of the system in normal progressive order and out of normal progressive order.

3.2 Component Test description

Below we summarize the modules, functionality, testing strategy, and test sets that we will be employing:

| Design Entities | Functionality/Tasks | Testing Strategy | Test sets |
|---|--|-------------------------------------|---|
| 1. App entity - Main | - Init task that initializes the MCU pinout, variables (b) (4) - ERROR handler that is used to respond to errors. - IDLE task that is used as - Capacitor Voltage handler that is used to manage the capacitor voltage. (b) (4) | Unit testing | (b) (4) MPLAB ICD 2 emulator and developer environment |
| 2. App entity - ISR | ISR task is called every 50ms and is used for timing purpose | Code Verification | Developer constructed debugger environment |
| 3. RFID entity – I2C:communication protocol | Receive; Send; Transmit; Stop; | Communication protocol verification | I2C verification tools |

(b) (4)

| Design Entities | Functionality/Tasks | Testing Strategy | Test sets |
|--------------------------------|---|------------------|---|
| (b) (4) | ;Read ; Write ; Deactivate | | MPLAB ICD 2 emulator and developer environment |
| 6. Capacitor Voltage entity | - Measure capacitor voltage - Charge Capacitors | Unit testing | MPLAB ICD 2 emulator and developer environment |
| 7. Code Self Check entity | Firmware code memory verification | Unit testing | Developer constructed debugger environment |
| 8. Misc entity | Create short delay time | Unit testing | Developer constructed debugger environment |
| 9. MMI entity | LEDs Indication. Buttons polling, New system Level Limitation | Unit testing | MPLAB ICD 2 emulator and developer environment |
| 10. Error entity | Check Errors | Unit testing | MPLAB ICD 2 emulator and developer environment |
| (b) (4) | | | |

3.3 Integration test description

3.3.1 Initialization

| | |
|--------------------------------|---|
| Test case group identification | Flow of initialization in the init state |
| Features to be tested | Reset and initialization |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.1.1 Initialization test case #1 – Normal flow

| | |
|-----------------------------|---|
| Test case identifier | System Initialization with correct firmware. |
| Input | Firmware is ready for integration testing |
| Expected output | Software initialized and pulse emission is enabled |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | The test case steps should follow the software detailed design sec 4.2.3 |
| References | Software System test plan and results document |

3.3.1.2 Initialization test case #2 –out of Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Initialization with corrupted firmware or other hardware and software errors. |
| Input | Firmware is ready for integration testing |
| Expected output | Software initialized and Error LED displayed with pulse emission is disabled. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | Make use of special code with wrong check sum of the firmware. |
| Precedence and dependencies | NA |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |

3.3.2 Set Energy Level

| | |
|--------------------------------|---|
| Test case group identification | Flow of set energy level |
| Features to be tested | Enables the capacitors charging and checks when the charged capacitors have reached to the required level. |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.2.1 Set Energy Level test case #1 – Normal flow

| | |
|-----------------------------|---|
| Test case identifier | System set energy level (General actions / Treatment) |
| Input | Firmware is ready for integration testing |
| Expected output | The capacitors charged to the required level |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Initialization test case #1 – Normal flow completed successfully. |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |

3.3.2.2 Set Energy Level test case #2 – out of Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Error handling of set energy level functionality(e.g. hardware will prevent the software from exceeding a maximal allowed charge voltage) |
| Input | Firmware is ready for integration testing |
| Expected output | All errors described in the software design document are performed. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Set Energy Level test case #1 – Normal flow completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |

3.3.3 Pulse Handler

| | |
|--------------------------------|---|
| Test case group identification | Pulse Handler functionality |
| Features to be tested | Emission of pulse upon successful fulfillment of required conditions |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.3.1 Pulse Handler test case #1 – Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Verify generate a pulse |
| Input | Firmware is ready for integration testing |
| Expected output | Power meter display actual pulse emitted value. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Set Energy Level test case #1 – Normal flow completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |

3.3.3.2 Pulse Handler test case #2 – out of Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Error handling of generates a Pulse. |
| Input | Firmware is ready for integration testing. |
| Expected output | No measurement value on Power meter display 9(Pulse is not generated). |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Set Energy Level test case #1 – Normal flow completed successfully and Pulse Handler test case #1 completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |

3.3.4 Monitor the capacitors voltage

| | |
|--------------------------------|---|
| Test case group identification | Monitor the capacitors voltage functionality |
| Features to be tested | Charge or discharge according to the Energy Level |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.4.1 Monitor the capacitors voltage test case #1 – Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Capacitors charging |
| Input | Firmware is ready for integration testing |
| Expected output | Capacitors charged and LED indicators turn on. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Set Energy Level test case #1 – Normal flow completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

3.3.4.2 Monitor the capacitors voltage test case #2 – out of Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Error handling of capacitors charging |
| Input | Firmware is ready for integration testing |
| Expected output | Capacitors can not be charged and Error LED indicators turn on. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Monitor the capacitors voltage test case #1 – Normal flow completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

(b) (4)

| | |
|--------------------------------|---|
| Test case group identification | (b) (4) |
| Features to be tested | (b) (4) |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.5.1

(b) (4)

| | |
|-----------------------------|---|
| Test case identifier | (b) (4) |
| Input | Firmware is ready for integration testing |
| Expected output | Capacitors charged and LED indicators turn on. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Turn on system, Use color Panel I (Go) |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

3.3.5.2

(b) (4)

| | |
|-----------------------------|---|
| Test case identifier | (b) (4) |
| Input | Firmware is ready for integration testing |
| Expected output | Capacitors can not be charged and Error LED indicators turn on. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | Turn on system, Use color Panel II (No go) |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

3.3.6 Hardware Check

| | |
|--------------------------------|---|
| Test case group identification | Hardware and software recovery |
| Features to be tested | Monitor hardware for abnormal or fault condition |
| Testing Approach | Behavioral and functional testing |
| Pass/Fail criteria | Refer to software system test plan for detailed criteria that the product must satisfy in order to pass this test case group. |

3.3.6.1 Hardware Check test case #1 – Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Verify system recovery |
| Input | Firmware is ready for integration testing |
| Expected output | Error LED is turned on and the firmware is stuck in an infinite loop. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | NA |
| Precedence and dependencies | All normal functional test cases are completed successfully |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

3.3.6.2 Hardware Check test case #2 – out of Normal flow

| | |
|-----------------------------|---|
| Test case identifier | Verify watchdog functionality within out of normal flow of software |
| Input | Special test firmware is ready for integration testing |
| Expected output | The system will repeatedly sound a beep for ready, delay then reset. |
| Environment | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring equipment. |
| Special Procedures | Add an infinite loop after system ready sounds |
| Precedence and dependencies | NA |
| References | Refer to software system test plan for detailed test environment requirements, including software, hardware and monitoring |

4 System Test Plans and Procedures

4.1 System test Strategy Overview

System testing will be based on system level scenarios and system level requirements such as safety, regulatory, performance and capacity..

4.1.1 Scenarios

Each system level scenario should have an associated system level test case. This provides traceability and validates/verifies that the scenario has been implemented. The collection of "system" level scenarios tests could make up the majority of the tests called out in the software system test plan. Scenario test cases could be developed in parallel with Flash N' Go software.

4.1.2 Safety and Regulatory

Safety and regulatory requirements described in the Flash N' Go System SRS and Risk Analysis must have test cases that verify Flash N' Go software compliance.

4.1.3 Performance

Performance requirements such as time for capacitor to be charged, need to be verified for all life cycle phases and should be part of the core system tests.

4.1.4 Regression

Testing that is used to confirm the system functionality has not been comprised by incremental code upgrades, improvements, and bug fixes. These tests may also have a specific focus on SPR resolution. A regression test should be developed to confirm that a SPR has been resolved. This test should remain in the regression test suite throughout product life.

4.1.5 Reliability/Stress/Simultaneity

Part of the system testing should include the means necessary to measure software reliability. Stress and Simultaneity testing should be designed to provide maximum functionality coverage.

4.2 System Test Description

Refer to Software System test plan and results

5 Defect tracking plans

This defect tracking document will include severity of the error (critical, high, medium, trivial), the phase in which the error is inserted and detected (requirements, coding, design, alpha test, beta test or acceptance test) and the user responsible for following up on the error. An appropriate description of the error detected should be provided as an abstract.

as a tester the user will record the software fault detected. The user will then notify the development team that a defect has been opened. The developer will resolve the defect and give the tester a new fixed version.

After the defect is verified the initial tester will run through the same steps to ensure that the original problem is corrected. The error will be closed after ensuring the problem is fixed else reopened for revaluation. As a developer, the user will monitor the open defects. If the user comes across a defect that was opened against code that he/she is responsible for, the user will view the defects details. After proper understanding of the defect the user will correct the code and put it into the verify state.

6 Traceability from SRS to SDD

Refer to the traceability matrix document.

7 Test requirements cross-reference Matrix

Refer to the traceability matrix document.

8 Acceptance Test and Preparation for delivery

8.1 *Procedure by which the software product will be acceptance tested*

The Flash N' Go firmware will be tested by the team, the client, a teacher, and a group of 8th grade students.

- The team will test all the test cases defined in the software system test plan to ensure they work as they are expected.

8.2 *Specific acceptance criteria*

Refer to the acceptance criteria of the system V&V.

8.3 Scenario by which the software product will be installed.

The system will be delivered with the firmware download to the device at the production.

Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go
Unresolved Software Anomalies

DO1009171

| | Role | Name | Signature | Date |
|---------|------|------|-----------|------|
| (b) (4) | | | | |

Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Table of Contents

| | | |
|-----|--|---|
| 1 | Introduction | 4 |
| 1.1 | Purpose | 4 |
| 1.2 | Scope | 4 |
| 1.3 | Definitions | 4 |
| 1.4 | References | 4 |
| 2 | Unresolved Anomalies – Software Ver. 1.0.0 | 5 |

1 Introduction

1.1 Purpose

The purpose of this document is to present a list of all unresolved software anomalies.

1.2 Scope

This document is used to specify software anomalies as specified in the "Guidance for the Content of Pre-market Submissions for Software Contained in Medical Devices"

1.3 Definitions

The following table defines terms and definitions used:

| Term | Definition |
|------|-------------------------|
| SPR | Software Problem Report |
| | |

1.4 References

The following table lists all documents which supply inputs to this document:

| Document | Doc No. | Revision | Location |
|--|-----------|----------|----------|
| Software System Verification and Validation Plan | DO1009131 | 1 | DHF |
| | | | |

2 Unresolved Anomalies – Software Ver. 1.0.0

| Bug No. | Problem | Impact on device performance | Timeframes for correcting the problem (where appropriate). |
|---------|---------|------------------------------|--|
| | | | |
| | | | |
| | | | |
| | | | |

Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go Software Revision History

DO1009161

| Authorization | | | | |
|---------------|------|------|-----------|------|
| | Role | Name | Signature | Date |

(b) (4)




Skinovations

Home Skinovations Ltd.
Apolo Building,
Soltam Industrial Zone.
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |



Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

1. General

A software history revision log is provided below. Clearance is being thought for version [1.2.0], which is the most recent version of the software.

2. Revision Level History

| Revision number | Creation date | Notes |
|-----------------|---------------|-------|
| (b) (4) | | |
| | | |
| | | |
| | | |
| | | |

166

1. Traceability Analysis

The following table gives information on the traceability between software requirements and detailed design and the corresponding verification and validation test that addresses that requirement.

| Software Requirement Identifier | Software Design Identifier | V&V Test |
|---------------------------------|-------------------------------|--|
| 2.2 (a) | 2.1.5 | 9.1.1 |
| 2.2 (b) | 2.1.1 | 9.2.1;9.2.2 (safety related) |
| 2.2 (c) | 2.1.5 | 9.4.2 |
| 2.2 (d) | 2.1.3 | 9.1.1;9.2.3 |
| 2.2 (e) | 2.1.2 | 9.1.4; 9.1.6; 9.2.3(4); 9.2.4(6); 9.3.3 |
| 2.2 (f) | 2.1.7 | 9.1.2;9.1.3 |
| 2.2 (g) | 2.1.4 | 9.3.1; 9.3.2; 9.3.3; 9.4.3 |
| 3.1.1.1 | HW Design | 9.1.1 (1) |
| 3.1.1.2 | HW Design | 9.1.1 (1) |
| 3.1.1.3 | 4.6.1.2 | 9.1.1(2) |
| 3.1.1.4 | 4.6.2.4 | 9.4.2(3) |
| 3.1.1.5 | 4.6.2.1;4.6.2.2 | 9.1.1(3) |
| 3.1.1.6 | 4.6.1.1 | 9.1.1(3) |
| 3.1.1.7 | 4.6.2.3;4.6.3 | 9.1.1(3) |
| 3.1.1.8 | 4.6.1.5 | 9.2.4(8) |
| 3.1.1.9 | 4.6.1.6 | 9.2.4(10) |
| 3.1.1.10 | 4.6.1.3 | 9.2.2(1) |
| 3.1.2.1 | 2.2.1 | 9.3.3(1-2) |
| 3.1.2.2 | 2.2.2 | 9.2.3(4) |
| 3.1.2.3 | 2.2.3 | 9.2.3;9.2.4 |
| 3.1.2.4 | 2.2.4 | 9.1.1(3) |
| 3.1.2.5 | 2.2.5 | 9.1.1(3) |
| 3.1.2.6 | 2.2.6 | 9.3.1 |
| 3.1.2.7 | 2.2.7 | 9.2.7, 9.3.7 |
| 3.1.2.8 | 2.2.8 | 9.2.7, 9.3.7 |
| 3.1.3 | 2.2 | 9.1.5 (1-4) |
| 3.1.4.1 | 5.2.2.1 | 9.1.6 (1-2) |
| 3.2.1 | 4.2 | 9.2.1(1-2);9.2.2(1-2) |
| 3.2.1.3.1 | 4.2.1 | 9.2.3 |
| 3.2.1.3.2 | 5.2.1.1 | 9.2.1 |
| 3.2.1.3.3 | 4.2.3 | 9.2.1 |
| 3.2.1.3.4 | 4.7.2;4.7.3;4.7.4;4.7.5;4.7.6 | 9.1.1 (2) |
| 3.2.1.3.5 | 5.2.2.3 | 9.1.1 (4); 9.1.4 (2) |
| 3.2.1.3.6 | 4.6.1.1;5.2.6.1 | 9.2.3 |
| 3.2.1.3.7 | 5.2.3 | 9.2.3(4) |
| 3.2.1.3.8 | 4.6.1.4;4.6.3 | 9.1.1 (3; step 3) |
| 3.2.1.3.9 | 4.6.1.2;5.2.6.1 | 9.1.1 (6) |
| 3.2.2.3.1 | 5.2.6.2 | 9.1.1(3;step1 and 2) |

Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

[illegible]

Home Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go
Traceability Analysis

DO1009151

(b) (4)



Skinovations

Home Skinovations Ltd.
Apolo Building,
Soltam Industrial Zone,
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Home Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' GO

Software Detailed Design Spcification

DO1008811

| | Role | Name | Signature | Date |
|---------|------|------|-----------|------|
| (b) (4) | | | | |

Home Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Table of Contents

| | | |
|-----|--|----|
| 1. | Introduction | 4 |
| 1.1 | Purpose | 4 |
| 1.2 | Scope | 4 |
| 1.3 | Definitions and acronyms | 4 |
| 1.4 | References | 4 |
| 1.5 | Design Conventions | 4 |
| 2 | Design Overview | 5 |
| 2.1 | Software Context Diagram | 5 |
| 2.2 | Hardware Software Interfaces | 6 |
| 2.3 | Theory of Operation | 10 |
| 3 | Design Rationale | 10 |
| 4 | Architecture | 10 |
| 4.1 | Operating System | 10 |
| 4.2 | Resets and Initialization | 10 |
| 4.3 | Interrupt service routine | 11 |
| 4.4 | Tasks | 11 |
| 4.5 | State machine Description | 12 |
| 4.6 | Presentation view | 13 |
| 4.7 | Error handling | 14 |
| 5 | Detailed Design Description | 15 |
| 5.1 | External Interface Data Elements | 15 |
| 5.2 | Design entities | 15 |

2.1.5 MMI

The MMI is used to manage turn ON/OFF LEDs check the keys and turn ON/OFF the buzzer, Limit new system levels

2.1.6 Light Pulse handler

The light pulse handler is used to set the Light enable signal if conditions are met.


2.1.7 Timer ISR

The Timer ISR is used to measure time, handle timer variables, and check for time outs

(b) (4)


2.2 Hardware Software Interfaces

(b) (4)




2.2.2 Capacitor interfacing

(b) (4)




2.2.3 Indication lights


(b) (4)




2.2.4 KEY polling




2.2.5 Buzzer interface



2.2.6 Fan status




2.2.7 Color detection (b) (4)



2.2.8 Color detection (b) (4)





2.3 Theory of Operation

The software will first initialize the hardware. Then initialize the software flags and variables. The software then will start a normal operation of charging capacitors waiting for a trigger button and generating a light pulse. During this time the software will check for any problems and will halt if any problem occurs.

3 Design Rationale

The software structure is based on the need for a simple, small sized customized to the system hardware design.

4 Architecture

4.1 Operating System

The operating system is built from the ground up using a simple design. The software operating system consists of an executive loop that invokes a series of functions in a round-robin fashion along with interrupts that provide pre-emptive processing. The functions invoked by the executive loop are called tasks. Each task provides a unique set of functionality. Some of the tasks that make up the software include the, high-voltage ramping task, and count rate protection task, among others.

4.2 Resets and Initialization

4.2.1 Power-On Reset

Upon POR the SW will wait for (b) (4), then will start the initialization process.

4.2.2 Watch-Dog

A WD reset occurs if the software fails to clear the WD timeout counter.


4.2.3 Complete initialization sequence

A complete initialization sequence is performed (b) (4)

(b) (4)

- Configures MCU pins
- Initialize the I/Os.
- Setup the LED indications.
- Check the code corruption by comparing the checksum of the ROM; if the checksum failed error LED is turned on and the software code is disabled.
- Initialize flags.

(b) (4)



4.3 Interrupt service routine

Interrupt service routines (b) (4) pre-emptive processing in the software. Since the operating system consists of a round-robin loop that executes tasks in a pre-determined order, only interrupts can provide pre-emptive processing.

This section includes the WD timer and the Timer interrupts.

4.3.1 WD timer interrupt

The WD timer will generate a non mask-able interrupt that will restart the MCU the same as POR.

4.3.2 Timer interrupt

The timer interrupt is set to occur every 50ms, timers are decremented down to zero and, Keys are sampled and checked.

4.4 Tasks

Tasks are invoked in round-robin fashion by the executive loop. The executive loop does manage the switching in and out of tasks, the tasks, therefore, must be capable of retaining their own state information. To ensure timely progress through the executive loop, tasks must execute quickly and return control to the executive loop. To accomplish this, each task runs to completion.


4.4.1 Normal Tasks

The normal tasks of the software include interfacing with the (b) (4), capacitor voltage management, MMI, (b) (4) error handling

4.4.2 Interrupt Tasks

Interrupt tasks include time based task including timing, key polling

4.5 State machine Description



4.5.1 Init State

In this state the Firmware initiates MCU registers, verifies Firmware code (b) (4), initiates variables and flags (b) (4)

4.5.2 Idle state

In this state the firmware checks the keys (b) (4)

4.5.3 Charge state

In this state the Firmware enables the capacitor charging (b) (4)

(b) (4)

4.5.4 Pulse state

In this state the Firmware generates a light pulse.

4.5.5 Error state

In this state the Firmware will handle errors were an error LED is turned on and the Firmware is stuck in an infinite loop.

4.6 Presentation view.

4.6.1 LEDs

4.6.1.1 Level Indication

There are 5 Level Indication LEDs that indicate the flash level. During Error the Level indication shows the error code. Upon Hand piece push button press Level indication will indicate Miss-matched color

4.6.1.2 Ready Indication

The Ready LED indicates that the unit is ready for next pulse.

4.6.1.3 Error Indication

The Error LED indicates that a fault occurred and the unit is in a fault protection.

4.6.1.4 Buzzer Enabled Indicator

The Buzzer LED indicates that the buzzer is enabled and will sound when required.

(b) (4)

4.6.2 Keys

4.6.2.1 Energy Plus

The Energy Plus key is used to increase the light pulse level. (b) (4)

(b) (4)

4.6.2.2 Energy Minus

The Energy Minus key is used decrease the light pulse level. (b) (4)

(b) (4)

4.6.2.3 Buzzer Enable/Disable

The Buzzer key, switches the buzzer mode between On and Off.

4.6.2.4 Hand piece button

The hand piece button is used to fire a single pulse of light.

4.6.3 Buzzer

The Buzzer is used to indicate to the user that the next pulse can be fired.

4.7 Error handling

Upon any error turn the system OFF, wait for a few minutes, turn the system ON. If the error re-occurs turn the system OFF, and contact service.

4.7.1 Corrupt Firmware code

Corrupt firm code will be handled by setting up the error LED and going into an endless loop.

(b) (4)



5 Detailed Design Description

5.1 External Interface Data Elements

| Interface | Type | Direction | Data Elements Description |
|---------------------|---------|----------------|---------------------------|
| KEYS | Digital | Input | Key event handler |
| Hand Piece PB | Digital | Input | Fire light pulse |
| LEDS | Digital | Outputs | User indications |
| Buzzer | Digital | Output | |
| (b) (4) | Digital | Input / Output | Cartridge interface |
| IPL | Digital | Output | Generate light pulse |
| Fan OK | Digital | Input | Fan status |
| (b) (4) | | | Temp temperature |
| Color Detection LED | Digital | Output | Color LED |
| (b) (4) | | | |

5.2 Design entities

5.2.1 APP entity

5.2.1.1 Main

The Main includes several tasks

- Init task that initializes the MCU pinout variables.

(b) (4)

- ERROR handler that is used to respond to errors.
- IDLE task that is used as
- Capacitor Voltage handler that is used to manage the capacitor voltage.
- Pulse task that is used enable a light pulse when pulse conditions are met

5.2.1.2 Isr

ISR task is called every 50ms and is used for timing purpose

(b) (4)

(b) (4)



5.2.3 Capacitor Voltage

5.2.3.1 Mesure_capacitor_voltage

The mesure_capacitor_voltage is used to sample the voltage using the MCUs built in ADC and calculate the voltage

5.2.3.2 Charge_capacitors

The charge_capacitors is used to enable the capacitor charging

5.2.4 Code Self Check

5.2.4.1 Check_program_memory

The check_program_memory is used to validate the Firmware code

5.2.5 Misc

5.2.5.1 Short_delay

The short_delay is used to create a short delay time

5.2.6 MMI

5.2.6.1 Mode_status

The mode_status is used to update the mode LEDs to the current mode.

5.2.6.2 Buttons_status

The buttons_status is used to poll the system keys and respond to the valid combinations.

5.2.6.3 Level_limitation

Level limitation is used to limit the energy level 4 and 5 for a new system, only level 1 to 3 are allowed for 250 first pulses

5.2.6.4 Turning on/off Level_limitation

Level limitation can be turned on/off by pressing the (=) and (-) keys simultaneously 3 short beeps indicate limitation on (pulse counter set to 250) , 3 short beeps indicate limitation off

5.2.7 Error

5.2.7.1 Check_error_status

The check_error_status is used to handle a condition were the pulse counter reaches zero.

(b) (4)



Home Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go

Software Requirements Specification

DO1009141

| Written by | Role | Name | Signature | Date |
|------------|------|------|-----------|------|
|------------|------|------|-----------|------|

(b) (4)

Home Skinovations

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | | |
|---------|------|-------------|--|--|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Table of Contents

| | |
|---|----|
| 1. INTRODUCTION | 4 |
| 1.1 PURPOSE | 4 |
| 1.2 SCOPE | 4 |
| 1.3 DEFINITIONS | 4 |
| 1.4 REFERENCES | 4 |
| 1.5 SRS ORGANIZATION | 5 |
| 1.6 OVERVIEW | 5 |
| 1.7 REQUIREMENT STATEMENT IDENTIFICATION CONVENTION | 6 |
| 2. OVERALL DESCRIPTION | 6 |
| 2.1 PRODUCT PERSPECTIVE | 6 |
| 2.2 PRODUCT FUNCTIONS | 7 |
| 2.3 USER CHARACTERISTICS | 8 |
| 2.4 GENERAL CONSTRAINS | 8 |
| 2.5 ASSUMPTIONS AND DEPENDENCIES | 8 |
| 3. REQUIREMENT SPECIFICATIONS | 8 |
| 3.1 INTERFACE REQUIREMENTS | 8 |
| 3.2 SYSTEM FEATURES | 11 |
| 3.3 PERFORMANCE REQUIREMENTS | 14 |
| 3.4 SAFETY REQUIREMENTS | 15 |
| 3.5 LOGICAL DATABASE REQUIREMENTS | 15 |
| 3.6 DESIGN CONSTRAINS | 16 |
| 3.7 SOFTWARE SYSTEM ATTRIBUTES | 16 |
| 3.8 TESTING REQUIREMENTS | 17 |
| 3.9 OTHER REQUIREMENTS | 17 |

1. Introduction

1.1 Purpose

The purpose of this document is to provide a consistent and complete description of the requirements for the software of Flash N' Go. The requirements will be presented using textual descriptions to explain concepts, different types of diagrams to illustrate complicated interactions, and tables to relate relevant information.

The intended audience of this document is all of the managers and engineering involving the development of the MCU firmware software. This includes, but is not limited to, software developers, project managers, quality assurance personnel, and customers.

1.2 Scope

The firmware of the Slik'n is responsible for the normal and efficient operation of all of the other components within the Flash N' Go device. The controller's main goal is essentially to handle input requests from other components and respond accordingly with generating pulses.

It is important to note that this SwRS document only pertains to the requirements of the software of the MCU. It does not include requirements for the hardware that it will be deployed on. It also does not include requirements for the other components of the Flash N' Go device. Other components are mentioned because the requirements of the controller are indeed based on the interactions between it and the other components of the system; however, the requirements of the other components, as well as the overall Flash N' Go device itself, is outside the scope of this document and is therefore not included.

1.3 Definitions

The following table defines terms and definitions used in this SRS:

| Term | Definition |
|---------------|--|
| SwRS | Software Requirements Specification |
| V&V | Verification and Validation |
| MMI | Man Machine Interface |
| CPU | Central Processing Unit |
| I/O | Input/ Output |
| Control Panel | Two operation button and LED Display |
| C | Program language for microcontrollers. |
| PS | Power Supply |
| RAM | Random Access Memory |

| Term | Definition |
|--------------------------------|-----------------------------|
| I ² C communication | I2C communications standard |
| Cartridge Memory | EEPROM |
| PB | Push Button |
| HP | Hand Piece |

1.4 References

The following table lists all documents which supply inputs to this SRS:

- 1) Flash N' Go operation Manual.
- 2) Flash N' Go System Requirements Specification
- 3) IEC 60601-1-4 General requirements for safety programmable electrical medical system.
- 4) IEEE 830-1998 standard for creating SRS documents is IEEE Recommended practice for software requirements specification by IEEE-SA computer society (IEEE-SA Standards Boards, 1998)

1.5 SRS Organization

In section 2 we will give a general overview of the functionality of Flash N' Go firmware. It is used to establish a context for both the informal requirements definition and for technical requirement specifications in the next chapter. This section consists of six subsections, as follows:

- Product perspective.
- Product functions.
- User characteristics.
- Constrains.
- Assumptions and dependencies.
- Apportioning and requirements.

In section 3 we will give a detailed overview of the functional and non-functional requirements of the Flash N' Go firmware. This section consists of six subsections, as follows:

- External interfaces.
- System features .
- Performance requirements.
- Design constrains.
- Software system attributes.

1.6 Overview

The purpose of the software in the Flash N' Go device:

- Control user interacts with the device (Power On/Off switch; Energy level settings buttons; Buzzer enable/disable and Pulse Button).
- Check the Hardware status and inform user through the indicators light.

1.7 Requirement Statement Identification Convention

The requirement statement identification convention is based on section numbering.

2. Overall description

This section describes the general factors that affect the Flash N' Go firmware software requirements. This section does not state specific requirements; instead it provides a background for those requirements, which are defined in detail in Section 3.

The Flash N' Go firmware will allow an operator to view system status and generate pulses.

The major functions of the Flash N' Go firmware are: MMI, Setup and Initialization; Error Handler; Pulse Handler; Capacitor voltage handler; RFID Handler; ISR Timer;

2.1 Product perspective

The Flash N' Go firmware directly interacts with the following other components of the entire device:

Indicators Light:

Desired parameters.


Pulse Button:

Desired parameters.

Energy Level Settings Buttons:

Desired parameters.

Context Diagram



2.2 Product Functions

The firmware is responsible for the following functions:

- a) Manage Man Machine Interface.
- b) Manage setup and initialization.
- c) Pulse Handler.
- d) Capacitor voltage handler.

(b) (4)

- f) ISR Timer.
- g) Error handler.

(b) (4)

2.3 User characteristics

Users are required to be familiar with the Operation Manual and the use of the Flash N' Go device.

The users of the firmware are essentially the other components of the device that interact with it. The integral characteristic that these other components must have is being able to process and react to the request sent to it by the controller in a timely fashion.

2.4 General constrains

The following constraints were factors in the development of the detailed requirements for the Flash N' Go firmware:

- a) Regulatory policies – according to the software life cycle procedure.
- b) Hardware limitations – The Flash N' Go firmware will reside on hardware MCU PIC16F887.

2.5 Assumptions and dependencies

The software is strongly hardware dependent;

3. Requirement Specifications

This section of the SwRS contains the Flash N' Go firmware requirements to a level of detail sufficient to enable designers to design the system, and acceptance testers to test that the device satisfies the requirements.

These software requirements will be revised by the system engineer to reflect the final approved Flash N' Go architecture design. This will be done before high-level software design phase. Every stated requirement is externally perceivable by users, operators, or other external systems.

This section describes each input into the software, every output from the software, and all functions performed by the software.

3.1 Interface Requirements

3.1.1 User Interfaces

The MMI will be implemented to allow the user to communicate with the device:

3.1.1.1 Power On/Off switch

The power on button shall activate the HW and the software will be activated as a result.

3.1.1.2 Power Indicator Light

The power indicator lights shall turn on after pressing the Power On/Off switch.

3.1.1.3 Ready Indicator Light

3.1.1.3.1 The ready indicator light shall turn on after pressing the Power On/Off switch.

3.1.1.3.2 The Ready Indicator Light shall turn on each time after a pulse has been fired and a time out had passed.

3.1.1.4 Pulse Button

The user shall handle pulses by pressing the Pulse Button on the Applicator unit.

3.1.1.5 +/- Energy Level Settings Buttons

The user shall select the required energy level using the Energy Level Settings Buttons. Pressing both (+) and (-) simultaneously shall toggle level limitation on / off.

3.1.1.6 Energy Level Indicator Light

The user shall observe the control panel for the proper energy level indicator lights on. The audio on/off switch used to enable/disable audible sound during normal operation.

3.1.1.7 Audio On/Off switch

The user shall hear a sound if (Audio is on) after the Ready indicator Light turned on.

(b) (4)



3.1.1.10 System Warning Indicator Light

The software shall turn on the System Warning Indicator once an error occurred.

3.1.2 Hardware Interface

The detailed Flash N' Go hardware to software interfaces are described below to illustrate the key interfaces to the Flash N' Go software.

(b) (4)



3.1.2.2 Capacitor Interfacing

Refer to the Software Detailed Design Specification.

3.1.2.3 Indication Lights

Refer to the Software Detailed Design Specification.

3.1.2.4 Key Polling

Refer to the Software Detailed Design Specification.

3.1.2.5 Buzzer Interface

Refer to the Software Detailed Design Specification.

3.1.2.6 Fan status

(b) (4)



3.1.3 Software Interfaces

Flash N' Go software is using:

Microchip: MPLAB ICD 2 emulator and developer environment.

The interface requirements for this product are defined by the MPLAB ICD 2 Reference Manual.

Flash N' Go software is planning on using the following capabilities of MPLAB ICD 2 emulator:

- During development using the Debug.
- Down load function.

(b) (4)

3.2 System Features

3.2.1 Initialization

3.2.1.1 Purpose

The Flash N' Go shall initialize the software and Hardware for treatment start.

3.2.1.2 Stimulus/Response sequence

Press the power On/Off switch; the initialization process will be activated; the Ready Indicator Light will turn on, accompanied by a sound to indicate that your device is ready for its first pulse.

(b) (4)

3.2.1.3 Associated functional requirements

3.2.1.3.1 The On/Off Power Switch shall be used to turn the device on/off.

3.2.1.3.2 The MCU internal registers and pin-outs shall be initialized.

3.2.1.3.3 The Software variables shall be initialized.

3.2.1.3.4 The software shall check the hardware.

3.2.1.3.5 The software shall read data from the RFID tag.

3.2.1.3.6 The level LED's shall indicate the current level.

3.2.1.3.7 The software shall activate the charger according to the required level (b) (4)

3.2.1.3.8 The Buzzer shall sound to indicate your device is ready for its first use (b) (4)

3.2.1.3.9 The Ready LED shall be turned on to indicate the device is ready for its first pulse.

3.2.2 Set Energy Level

3.2.2.1 Purpose

The device shall have an energy setting buttons for (b) (4)

3.2.2.2 Stimulus/Response sequence

Press the +/- Energy Level Setting buttons to switch between treatment energy levels.

3.2.2.3 Associated functional requirements

3.2.2.3.1 The Energy Level Setting Buttons shall be used to select one of the (b) (4) treatment energy

3.2.2.3.2 The Energy Level Indicator shall be turned on after selecting the required energy level.

3.2.2.3.3 The Energy setting buttons shall be used to turn off the Ready Indicator Light.

3.2.2.3.4 (b) (4)

3.2.2.3.5

3.2.3 Pulse Handler

3.2.3.1 Purpose

The device shall have a pulse handle to enable a pulse emission.

3.2.3.2 Stimulus/Response sequence

Press the Pulse button; Light pulse emitted only after software check required conditions.

3.2.3.3 *Associated functional requirements*

3.2.3.3.1 The pulse handler shall emit pulse only upon the completion of the following conditions:

3.2.3.3.1.1 *Pulse Counter is greater than zero.*

3.2.3.3.1.2 *Pulse Button and Lamp cartridge should be pressed.*

3.2.3.3.1.3 *Capacitor voltage has reached the required level.*

3.2.3.3.1.4 *Minimal time between two pulses should be TBD seconds.*

(b) (4)

3.2.3.3.2 The Pulse handler shall decrease the Pulse Counter by one.

3.2.4 *Monitoring the Capacitors voltage*

3.2.4.1 *Purpose*

This feature shall monitor the capacitors voltage.

3.2.4.2 *Stimulus/Response sequence*

If the software get event from the Energy Level Settings Buttons the capacitors will be charged or discharged according to the current Energy Level and the new selected energy level and according to the calibration factor.

3.2.4.3 *Associated functional requirements*

3.2.4.3.1 Monitoring the Capacitors Voltage shall be used to measure the capacitor voltage.

3.2.4.3.2 The Energy setting buttons shall be used as an event to Charge the Capacitors.

3.2.4.3.3 The capacitor voltage will vary according to the calibration factor.

3.2.5 Hardware Check

3.2.5.1 Purpose

The Flash N' Go software shall monitor the hardware for abnormal or fault condition

3.2.5.2 Stimulus/Response sequence

Upon any fault or miss function of the hardware the software shall handle this as an error condition refer to the Error Detection and Handling section in this document

3.2.5.3 Associated functional requirements

3.2.5.3.1 Fan: The software shall monitor the Fan for miss function

3.2.5.3.2 Lamp too hot: The lamp temperature shall be checked to protect against over heated lamp

(b) (4)

3.3 Performance Requirements

3.3.1 The software has the capability of responding to user inputs within a reasonable period of time.

3.3.2 Detects abnormalities and corrects them at least every 50 millisecond (The Timer interrupt is set to occur every 50ms), for example, monitoring of Capacitor voltage, Keypad.

3.3.3 Detects critical faults once per second, (b) (4) and stop the system if necessary.

3.3.4 Within 5 seconds (b) (4) capacitor should be charged up to required level.

3.4 Safety Requirements

- 3.4.1 The software stop pulse emitted if the Applicator Tip is disconnected.
- 3.4.2 The corrupted firmware will be stopped by watchdog.
- 3.4.3 The software shall check the hardware and if error occurred; the system shall be shutdown.
- 3.4.4 The software shall perform CRC and checksum for communication error.

3.5 Logical Database Requirements

NA

3.6 Design Constrains

3.6.1 Capacitors charge fault protection

A maximal time out of 5 seconds used in protection against fault in the capacitor charging

3.6.2 Charge level

Hardware will prevent the software from exceeding a maximal allowed charge voltage.

(b) (4)



3.7 Software system attributes

3.7.1 Reliability

3.7.1.1 None the less we expect our system shall experience no severity 1 or 2.

3.7.2 Availability

3.7.2.1 In case of device failure, refer to operation manual section troubleshooting.

(b) (4)



3.7.4 Error Detection and Handling requirements

3.7.4.1 General

3.7.4.1.1 Firmware must be able to handle the errors produced by software code.

3.7.4.2 User Error

3.7.4.2.1 The firmware shall indicate system error by turn on the Error Indicator Light.

3.7.5 Maintainability

Refer to operation manual.

3.7.6 Portability

NA

3.8 Testing Requirements

All files defined in the firmware shall be subject to the testing methodology defined in the V&V plan.

3.9 Other Requirements

3.9.1 Fault Detection and Isolation

Fault detection and isolation requirements specified by the MCU are met by the Diagnostic, Memory, Debugging features described in the "Flash N' Go Software Detailed Design" document.

3.9.2 Self Test

The firmware software shall perform self-test operations as an integral part of its normal activities. Refer to Software detailed design document for further information.

3.9.3 Redundancy Management and Hardware Re-configuration

MCU redundancy requirements are not handled by the Flash N' Go Software, but are instead covered by independently hardware subsections.

3.9.4 Watchdog Reset and crash recovery

If the MCU re-boots due to a fatal error or watchdog reset, the MCU start-up software shall not install any patches. The software will idle.

3.9.5 Bad Parameter Block Corruption Handling.

In the event that the firmware receives a parameter block whose computed integrity check code (CRC or Checksum) does not match the one contained within the block, the software will indicate the error in the LED and goes to endless loops.


HomeSkinovations

Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
www.homeskinovations.com

Appendix 2
From 3674



DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug Administration

Certification of Compliance, under 42 U.S.C. § 282(j)(5)(B), with
Requirements of ClinicalTrials.gov Data Bank (42 U.S.C. § 282(j))

(For submission with an application/submission, including amendments, supplements, and resubmissions, under §§ 505, 515, 520(m), or 510(k) of the Federal Food, Drug, and Cosmetic Act or § 351 of the Public Health Service Act.)

SPONSOR / APPLICANT / SUBMITTER INFORMATION

| | |
|---|--|
| 1. NAME OF SPONSOR/APPLICANT/SUBMITTER Home Skinnovations Ltd. | 2. DATE OF THE APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES (b) (4) |
| 3. ADDRESS (Number, Street, State, and ZIP Code) Yokneam Industrial area P.O.B 533 20692 ISRAEL | 4. TELEPHONE AND FAX NUMBER (Include Area Code) (Tel.) +972 547 404 616 (Fax) +972 4 909 7471 |

PRODUCT INFORMATION

5. FOR DRUGS/BIOLOGICS: Include Any/All Available Established, Proprietary and/or Chemical/Biochemical/Blood/Cellular/Gene Therapy Product Name(s)
FOR DEVICES: Include Any/All Common or Usual Name(s), Classification, Trade or Proprietary or Model Name(s) and/or Model Number(s)
(Attach extra pages as necessary)

~~Flash N Go~~
Flash N Go

APPLICATION / SUBMISSION INFORMATION

| |
|---|
| 6. TYPE OF APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES <input type="checkbox"/> IND <input type="checkbox"/> NDA <input type="checkbox"/> ANDA <input type="checkbox"/> BLA <input type="checkbox"/> PMA <input type="checkbox"/> HDE <input checked="" type="checkbox"/> 510(k) <input type="checkbox"/> PDP <input type="checkbox"/> Other |
| 7. INCLUDE IND/NDA/ANDA/BLA/PMA/HDE/510(k)/PDP/OTHER NUMBER (If number previously assigned) K082298 |
| 8. SERIAL NUMBER ASSIGNED TO APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES |

CERTIFICATION STATEMENT / INFORMATION

9. CHECK ONLY ONE OF THE FOLLOWING BOXES (See instructions for additional information and explanation)
- ☐ A. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, do not apply because the application/submission which this certification accompanies does not reference any clinical trial.
- ☒ B. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, do not apply to any clinical trial referenced in the application/submission which this certification accompanies.
- ☐ C. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, apply to one or more of the clinical trials referenced in the application/submission which this certification accompanies and that those requirements have been met.
10. IF YOU CHECKED BOX C, IN NUMBER 9, PROVIDE THE NATIONAL CLINICAL TRIAL (NCT) NUMBER(S) FOR ANY "APPLICABLE CLINICAL TRIAL(S)," UNDER 42 U.S.C. § 282(j)(1)(A)(i), SECTION 402(j)(1)(A)(i) OF THE PUBLIC HEALTH SERVICE ACT, REFERENCED IN THE APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES (Attach extra pages as necessary)
NCT Number(s):

The undersigned declares, to the best of her/his knowledge, that this is an accurate, true, and complete submission of information. I understand that the failure to submit the certification required by 42 U.S.C. § 282(j)(5)(B), section 402(j)(5)(B) of the Public Health Service Act, and the knowing submission of a false certification under such section are prohibited acts under 21 U.S.C. § 331, section 301 of the Federal Food, Drug, and Cosmetic Act.
Warning: A willfully and knowingly false statement is a criminal offense, U.S. Code, title 18, section 1001.

| | | |
|--|---|---|
| 11. SIGNATURE OF SPONSOR/APPLICANT/SUBMITTER OR AN AUTHORIZED REPRESENTATIVE (Sign) Amir Waldman | 12. NAME AND TITLE OF THE PERSON WHO SIGNED IN NO. 11 (Name) AMIR WALDMAN (Title) VP Regulation | |
| 13. ADDRESS (Number, Street, State, and ZIP Code) (of person identified in No. 11 and 12) Yokneam Industrial area POB 533 20692 ISRAEL | 14. TELEPHONE AND FAX NUMBER (Include Area Code) (Tel.) +972 547 404 616 (Fax) +972 4 909 7471 | 15. DATE OF CERTIFICATION 10/2/08 |

HomeSkinovations

Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
www.homeskinovations.com

Appendix 3 Revised user manual



Flash N' Go User manual

Copyright © 2008 Home Skinovations Ltd. All rights reserved.

Print date: February 2008

Home Skinovations Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd. is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd.

Data is subject to change without notification.

Home Skinovations Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd.

Specifications are subject to change without notice.

Home Skinovations Ltd.
Apolo Building, Shaar Yokneam
POB 533
Yokneam 20692
ISRAEL

Website: www.silkn.com
e-mail: info@silkn.com

Table of Contents

| | |
|---|-------------------------------------|
| 1. Before You Start | <i>Error! Bookmark not defined.</i> |
| 1.1. What is Flash N' Go™? | <i>Error! Bookmark not defined.</i> |
| 1.2. Flash N' Go™ intended use | <i>Error! Bookmark not defined.</i> |
| 2. Safety with Flash N' Go™ | 5 |
| 2.1. With Flash N' Go™ Safety Comes First | 5 |
| 2.2. Is Flash N' Go™ for you? | 5 |
| 2.3. When to avoid using Flash N' Go™? | 7 |
| 2.4. Precautions -- How to use Flash N' Go™ Safely | 8 |
| 2.5. Reducing the risk of injury | 9 |
| 2.6. Possible Side Effects using Flash N' Go™ | 10 |
| 3. Long Term Hair Removal the Flash N' Go™ Way | 13 |
| 3.1. How does light remove hair? | 13 |
| 3.2. How does the hair growth cycle impact light-based hair removal? | 13 |
| 3.3. What is different about Home Pulsed Light™ technology in Flash N' Go™? | <i>Error! Bookmark not defined.</i> |
| 3.4. Plan your Flash N' Go™ hair removal for best results | 14 |
| 4. First Steps with Flash N' Go™ | 16 |
| 4.1. Flash N' Go™ device description | <i>Error! Bookmark not defined.</i> |
| 4.2. Package contents | <i>Error! Bookmark not defined.</i> |
| 4.3. Flash N' Go™ Setup | 16 |
| 4.4. Treating with Flash N' Go™ for the first time | 16 |
| 4.5. What to Expect when treating with Flash N' Go™? | 18 |
| 4.6. After treating with Flash N' Go™ | 19 |
| 5. Maintenance of Flash N' Go™ | 20 |
| 5.1. Cleaning Flash N' Go™ | 20 |
| 5.2. Replacing the Lamp Cartridge | 20 |
| 6. Troubleshooting | 22 |
| 6.1. "My Flash N' Go™ does not start." | 22 |
| 6.2. "A light pulse is not emitted when I press the PULSE SWITCH." | 22 |
| 7. Customer Service | 23 |
| 8. Frequently Asked Questions | 24 |
| 9. Specifications | 26 |
| 10. Labeling | 27 |
| EC Declaration of Conformity | 28 |

Before using Flash N' Go™ for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures.

We recommend you re-familiarize yourself with this User Manual before each use of Flash N' Go™.



Flash N' Go™ is a powerful electrical device. As such, it should be used with special attention to safety.



Please read all warnings and safety precautions before use, and strictly follow them when using Flash N' Go™

1. What is Flash N' Go™?

Flash N' Go™ is a light-based device for long-term hair removal designed for home-use.

1.1. Flash N' Go™ device description

Flash N' Go consists of a Base Unit, APPLICATOR and Disposable Lamp Cartridges.



1. BASE UNIT
2. APPLICATOR
3. PULSE BUTTON
4. DISPOSABLE LAMP CARTRIDGE
5. SKIN COLOR SENSOR
6. CONTROL PANEL
7. POWER ON/OFF SWITCH
8. POWER INDICATOR LIGHT
9. ENERGY LEVEL SETTING BUTTONS
10. ENERGY LEVEL INDICATOR LIGHTS
11. READY INDICATOR LIGHT
12. AUDIO ON/OFF SWITCH
13. CARTRIDGE 90% USAGE INDICATOR LIGHT
14. CARTRIDGE 100% USAGE INDICATOR LIGHT
15. SYSTEM WARNING INDICATOR LIGHT

Your Flash N' Go™ consists of a **BASE UNIT** and an **APPLICATOR**.

On the **BASE UNIT** you can find the Flash N' Go™ **CONTROL PANEL** including **POWER ON/OFF SWITCH** and **POWER INDICATOR LIGHT**, the **ENERGY LEVEL SETTING BUTTONS** and **ENERGY LEVEL INDICATOR LIGHTS**, the **AUDIO ON/OFF SWITCH**, **READY INDICATOR**, **SYSTEM WARNING INDICATOR LIGHT** and **CARTRIDGE 90% and 100% USAGE INDICATOR LIGHTS**.

On the Flash N' Go™ **APPLICATOR** you can find the **PULSE BUTTON**. The **DISPOSABLE LAMP CARTRIDGE** is located at the **APPLICATOR TIP**.

1.2. Package contents

Upon opening the Flash N' Go™ package, you will find the following parts:

- Flash N' Go™ **BASE UNIT** and **APPLICATOR**
- An AC cord
- A second **DISPOSABLE LAMP CARTRIDGE**
- This User's Manual and a Quick Start Guide leaflet
- An Instructional DVD

1.3. The Disposable Lamp Cartridge

The Flash N' Go™ **DISPOSABLE LAMP CARTRIDGE** can fire 750 light pulses (regardless of the energy level of these pulses). For average body size 750 light pulses could cover for example:

One treatment sessions of 2 Legs, 2 Arms, 2 Underarms, and Bikini line

Or

Two treatments sessions of 2 legs and 2 Underarms

Or

Ten treatments sessions of 2 Underarms and Bikini line

Depending on the body areas you wish to treat, you will typically need to use from one to seven **DISPOSABLE LAMP CARTRIDGES** during the first year of using Flash N' Go.

2. Flash N' Go™ Intended Use

Flash N' Go™ is intended for removal of unwanted hair.

Flash N' Go™ may be used to remove unwanted body hair. Ideal body areas for Flash N' Go™ use includes the underarms, bikini line, arms and legs.

3. Safety with Flash N' Go™

3.1. With Flash N' Go™ Safety Comes First

- **HPL™ technology in Flash N' Go™ - Superior safety with lower energy level**
Home Pulsed Light™ technology is able to achieve long-term hair removal results at a fraction of the energy level used in other light-based hair removal equipment. The low energy used in Flash N' Go™ reduces its potential to cause harm or complications, and contributes to your overall safety.

- **Flash N' Go™ protects your eyes**

The Flash N' Go™ *APPLICATOR* has a **built-in safety feature** for eye protection. It has been designed so that a light pulse can not be emitted when the *APPLICATOR* is facing open air. The safety switch is activated only when the *APPLICATOR TIP* is in full contact with the tissue and pressed.

- **Flash N' Go™ protects your skin**

Flash N' Go™ comes with a *SKIN COLOR SENSOR* enabling use only on lighter skin complexions to ensure skin safety.

Furthermore, the 2X3cm² *OPTICAL LENS* through which pulses of light are delivered is recessed inside the *LAMP CARTRIDGE* at the *APPLICATOR TIP*. This enables Flash N' Go™ to protect your skin by avoiding direct contact between the *OPTICAL LENS* and the skin.

4. Contraindications

Important Safety Information – Read Before Use!


Flash N' Go™ is not designed for everyone. Please read and consider the information in the following section before use. For further information and personalized advice you may also visit www.silk'n.com or your local Flash N' Go™ domain.

If you are unsure if Flash N' Go™ is safe for you to use, please ask your doctor or dermatologist!

- **DO NOT USE Flash N' Go™ on naturally dark skin complexion!**

Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™.

Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).



DO NOT USE Flash N' Go™ on naturally dark skin!

A unique *SKIN COLOR SENSOR* is embedded in Flash N' Go™ to measure the treated skin complexion at the beginning of each session and occasionally during the session. *SKIN COLOR SENSOR* ensures that pulse will be emitted on suitable skin tones.

- **DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure!**

Tanned skin particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk when using Flash N' Go™.

Using Flash N' Go™ on skin that has been recently exposed to sunlight can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).



DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure! Such use can cause serious burns or skin injury. Avoid exposure to the sun for 4 weeks before your Flash N' Go™ treatment!

The unique Flash N' Go™ *SKIN COLOR SENSOR* will also help you avoid treating sun exposed skin.

- **DO NOT USE Flash N' Go™ on the face or neck**

Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair.

Using Flash N' Go™ to remove facial hair may cause serious eye injury and may stimulate facial hair growth.



Flash N' Go™ is not recommended for use on the face or neck.

- **NOTE!** Flash N' Go™ is not effective on light hair

The Flash N' Go™ device is not effective on naturally white, grey, blond and red body hair. If your body hair is of these colors, Flash N' Go™ will not work on you.

5. When to avoid using Flash N' Go™?

Certain conditions may prevent the use of Flash N' Go™ temporarily. **DO NOT USE Flash N' Go™** if any of the following currently apply to you:

- If you are pregnant or nursing (lactating).
- If you were exposed to strong sunlight or an artificial tanning machine during the past 28 days.
- If you have a tattoo or permanent makeup on the area to be treated.
- If you have dark brown or black spots, such as large freckles, birth marks, moles or warts on the area to be treated.
- If you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sore in the area to be treated. Wait for the effected area to heal before using Flash N' Go™.
- If you have a history of keloidal scar formation, a known sensitivity to light (photosensitivity) or are taking medication that makes the skin more sensitive to light, including non-steroidal anti-inflammatory agents, (e.g., aspirins, ibuprofens, acetaminophen), tetracyclines, phenothiazines, thiazide, diuretics, sulfonyluraes, sulfonamides, DTIC, fluorouracil, vinblastine, griseofulvin, Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, Accutane® and/or topical retinoids.
- If you have abnormal skin conditions caused by diabetes, for example, or other systemic or metabolic diseases
- If you are currently or have recently been treated with Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, topical retinoids or azelaic acid.

- If you have been treated with Accutane® (isotretinoin) within the past 6 months.
- If you have been on a steroid regimen within the past 3 months.
- If you have a history of herpes outbreaks in the area of treatment, unless you have consulted your physician and received preventative treatment before using Flash N' Go™.
- If you suffer from epilepsy.
- If you have an active implant, such as a pacemaker, incontinence device, insulin pump, etc.
- If you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.
- If you have a history of skin cancer or areas of potential skin malignancies.
- If you have received radiation therapy or chemotherapy treatments within the past 3 months.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.

If you are unsure if Flash N' Go™ is safe for you to use, please consult with your doctor or dermatologist.

6. Precautions – How to use Flash N' Go™ Safely

- **Choose your energy levels CAREFULLY!**

Energy level refers to the intensity of the light pulse that is projected on your skin during use, from the lowest level (-) to the highest level (+). *INDICATOR LIGHTS* on the *CONTROL PANEL* illustrate the energy level at which the machine is set. As the energy level increases, so do the results of Flash N' Go as well as the risk of side effects (see "Possible Side Effects" below).

Always begin your first use of Flash N' Go™ at the lowest energy setting (one light at "-")!

Only if you experience little or no discomfort during and after use of Flash N' Go™ at the lowest energy setting, raise the energy level by one *INDICATOR LIGHT* the next time you use Flash N' Go™, and so on for each subsequent hair removal session.

For detailed instructions on energy level setting see "Energy Level" box in section 4.4 - "Treating with Flash N' Go™ for the first time".

- **Avoid adverse effects!**

Do not treat the same area of skin more than once per hair removal session!

Try to avoid overlapping pulses!

If your skin blisters or burns, **STOP USE IMMEDIATELY!**

- **Avoid complications after use of Flash N' Go™!**

Do not expose treated areas of skin to the sun!

Sun exposure includes unprotected exposure to direct sunlight of over 15 minutes constantly, or unprotected exposure to diffuse sunlight of over 1 hour constantly.

To protect recently treated skin when exposed to sunlight, be sure to thoroughly apply sunscreen SPF 30 or higher, for 2 weeks after each hair removal session.

- Always shave the area to be treated and make sure that the skin is clean and dry before using the Flash N' Go™.
- Cover birthmarks and tattoos before Flash N' Go™ application.
- Cover dark brown or black spots, such as large freckles, birth marks, moles or warts before Flash N' Go™ application.
- Never look directly at the light coming from the Flash N' Go™ *APPLICATOR* and *LAMP CARTRIDGE*.
- Do not use Flash N' Go™ on nipples and genitals (male or female).
- Do not use Flash N' Go™ on any body site where you might later want hair.
- Do not use Flash N' Go™ for any purpose other than hair removal.
- Never point the Flash N' Go™ *APPLICATOR* in an attempt to emit a light pulse into open space. Always make sure that the *APPLICATOR* is pointed at, and in full contact with the skin during application.
- Remove the Flash N' Go™ *APPLICATOR* from the skin if either the skin or the *APPLICATOR* is too hot.
- Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone to clean the skin before using Flash N' Go™.
- Use of Flash N' Go™ may cause temporary pigmentation changes (See "Possible Side Effects" below).
- Keep this device out of the reach of children. Do not use Flash N' Go™ on children or allow children to use it.

7. Reducing the risk of injury

As with any electrical device, certain precautions must be taken in order to ensure your safety when using Flash N' Go™.

- **Keep Flash N' Go™ away from water!**
Flash N' Go™ is an electrical device. As such it should always be kept away from water.
Do not place or store Flash N' Go™ where it can fall or be pushed into a tub, sink or any other vessel containing water. Do not place in, or drop into water or any other liquid.
This may cause severe electrocution.
Do not use Flash N' Go™ while bathing.
Do not use Flash N' Go™ if it becomes damp or wet.
Do not reach for Flash N' Go™ if it has fallen into water.
Unplug Flash N' Go™ immediately if it has fallen into water.



Keep Flash N' Go™ away from water!

- **Never open Flash N' Go™!**

Do not attempt to open or repair your Flash N' Go™ device. Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

- Trying to open Flash N' Go™ may also damage the device and will void your warranty.
Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.
- Use Flash N' Go™ only for its intended use and as described in its manual.
- Flash N' Go™ should never be left unattended when plugged into an outlet.
- Do not operate Flash N' Go™ if it has a damaged cord or plug and keep the power cord away from heated surfaces.
- Do not use Flash N' Go™ if you see or smell smoke when it is in use.
- Do not use Flash N' Go™ if it is not working properly or if it appears damaged.
- Do not use Flash N' Go™ if the fan vent in its *APPLICATOR* is cracked, coming off or missing.
- Do not use Flash N' Go™ if the *SKIN COLOR SENSOR* in its *APPLICATOR* is cracked, or broken.
- Do not use Flash N' Go™ if the outer shell is cracked or is coming apart.
- Do not use Flash N' Go™ with a damaged *DISPOSABLE LAMP CARTRIDGE*, or if its *OPTICAL LENS* is cracked, chipped or missing.
- Always unplug Flash N' Go™ from the electrical outlet immediately after use.
- Unplug Flash N' Go™ before cleaning.
- Do not use Flash N' Go™ with any attachments or accessories not recommended by Home Skinovations Ltd.

8. Risks of using Flash N' Go Flash N' Go™

When used according to the instructions, side effects and complications associated with use of Flash N' Go™ are uncommon. However every cosmetic procedure,

including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that can occur with pulsed light hair removal systems designed for home use.

| Adverse Event | Likelihood of the Adverse Event | Adverse Effect | Likelihood of the Adverse Effect |
|--|---------------------------------|-----------------------------------|----------------------------------|
| Stacking or Overlapping of multiple pulses on the same skin spot | Minor | Minor Skin Discomfort | Minor |
| | | Skin Redness | Minor |
| | | Increased Sensitivity of the Skin | Minor |
| | | Skin Wounds and Burns | Rare |
| | | Scarring | Rare |
| | | Pigment Changes | Rare |
| | | Excessive Redness and Swelling | Rare |
| | | Infection | Negligible |
| | | Bruising | Negligible |

- **Minor Skin Discomfort**

Although home pulsed light hair removal is generally very well-tolerated, most users do feel some mild discomfort during use, usually described as being a mild stinging sensation on the treated skin areas. The stinging sensation usually lasts during the time of the application itself or for a few minutes thereafter. Anything beyond this minor discomfort is abnormal and means that either you should not continue to use Flash N' Go™ because you are unable to tolerate the hair removal application, or that the energy level setting is too high.

- **Skin Redness**

Your skin may become red right after using Flash N' Go™ or within 24 hours of using Flash N' Go™. Redness generally clears up within 24 hours. See your doctor if redness does not go away within 2 to 3 days.

- **Increased Sensitivity of the Skin**

The skin of the treated area is more sensitive so you may encounter dryness or flaking of the skin.

- **Skin Wounds and Burns**

Very rarely, burns or wounds to the skin can occur following the application. The burn or wound can require a few weeks to heal and, extremely rarely, may leave a noticeable permanent scar.

- **Scarring**

Although very rare, permanent scarring may occur. Usually when scarring occurs it is in the form of a flat and white lesion on the skin (hypotrophic). However, it can be large and red (hypertrophic) or large and extend beyond the margins of the injury itself (keloid). Subsequent aesthetic treatments may be required to improve the appearance of the scar.

- **Pigment Changes**

Flash N' Go™ targets the hair shaft, in particular the pigmented cells in the hair follicle and the hair follicle itself. Nevertheless there is risk of temporary hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. This risk of changes in skin pigmentation is higher for people with darker skin tones.

Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.

- **Excessive Redness and Swelling**

In rare cases treated skin may become very red and swollen. This is more common in sensitive areas of the body. The redness and swelling should subside within 2 to 7 days and should be treated with frequent applications of ice. Gentle cleansing is OK, but one should avoid exposure to sun.

- **Infection**

Infection of the skin is exceedingly rare but is still a possible risk following a skin burn or wound caused by Flash N' Go™.

- **Bruising**


Very rarely, use of Flash N' Go™ may cause blue-purple bruising that can last 5 to 10 days. As the bruising fades, there may be a rust-brown discoloration of the skin (hyperpigmentation) that can be permanent.

9. Long Term Hair Removal the Flash N' Go™ Way

Flash N' Go™ is a personal light-based device for long-term hair removal. The process of laser and light-based hair removal is well known and established. It has been proven in clinical use around the world for over 15 years as a safe and effective way to achieve long-term hair reduction.

9.1. How does light remove hair?

Light-based hair removal is based on the theory of selective photothermolysis in which optical energy is used to disable hair growth. In order to achieve such thermal effect the hair shaft needs to selectively absorb light energy and transform it into heat. This selectivity is achieved when optical energy that is delivered to the tissue is mostly absorbed by hair shaft pigment, while the skin and the surrounding tissue stays cool. Melanin is the pigment in the hair shaft that is responsible for the absorption of the light, which generates the heat that eventually disables hair growth. Therefore the more melanin present in the hair (i.e. the darker the hair) the more light that can be absorbed and the more effective light can be at removing hair.




Before Flash N' Go™ Application
Go™ Application

During Flash N' Go™ Pulse After Flash N'

9.2. How does the hair growth cycle impact light-based hair removal?


Every hair in our body goes through the three phases of the hair growth cycle: Anagen, Catagen and Telogen. These phases have an important impact on how the process of light-based hair removal works.

Anagen is the hair growth phase while Catagen and Telogen, both are resting phases.



Anagen-Growth phase Catagen-Resting phase Telogen-Resting phase

The time it takes to complete a full hair growth cycle varies from person to person and the location of the hair on the body, but is typically 18-24 months. At any given moment the majority of the hair follicles in any skin area are in the resting phases. These resting hairs cannot be affected by Flash N' Go™. However, hairs in the growing Anagen phase will respond to Flash N' Go™ applications. It is important to understand that it may take a full hair growth cycle to realize complete hair removal results with Flash N' Go™.



9.3. Plan your Flash N' Go™ hair removal for best results

A typical full hair growth cycle may take 18-24 months. During this time multiple Flash N' Go™ sessions may be required in order to achieve long term hair removal.

The efficiency of hair removal varies from person to person according to body area, hair color, and how Flash N' Go™ is used.

Typical Flash N' Go™ hair removal plan during a full hair growth cycle:

- The first **3-4 hair removal sessions** with Flash N' Go™ will be approximately **two weeks apart**.
- Hair removal **sessions 5-7** with Flash N' Go™ will be approximately **four weeks apart**.
- After that you will typically use Flash N' Go™ again **from time to time** if and when needed, until long-term results are achieved.

This is the recommended treatment schedule that has proven to produce the best results, but you may plan your personal treatment schedule differently and still get satisfactory results.

Note that treating the same body part more often than once in two weeks, will not improve hair removal results.

Typical maintenance with Flash N' Go™ after achieving long-term hair removal:

Due to hormonal or other physiological changes dormant hair follicles may become active. Maintenance hair removal sessions with Flash N' Go™ may be required from time to time.

10. How To Use Flash N' Go™

10.1. Flash N' Go™ Setup

1. Remove Flash N' Go™ *BASE UNIT*, *APPLICATOR* and other components from box.
 2. Insert the *APPLICATOR* into its cradle in the Flash N' Go™ *BASE UNIT*.
 3. Verify that the *DISPOSABLE LAMP CARTRIDGE* is inserted correctly into the *APPLICATOR*.
 4. Plug the power cord into the Flash N' Go™ *BASE UNIT* socket.
 5. Plug the other end of the power cord into an electrical outlet.
- Your Flash N' Go™ is now ready to start.

10.2. Treating with Flash N' Go™ for the first time

The skin should be **shaved, clean, dry** and free of any powders, antiperspirants or deodorants.

1. Press the *POWER ON/OFF SWITCH*. The *POWER INDICATOR LIGHT* will turn on and a fan sound (similar to the sound of a hairdryer) will start.
 2. Approximately 3.5 seconds after pressing the *POWER ON/OFF SWITCH*, the *READY INDICATOR LIGHT* will turn on and the device beeps. The device sets itself to deliver the lowest energy level pulses. The device is then ready for you to trigger the first pulse.
 3. If the energy level should be higher than the lowest, press the *ENERGY LEVEL SETTING BUTTONS* using the "-" or "+" to respectively decrease or increase the energy level, until the desired energy level is set. *ENERGY LEVEL INDICATOR LIGHTS* will coincide with the energy level setting. (See "Energy Level" frame).
 4. Apply the *APPLICATOR* firmly to the skin, making sure the skin is spread evenly and smoothly. Make sure the *APPLICATOR TIP* is in full contact with the skin.
 5. Press the *PULSE BUTTON*. The device will first determine the color of your skin.
If the color of the skin is light enough for safe application, the device will flash a pulse of light onto your skin. You will see a bright flash of light and simultaneously hear a subtle pop sound, which is a normal noise for the device. You may feel a mild sensation of warmth and tingling.
Flash N' Go™ will immediately recharge for the next pulse. After 3.5 seconds the *READY INDICATOR LIGHT* will turn on again and the device will beep.
If you see no light pulse, and the 5 green energy level LEDs blink for 3 seconds, your skin tone is too dark for safe application. Try using the device on a different body part or contact the Flash N' Go support.
 6. Remove the *APPLICATOR TIP* from the treated area of skin.
 7. Move the *APPLICATOR TIP* to another area of skin.
Use the pressure marks the device just made in your skin to guide you for proper positioning of the next pulse, avoiding both gaps and overlaps between pulses (See "Covering skin areas" frame).
- Warning: Do not treat the same area of skin more than once per hair removal

session! Treating the same area of skin more than once per session increases the likelihood of adverse effects.

8. Repeat the process starting again with Step 5.

Energy Level

Energy level determines the intensity of the Flash N' Go™ light pulse delivered to your skin, from the lowest level (-) to the highest level (+). Corresponding *INDICATOR LIGHTS* on the *CONTROL PANEL* represent the increases in energy. As energy level increases, so does hair removal results as well as the risk of possible side effects and complications.

Always start your first Flash N' Go™ hair removal session at the lowest energy setting.

If you experience little or no discomfort during and after the hair removal session using the lowest energy level, raise the energy level by one *INDICATOR LIGHT* at the next session, and so on for each subsequent session.



Do not raise the energy level if you experience abnormal discomfort during or after treating with Flash N' Go™ (See section 2.6 - "Possible Side Effects").

Do not raise the energy level during hair removal session even if you experience no discomfort. Discomfort may also appear some time after the session.

Whenever Flash N' Go™ is turned on its energy level will automatically be reset to the lowest energy level. Only one *ENERGY LEVEL INDICATOR LIGHT* will be on.


To set the energy level, press the *ENERGY LEVEL SETTING BUTTONS* using "-" or "+" to respectively decrease or increase the energy level. The number of *ENERGY LEVEL INDICATOR LIGHTS* will coincide with the change in energy level.

For your safety, when using Flash N' Go™ for the very first time, the system is automatically set to deliver the first 50 pulses at the lowest energy level, and the next 200 pulses at up to level 3.

To cancel these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 2 consecutive "beep" sounds. To restore these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 3 consecutive "beep" sounds.

Covering skin areas

Flash N' Go™ pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas



When applied to the skin, the Flash N' Go™ *APPLICATOR TIP* is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.



Try to avoid overlapping pulses!

Do not treat the same area of skin more than once per hair removal session!



If your skin blisters or burns, STOP USE IMMEDIATELY!

10.3. What to Expect when treating with Flash N' Go™?

For many people, using Flash N' Go™ may be their first experience with a light-based device designed for home use. Flash N' Go™ is simple to use, and hair removal sessions go by quickly. During a Flash N' Go™ session it is **normal** to experience and feel:

- **A Flash of Light** – The bright light of Flash N' Go™ will not harm the eyes when applied to non-facial sites, and special eye-protection is not needed when using Flash N' Go™.
- **A Fan Noise** – The cooling fan in Flash N' Go™ makes noise similar to a hairdryer. This is normal.
- **A Pop Sound with Each Pulse** – When a pulse of HPL™ light is activated, it is normal to hear a subtle pop sound simultaneously with the flash of light.
- **Moderate Pressure of the *APPLICATOR*** – This is necessary and helpful for placement of adjacent pulses of light, and is part of the unique safety feature of Flash N' Go™.
- **A Sensation of Warmth and Tingling** – During each pulse of light it is normal to feel a mild sensation of warmth and tingling from the light energy. Remember it is important to always use low energy settings for initial hair removal sessions. You may feel some warmth for up to an hour after your Flash N' Go™ session.
- **Some Mild Red or Pink Color** – During and just after your Flash N' Go™ session it is not uncommon to see some very mild, pink-like color of the skin. This is usually most noticeable around the hairs themselves. However if you see full redness of the skin, blistering or burns stop use of Flash N' Go™ immediately.

10.4. After treating with Flash N' Go™

- When Flash N' Go™ session has been completed turn Flash N' Go™ off by pressing the *POWER ON/OFF SWITCH*. (Be sure to remember the last energy level setting you used, as it will not be restored when turning Flash N' Go™ on again.)
- Unplug the power cord from the electrical outlet.
- After each hair removal session it is recommended that you clean your Flash N' Go™ device, especially the *APPLICATOR TIP* (See section 5.1: "Cleaning Flash N' Go™").
- After cleaning, it is recommended to store your Flash N' Go™ device in its original box, and keep it away from water.



Skin care following hair removal session

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last Flash N' Go™ session.

Side effects and complications

Some patients may experience pigmentation changes resulting from treating with Flash N' Go™. These effects, if they occur, are generally mild and transient.



In case you experience any complication (See section 2.6 Possible Side Effects using Flash N' Go™) ***please contact your physician immediately.***

11. Maintenance of Flash N' Go™

11.1. Cleaning Flash N' Go™

After each hair removal session, it is recommended to clean your Flash N' Go™ device, and especially the *APPLICATOR TIP*.

- Unplug Flash N' Go™ before cleaning.
- Use a dry, clean cloth and a specially formulated cleaner for electronic equipment to gently wipe Flash N' Go™ surface, and especially the *APPLICATOR TIP*.
Never immerse Flash N' Go™ or any of its parts in water!

11.2. Replacing the LAMP CARTRIDGE


The Flash N' Go™ *DISPOSABLE LAMP CARTRIDGE* can fire 750 light pulses that would typically cover 2 Legs, Arms, Under arms, and Bikini line. Pulse intensity is determined only according to the energy level setting of the device. There is no decrease of energy during the usable lifetime of the *LAMP CARTRIDGE*.

When a *DISPOSABLE LAMP CARTRIDGE* has reached 90% of its possible lifetime the *CARTRIDGE 90% USAGE INDICATOR LIGHT* on the *CONTROL PANEL* will be activated, indicating that the *DISPOSABLE LAMP CARTRIDGE* should soon be replaced.

Once all 750 pulses in a *DISPOSABLE LAMP CARTRIDGE* have been used the *CARTRIDGE 100% USAGE INDICATOR LIGHT* will be activated, indicating that further pulses can not be emitted.

In order to continue the hair removal session, it will now be necessary to replace the *DISPOSABLE LAMP CARTRIDGE*.

Note: A *DISPOSABLE LAMP CARTRIDGE* should also be replaced if large spots appear inside it or if its *LENS* is broken.



To replace a *DISPOSABLE LAMP CARTRIDGE*:

1. Press the *POWER ON/OFF SWITCH* to turn Flash N' Go™ off.
2. Unplug the power cord from the electrical outlet.
3. Grasp the used *DISPOSABLE LAMP CARTRIDGE* on both sides, pull it out of the

12. Troubleshooting

12.1. "My Flash N' Go™ does not start."

- Make sure the power cord is properly connected to the Flash N' Go™ device.
- Make sure the power cord is plugged into an electrical outlet on the wall.

12.2. "A light pulse is not emitted when I press the PULSE SWITCH."

- Make sure that you have good contact with the skin and that the *APPLICATOR TIP* is evenly and firmly pressed to the skin. For your safety, the *PULSE BUTTON* will activate a pulse only if the *APPLICATOR TIP* is firmly pressed against the skin.
- Skin color sensor stopped the light pulse, your skin may be too dark for using the Flash N' Go. Try to apply Flash N' Go™ again.
- Check the *CARTRIDGE 100% USAGE INDICATOR*. If it is on, disconnect Flash N' Go™ and replace the *DISPOSABLE LAMP CARTRIDGE*.
- Make sure that the *READY INDICATOR* is on.
 - a. If within 10 seconds the *READY INDICATOR* remains off turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - b. If the problem persists, contact your local Flash N' Go™ Customer Service Center.
- Check the *SYSTEM WARNING INDICATOR LIGHT*.
 - c. If the light is on turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - d. If the problem persists, contact your local Flash N' Go™ Customer Service Center.

If these problems persist, contact your local Flash N' Go™ Customer Service Center.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.

Trying to open Flash N' Go™ may also damage the device and will void your warranty.

Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.

13. Customer Service

For more information about Flash N' Go™ please enter www.silkn.com.

If your Flash N' Go™ is broken, damaged, in need of repair, or for any other Flash N' Go™ user assistance, please contact Flash N' Go™ Customer Service:

1-877-DO-SILKN / 1-877-367-4556

contact@silkn.com

14. Frequently Asked Questions


1. **Does Flash N' Go™ really work?**
Yes. In clinical trials held by physicians, Flash N' Go™ was proven to safely achieve long-term hair removal results.
2. **Where on my body can I use Flash N' Go™?**
The Flash N' Go™ device has been designed for body hair removal anywhere below the neck. The most common areas treated with Flash N' Go™ are: legs, underarms, arms and bikini line. It is not recommended to use Flash N' Go™ on the face. Flash N' Go™ is not suitable for everyone. We recommend that you read all the Flash N' Go™ warnings and contraindications in this User's Manual.
3. **How long does a Flash N' Go™ treatment session take?**
The time can vary depending on the area of the body treated. A full leg can take up to 30 minutes, or two underarms could take 10 minutes. Because Flash N' Go™ runs on regular electric power it can be used for as long as needed to complete a full hair removal session of the desired body part(s).
4. **Is Flash N' Go™ safe?**
Flash N' Go™ has been designed with your safety in mind, and tested and approved by top dermatologists and plastic surgeons to meet their safety standards for a home-use device. But like any skin product or electronic device, one must use according to the operating instructions and user precautions.
5. **Will Flash N' Go™ hurt?**
When used properly most users of Flash N' Go™ report feeling a slight sensation of heat at the time of the pulse of light. Users with thicker and darker hairs may feel slightly more discomfort, but this discomfort subsides once the hair removal session is completed. For your convenience Flash N' Go™ has five setting levels that can be used according to your sensitivity.
6. **How often should I use Flash N' Go™?**
Hair removal sessions with Flash N' Go™ should be spaced every two weeks for the first three to four sessions. After that hair removal sessions should be done if hairs have grown back, until the desired results are achieved.
7. **Is Flash N' Go™ effective on white, grey or blonde hairs?**
Flash N' Go™ works best on darker hair types, or hair that contains more melanin. Melanin is the pigment that gives hair and skin their color, and will absorb light energy. Black and dark brown respond the best, although brown and light brown hairs will also respond but typically require more hair removal sessions. Red may show some response. White, grey or blonde hairs usually don't respond to Flash N' Go™ though some users have noted results after multiple hair removal sessions.
8. **Can I use Flash N' Go™ on brown or black skin**
Do not use Flash N' Go™ on naturally dark skin complexion! Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™. Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).
9. **When will I see results from Flash N' Go™?**
As with any light-based or laser hair removal device, results are not immediate, and in fact you may not think anything happened at all. Hair may sometimes appear to be growing back after a hair removal session, but typically after two weeks many of these hairs will simply fall out. Additionally hairs grow in three different stages and only hairs in an active growth stage will be affected by Flash N' Go™. This is one of the main reasons that multiple sessions are required to achieve the desired result.

10. **Can a man use Flash N' Go™?**
Though designed exclusively for women, Flash N' Go™ may be suitable for use by men. However, hairs on men, typically those on the chest, will require more hair removal sessions than that of women to get the desired results. As with women, using Flash N' Go™ on men's facial hair such as beard and mustache, is not recommended.
11. **Why is my hair growing, even though I treated it a week ago?**
It is quite common for hair to appear as if it is still growing up to two weeks after a hair removal session with Flash N' Go™. This process is known as "ejection" and at around two-weeks you'll see that these hairs simply fall out or slide out with a slight tug. (We don't however recommend pulling on the hairs – just let them come out naturally.) It is also possible that some hairs, due to missed application or different stages of growth, were not affected by Flash N' Go™. These hairs will be treated in follow-up sessions, and hence the reason multiple hair removal sessions are needed to get the best result with Flash N' Go™.
12. **I've heard that some hairs grow back lighter and finer after light based hair removal?**
This phenomenon is well documented amongst aestheticians and doctors using light and laser devices for hair removal. It is possible that some hairs will grow back lighter and finer after hair removal with Flash N' Go™. Usually these hairs are a fraction of what was originally there, and continued use may have a desirable effect on them.
13. **Why can't I treat myself if I have an "active" suntan?**
Do not use Flash N' Go™ on tanned skin or after sun exposure! Tanned skin and particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk of adverse effects when using Flash N' Go™ including burns, blisters, and skin color changes (hyper- or hypopigmentation).
14. **Is long-term use of Flash N' Go™ dangerous for my skin?**
The use of light and laser energy in aesthetic medicine has been well documented for over 15 years in professional peer-reviewed journals, and by well respected institutions like the Mayo Clinic. These journals and institutions have not reported any side-effects or damage from long-term use of light and laser device.
15. **Can I use Flash N' Go™ to remove my chin hair or elsewhere on my face?**
It is not recommended to use Flash N' Go™ on the face or neck. Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair. In addition, using Flash N' Go™ to remove facial hair may cause serious eye injury.
16. **How long should I wait to treat with Flash N' Go™ after unprotected exposure to the sun?**
One should wait 4 weeks before using Flash N' Go™ after unprotected exposure to the sun. However, if there is ever any uncertainty about sun exposure please contact your physician or Home Skinovations customer support.
17. **Should I do anything before using Flash N' Go™?**
Before any Flash N' Go™ session it is important to avoid sun exposure on the treated area for at least four weeks. A high level UV Sun Screen (SPF 50+) will help, as will clothing covering the treated area. The area to be treated should also be cleaned with mild soap and water, and the hairs shaved down to skin level.
18. **How should I care for the treatment area after using Flash N' Go™?**
The area treated with Flash N' Go™ can be cleaned and maintained with standard skin care products. Special care must be taken to avoid unprotected sun exposure. Strong sunscreens (50+ SPF) and covering clothing are suitable for protection from the sun.
19. **Should I pull the hairs out after treatment?**
No, let the hairs gradually fall out on their own. This may take up to 2 weeks.

15. Specifications

| | |
|---------------------------------|--|
| Spot size | 2cm x 3cm [6cm ²] |
| Speed | 1 pulse every 3.5 second: 1.7 cm ² /sec |
| Technology | New Home Pulsed Light™ |
| Max Energy Level | Max 5J/cm ² |
| Wavelength | 475-1200nm |
| Charging time / Power source | Continuous operation |
| Electrical requirements | 100-240VAC, 2A |
| Time needed to treat lower legs | 30 minutes |
| Operation and Safety | SKIN COLOR SENSOR seamlessly ensures use only on appropriate skin types. |
| package size | Height 9 inch, Width 9 inch, Depth 5.2 inch |
| System weight | 4 pounds |
| | |

16. Labeling



Class II equipment

Degree of protection against electric shock: type BF applied part

WEEE - Waste Electrical and Electronic Equipment

CE Mark

Follow operating instructions

This device comply with part 15 of the (FCC) Federal Communications Commission.

CSA Mark for USA and Canada

Degree of protection against ingress water: ordinary

This device is not suitable for use in the presence of flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.



HomeSkinovations

EC Declaration of Conformity

We Home Skinovations Ltd. declare under our sole responsibility that the home Skin treatment products

Brand Name: **Home Skinovations**
Model Name: **Flash N' Go, Disposable Lamp Cartridge**
Models Number: **AS100001E, AS100674B, AS100756A**

are fully in conformity with the essential requirements of Council Directive 89/336/EEC and 73/23/EEC amended by Directives 93/68/EEC.

This declaration is based on the full compliance of the products with the following European standards:

EMC: EN 61000-3-2:2000; (A2:2005),
EN 61000-3-3:1995; (A1:2001)
EN 55014-1:2000; (A2:2002),
EN 55014-2:1997; (A1:2001)

Electrical safety: EN 60335-1:2002; (A1:2004; A12:2006),
EN 60335-2-23:2003,
EN 60335-2-27:2003

Human Exposure to electromagnetic field radiated: EN 50364:2001

Manufacturer:
Home Skinovations Ltd.
Apolo Building, Shaar Yokneam,
POB 533, Yokneam 20692, ISRAEL
Tel: +972-4-9097470
Fax: +972-4-9097471

Authorized Representative:
Home Skinovations GMBH
Dr. Kurt Huber Str. 6
D-82031 Grünwald, GERMANY
Tel. +49-89-64919530
Fax +49-89-64919531

Signed:

Signed:

Gabi Lavi
General Manager
Yokneam, February 4, 2008

Dr. Amir Waldman
VP Regulatory & Clinical Affairs
Yokneam, February 4, 2008

For Information: On the basis of this declaration, these products and packaging will bear the following mark:



Flash N' Go™ is a **physician recommended** and innovative light-based device for hair removal in the **privacy of your home**. **Safe, easy to use** and **cost effective**, Flash N' Go™ achieves **marked clinical results** while offering personal convenience.

Flash N' Go™ utilizes the highly sophisticated and new **Home Pulsed Light™ (HPL™)** technology developed by **Home Skinovations Ltd.**



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, MD 20850

510(k) Submission for Flash N' Go

Home Skinovations Ltd.
August 10, 2008



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

Flash N' Go 510(k) Notification

(The Flash N' Go is an over the counter device intended for the removal of unwanted hair)

Table of contents

Cover Letter

Principal of Operation Chapter 1

Device description Chapter 2

Substantial Equivalence and Predicate Devices Chapter 3

User's Manual & Professional Instruction For Use Chapter 4

Risk Analysis and Safety Validation Chapter 5

Software Validation Chapter 6

Clinical Information Chapter 7

Premarket Notification Truthful and Accuracy Statement Chapter 8

"Indication for use" Page Chapter 9

510(k) Summary of Safety and Effectiveness Chapter 10



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

August 10, 2008

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, MD 20850

Subject: 510(k) Submission for Flash N' Go

Dear Sir/Madam,

Home Skinovations Ltd. is submitting this premarket notification for the purpose of obtaining FDA marketing clearance for its product. A table of contents for this 510(k) application is presented immediately after this cover letter. The following is information required under 21 CFR 807.87:

This is to notify you of the intention of Home Skinovations Ltd. to manufacture and market the following device:

Please send all correspondence by fax (+972-4-9097471), in addition to mail

| | |
|--------------------------------|--|
| Trade/Proprietary name: | Flash N' Go |
| Common name: | Light based hair removal system |
| Classification name: | The subject of this application is a Laser surgical instrument for use in general and plastic surgery and in dermatology (21 CFR 878.4810) |
| Product codes: | GEX |
| Class of device: | The device should be considered a class II device, as are the predicate devices. |
| Reason for Application: | New device. |
| Name and address of | Home Skinovations Ltd. |
| Manufacturer: | Home Skinovations Ltd. |



Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

Contact person:

Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471
Dr. Amir Waldman
VP regulatory affairs
Tel. +972-547-404616, +972-4-9097440
Fax +972-4-9097471
E-mail - amirw@silkn.com

Applicable standards:

Action taken by the person required to register to comply with the requirements of the Act under section 514 for performance standards or section 513 for special controls: There are no mandatory performance standards applicable to this device. There are voluntary performance standards that the Flash N' Go complies with as discussed in chapter 4.

Labeling/Promotional material:

The labeling and promotional materials for the Flash N' Go, and for the substantially equivalent predicate devices, are presented in later sections of this application.

Substantial equivalence predicate devices:

Flash N' Go is substantially equivalent to the following devices:

- ABC hair removal system, Indication for use: "The ABC hair removal system is an over the counter device intended for adjunctive use with shaving for hair removal sustained with periodic treatments". Manufactured by Palomar Medical Technologies Inc. and subject of K060839.
- Spectra hair removal laser system, Indication for use: The Spectra hair removal laser system is an

HomeSkinovations

Apolo building
POB 533
Yokneam 20692
ISRAEL
Tel: +972(4)9097440
Fax: +972(4)9097471

over the counter device intended for adjunctive use with shaving for hair removal sustained with oeriodic treatments. Manufactured by SpectraGenics, Inc. and subject of K053527.

- Silk'n, Indication for use: The Silk'n is intended for removal of unwanted hair by using selective photothermal treatment. The device is generally indicated for dermatological use. The Silk'n is specifically indicated for patient removal of unwanted hair by using selective photothermal treatment under the direction of a physician, after training by a healthcare professional. Manufactured by Home Skinovations Ltd. and subject of K072906.

Confidentiality:

The material in this application is considered confidential and is not meant for public release.

**Truth and Accuracy
Statement:**

A truthful and accurate statement, according to 21CFR 807.87(k), is presented after a separate page with the indications for use.

**SMDA Safety and
Effectiveness Summary:**

Summary of Safety and is presented following the table of contents.

The therapeutic light output of the Flash N' Go (b) (4)

(b) (4)

We trust that the information included in this application will be adequate to allow for its prompt review. Please contact me if there are any questions.

Sincerely yours,

Dr. Amir Waldman
VP regulatory affairs
Home Skinovations Ltd.



Chapter 1 Principle of operation

Introduction to Flash N' GO

The Flash N' Go device uses light energy in order to remove hairs. The Flash N' Go device is an over the counter device indicated for the removal of unwanted hair.

The treatment with Flash N' Go is based on the principle of *selective thermolysis*. According to this principle, parameters of optical energy (spectrum, exposure duration and energy density) are chosen and optimized to selectively treat hair and follicle without damaging the surrounding tissues.

Background

Optical energy can be used for thermal destruction of hair follicles. Optical energy penetrates through the upper skin and is absorbed by the pigment of the hair shaft and converted into heat energy, raising the hair temperature. The hair follicle structure has no optical contrast with surrounding tissue and does not selectively absorb light. Heat can be conducted from the shaft to the hair follicle thus creating heating of the follicle to temperatures that are higher than the temperature of the surrounding tissue. When the follicle reaches a high enough temperature, it is damaged and this leads to delay in hair growth.

In order to remove hair safely, the light energy type has to be selected in such a way that the skin or other healthy tissue absorbs the smallest portion of energy while the largest portion is absorbed by the hair shaft. The schematic representation of skin with hair is shown in Figure 1.

The term selective photothermolysis was coined to describe heating of pigmented targets by selectively absorbed pulses of light radiation. Laser-assisted photothermolysis is based on applying single wavelength, selectively absorbed by target^[1]. This principle was later extended to multiple wavelengths^[2].




Figure 1. Schematic representation of skin with hair.

It was shown that different wavelengths could be absorbed by different parts of the target, providing more uniform heating of the treated object. However, in both cases the treated target that absorbs the light should have a higher absorption coefficient than the surrounding tissue. In the present device, this principle is applied by using low energy pulses of light.

The basic elements that are necessary to achieve selective Photo-thermolysis are:

- Energy that reaches and is preferentially absorbed by the desired target structures.
- Exposure duration less than, or equal to, the time necessary for the target structures to cool.

Selective Photo-thermolysis calls for the use of appropriate light energy parameters such as wavelength range and pulse duration to maximize the absorption of hair follicle within the irradiated area and minimize the heating of nearby tissue.

Selective Photo-thermolysis is successfully applied by light sources such as the Ruby laser, Alexandrite laser, Diode laser, Filtered Xenon Flash Lamp. The Flash N' Go uses the same principle of selective Photo-thermolysis; the skin is exposed to light

energy that penetrates through it and is absorbed in the hair shaft and follicle being treated. The parameters of the light energy generated by the Flash N' Go'n are optimized to achieve maximum selectivity and therefore are safe and effective.

The penetration of light into the skin and the absorption of the different components of the skin and the hair shaft are a function of the wavelength range of light used. Taking all skin chromophores into account, the visible, near-Infrared wavelengths are the most suitable for hair removal treatments. Indeed, most of the cleared devices, lasers and pulse light use spectral range of 400-1200 nm.

Light penetration depth as a function of wavelength is shown in Figure 2. Strong light absorption by blood limits the penetration in the spectral range of 400-600 nm. At higher wavelengths light penetration increases due, mostly, to decrease in the scattering coefficient. At wavelengths longer than 1300nm, the strong light absorption by water prevents its penetration into the dermis.




Figure 2. Light penetration depth into the skin as a function of wavelength.

The second requirement is the use of wavelengths in which light is preferably absorbed by the hair shaft compared to the surrounding tissue (dermis and epidermis). The absorption coefficients of epidermis and dermis are shown in Figure 3. The

graph, based on the data of Svaasand^[3] and Anderson and Parish^[4], is for skin with a 1% blood content.




Figure 3. Absorption and scattering coefficient of epidermis and dermis with 1% of blood content in the wavelength range of 400-1200nm.

In the dermis, absorption in the wavelength range of 400-1200nm is determined by the blood content. In the range of 980-1000nm, dermal absorption is mainly due to water absorption peak. In contrast to the dermis, the absorption of the epidermis is mainly controlled by melanin content. Figure 3 clearly demonstrates this point since it is seen that the absorption coefficient in the epidermis does not show the influence of water or hemoglobin absorption.

The behavior of the absorption coefficient of hair is similar to that of the epidermis since it is mainly controlled by melanin. Light absorption for all hair colors decreases as wavelength increases. The absorption coefficient is directly related to the darkness of hair. Absorption is much higher for black hair than for blond hair. Figure 3 shows the difference between dermal and epidermal (melanin) absorption which is maximal at the shorter wavelengths around 400.

The result of this analysis is that effective and safe hair removal can be achieved at the spectral range of 475 to 1200nm.

Pulse width

Thermal conductivity plays an important role in optimizing selective Photothermolysis. The cooling time (thermal relaxation time) of biological objects with temperatures higher than the surrounding tissues is a function of their size. The cooling time, t , of a cylindrical object (such as a hair shaft) is calculated as:

$$t = \frac{d^2}{16\alpha} \quad (3.4)$$

where α is the thermal diffusivity and d is the diameter of the cylindrical object^[6].

The cooling time of a hair follicle is a function of its diameter. Cooling time varies by a factor of 10 between a 100 μm diameter hair and a 300 μm diameter hair.

Diffusivity of hair is about $0.5 \cdot 10^{-7} \text{ m}^2/\text{sec}$. Dependence of hair cooling time as function of hair follicle diameter is shown in Figure 4.

Because our targets are both hair shafts and hair follicles, the pulse duration should be shorter than the thermal relaxation time of the hair shaft. Human hairs have diameters in the range of 100 μ to 300 μ . Thus, pulse duration should be shorter than 10 milliseconds.

Importance of spot size

The geometry of irradiation plays an important role in achieving good penetration and controlling the scattering in the dermis. When the light is applied to the skin with a small spot size, the scattering of photons diffuses the beam rapidly. The fluence decreases very quickly as a function of depth. Using a spot size that is less than the light penetration depth, the energy distribution can be approximated to that of a spherical decay with the fluence decreasing, through scattering, proportionately to depth squared. With a large spot size, light penetration is more efficient since the "source" of photons has an almost planar geometry^[2].





Figure 4 Cooling time of a hair follicle vs. its diameter

Safety

(b) (4)



Summary

It was shown that the choice of the proper optical energy parameters (b) (4) is necessary to achieve safe and effective hair removal on all potential Flash N' Go users.

The optimal range of parameters for safe and effective hair removal as is implemented with the Flash N' Go device are summarized below:

| | |
|-----------|------------------------|
| Spectrum | 475-1200nm |
| Fluence | 3-5 J /cm ² |
| Spot size | 2x3 cm ² |

References

1. R. Rox Anderson, "Laser-tissue Interactions", in Cutaneous Laser Surgery by Mitchel P. Goldman and R.E. Fitzpatrick (1994, Mosby p.9)
2. G. Lask, S. Eckhouse, M. Slatkine, A. Waldman, M. Kreindel, V. Gotfried, "The role of laser and intense light sources in photo-epilation: comparative evaluation". Cutaneous Laser Therapy, 1999, 1:3-13.
3. Svaasand, L.T. Norvang, E.J. Fiskerstrand, E.K.S. Stopps, M.W. Beams, J.S. Nelson, "Tissue Parameters determining the Visual Appearance of Normal Skin and Port-wine stains". Lasers in Medical Science 1995, 10:55-65
4. R.R. Anderson, J.A. Parrish, "The Optics of Human Skin." Journal of Investigative Dermatology 1981, 77:13-19
5. M.G.C. Van Gemert, A.J. Welch. Time constant in thermal laser medicine. Lasers Surg. Med., 1989, 9:405-21

Chapter 2 Device Description

Flash N' Go hair removal system is basically a low power IPL system. The device has two filtered Xenon flash lamps that are mounted inside of a hand held applicator which is connected to a main console with a cord. The on-off switch is located on the main console, along with a switch to set the desired energy level, as well as a set of LEDs that indicate the system's operation mode. A trigger push-button is located on the applicator to emit light pulse. (b) (4)

1. (b) (4)

2.

Device specifications:

| | |
|-----------------------|---|
| Wavelengths range | 475-1200 nm |
| Energy density | 3- 5 J/cm ² |
| Spot size | 2*3 cm ² |
| Pulse repetition rate | 0.3 Hz |
| Weight | Console 1.5 Kg (3 lbs) Applicator .25 Kg (1 lbs) |

Control

The control of the system is carried out through a combination of hardware controls and a real-time controller that monitors the status and parameters of the device subsystems. The micro controller translates the treatment parameters selected by the user into electrical signals that generate the proper energy out of the Applicator.

The interface of the real-time controller with the operator is carried out through a control panel that shows the status of the system and is also used to select the treatment parameters.

The control functions and the structure of the hardware element of the control system are shown in *Figure 1*.




Figure 1. Control system structure and functional relationship.

Turn on and turn off sequence

A push-button switch on the front panel is used to turn the machine on by the user. The controller gets its power from the power line through an isolation transformer. Check sum of the software is carried out whenever the device turns on. By turning the mains switch off, the user immediately disconnects the system from the power line and the system is turned off. Whenever the treatment mode is initiated by the user through the user interface, the controller initiates the new power.

Control during treatment

(b) (4)



Safety function monitoring

(b) (4)



During the entire time that the machine is on, the operational functionality of the cooling system is checked and in the event that there is a problem with this subsystem- the system stops its operation until the safety mode is over. (b) (4)

(b) (4) in the applicator, and in the event of over- heating beyond a predetermined point, the systems stops its operation until the safety mode is over.

Quality System Specification

CERTIFICATE

Number: 2109191

The management system of:

Home Skinovations Ltd.

4 Hacharoshet street
P.O. Box 7028
44640 Kfar Saba
Israel

including the implementation meets the requirements of the standard:

ISO 13485:2003

Scope:

Design, manufacture, marketing and sales of home use skin treatment devices


The file that forms the basis of this certificate:

2109191-QUA


This certificate is valid until: December 1, 2010

Issued for the first time: December 10, 2007

KEMA Quality B.V.



drs. G.J. Zoetbrood
Managing Director



dr. ir. G.W. Bos
Certification Manager

© KEMA Quality B.V.

Integral publication of this certificate and adjoining reports is allowed.



Chapter 3: Substantial equivalence

The Flash N' Go is substantially equivalent (SE) to the following devices:

- ABC hair removal system, Indication for use: "The ABC hair removal system is an over the counter device intended for adjunctive use with shaving for hair removal sustained with periodic treatments". Manufactured by Palomar Medical Technologies Inc. and subject of K060839.
- Spectra hair removal laser system, Indication for use: The Spectra hair removal laser system is an over the counter device intended for adjunctive use with shaving for hair removal sustained with poeriodic treatments. Manufactured by SpectraGenics, Inc. and subject of K053527.
- Silk'n, Indication for use: The Silk'n is intended for removal of unwanted hair by using selective photothermal treatment. The device is generally indicated for dermatological use. The Silk'n is specifically indicated for patient removal of unwanted hair by using selective photothremal treatment under the direction of a physician, after training by a healthcare professional. Manufactured by Home Skinovations Ltd. and subject of K072906.

Labeling and promotional material of the substantially equivalent devices are presented in Appendices 5-1, 5-2, 5-3 of this chapter.

Substantial equivalence is demonstrated according to the criteria outlined in the 510(k) "Substantial Equivalence" Decision Making Process, as published in the FDA 510(k) Guidance (FDA 90-4158).

A comparison table of the subject and predicate devices is included in this discussion.

Intended use equivalence

The Flash N' Go device is an over the counter device intended for the removal of unwanted hair.

The Flash N' Go employs pulses of low energy light to achieve delayed hair growth.

This intended use is similar to the intended use of the Silk'n with respect to the clinical result: "The Silk'n Photoepilation system is intended for removal of unwanted hair by using a selective photothermal treatment. The device is generally indicated for dermatological use. The Silk'n is specifically indicated for patient removal of unwanted hair by using a selective photothermal treatment under the direction of a physician, after training by healthcare professional". (b) (4)

(b) (4)

The Flash N' Go intended use is similar to the ABC and the Spectra devices as an over the counter device: "the ABC hair removal system is an over the counter device intended for adjunctive use with shaving for hair removal sustained with periodic treatments". (b) (4)

(b) (4)

Technological equivalence

The Flash N' Go has technological characteristics that are substantially equivalent to predicate devices. (b) (4). ABC and Spectra devices are laseres at 810 nm with maximum energy density of 12 J/cm², and 20 J/cm².

No new types of safety or effectiveness questions are raised and differences in the technological characteristics between the Flash N' Go and the predicate devices can be assessed for their effects on safety and effectiveness.

A comparison is provided in the table below:

| <u>Technological Characteristic</u> | <u>Flash N' Go</u> | Silk'n | Spectra | ABC device | Comments |
|---|--------------------------------------|--------------------------------------|--------------------|-------------------|---|
| Nature and Source of Energy Imparted to Patient | Pulsed light via external applicator | Pulsed light via external applicator | Pulsed Diode laser | Pulse Diode Laser | All devices use optical energy as the power source for hair removal treatment |

| | | | | | |
|---|----------|----------|-----|-----|--|
| Wavelength range (nm) | 475-1200 | 475-1200 | 810 | 810 | All devices use Visible / Near IR Wavelengths range to achieve reduced hair growth |
| Maximum Light Energy (J/cm ²) | (b) | (b) | 20 | 12 | Same energy like the Silk'n |

Safety and Effectiveness equivalence

The Flash N' Go does not present any new safety and effectiveness issues as compared with the predicate devices.

(b) (4)

As far as effectiveness, (b) (4)

(b) (4)

In light of the results, we conclude that the Flash N' Go device is safe and effective as an over the counter device indicated for the removal of unwanted hair.

Flash N' GoTM

Personal Hair Removal at Home
with Home Pulsed LightTM




Flash N' Go™ is a physician recommended and innovative light-based device for hair removal in the privacy of your home. Safe, easy to use and cost effective, Flash N' Go™ achieves marked clinical results while offering personal convenience.

THE HPL™ INNOVATION


The Flash N' Go™ device utilizes the highly sophisticated and new **Home Pulsed Light™ (HPL™)** technology developed by Home Skinovations Ltd. especially for safe use at home.

HPL™ Epilation is based on the theory of selective photothermolysis in which **optical** energy is used to disable hair growth. Epilation with Home Skinovations' proprietary HPL™ technology is further benefited by a unique **acoustic** effect that enhances the normal process of epilation by photothermolysis.


To safely achieve the desired clinical results the hair shaft needs to selectively absorb HPL™ light energy and transform it into heat. Melanin is the pigment in the hair shaft that is responsible for the absorption of light, which generates the heat at the follicle. Selectivity is achieved when the HPL™ optical energy delivered to the tissue is mostly absorbed by pigment in the hair shaft, while the epidermis and the surrounding tissue remain cool. Additionally HPL™ technology creates an acoustic effect that enhances the overall hair removal result. Together the optical and acoustic effects of HPL™ technology are what enable the Flash N' Go™ device to be safely used by patients for home based hair removal.



● Before using Flash N' Go™



● During Flash N' Go™ pulse



● After Flash N' Go™ application

Flash N' Go™ ADVANTAGES

Physician Recommended

Tested over 12 months in a multi-center clinical study involving over 150 female patients.

Small and Portable

Flash N' Go™ is about the size of a hairdryer and it runs off of regular electric power, so treatment can be done anywhere, anytime and non-stop for as long as necessary.

Safe

When used according to instructions, hair removal with Flash N' Go™ is safe and with minimal discomfort compared to waxing or laser.

Easy to Use

Flash N' Go™ is easy to set up, and treatments are simple, clean and very quick.

Cost Effective

Using Flash N' Go™ will save you money when compared to extended use of professional laser treatments, waxing or shaving.

Marked clinical Results

Although small and portable, the HPL™ technology in Flash N' Go™ captures all the clinical benefits of light-based hair removal.

TECHNICAL SPECIFICATIONS

| | |
|---------------------------------|--|
| Spot size | 2cm x 3cm [6cm ²] |
| Speed | 1 pulse every 3.5 second: 1.7 cm ² /sec |
| Technology | Home Pulsed Light™ |
| Max Energy Level | 5 J |
| Wavelength | 475-1200 nm |
| Charging time / Power source | Continuous operation |
| Time needed to treat lower legs | 30 minutes |
| Operation and safety | Safety sensor tip enables maximum safety & control. Easy to use. |

CONTACT INFORMATION

Home Skinovations Ltd.

Apollo Building, Shaar Yekneam, POB 533 Yekneam 20692 Israel

Tel: +972-1-9097470, Fax: +972-1-9097471

info@silkn.com www.silkn.com

ISO 13485





DISCOVER TRIA

BUY TRIA

MEDIA & STORIES

TRIA

BUY TRIA

STORE LOCATOR

COMING SOON

Smooth, sleek, silky skin at the touch of a button

The TRIA Personal Laser Hair Removal System is the first FDA-cleared laser hair removal product on the market that delivers clinical results in the privacy, comfort, and convenience of your own home.

Kit Includes:

- TRIA laser handpiece
- TRIA Skin Sensor
- Battery charger
- *Instructions for Use* booklet
- Instructional DVD
- Warranty card
- Storage case

Benefits:

- No more endless shaving
- No more painful waxing
- No more inconvenient office appointments
- Laser hair removal in the privacy, comfort, and convenience of your own home
- Beautiful, smooth skin that looks and feels as though the hair was never there

Price: \$995.00 (shipping & handling included)


To buy the TRIA System:

Call: 1-877-321-TRIA (8742) or

Email us at clientservices@triabeauty.com

Want to experience the TRIA laser in person?

[Click here](#) to find a TRIA MD near you.



INTERNATIONAL SITES | ABOUT US | WARRANTY | TERMS OF USE | PRIVACY POLICY | AFFILI
® 2008 TRIABEAUTY



Helping physicians enhance
and maximize their hair removal business





Flash N' Go User manual

Copyright © 2008 Home Skinovations Ltd. All rights reserved.

Print date: February 2008

Home Skinovations Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd. is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd.

Data is subject to change without notification.

Home Skinovations Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd.

Specifications are subject to change without notice.

Home Skinovations Ltd.
Apolo Building, Shaar Yokneam
POB 533
Yokneam 20692
ISRAEL

Website:

www.silkn.com

e-mail:

info@silkn.com

Table of Contents

| | |
|---|-----------|
| 1. Before You Start | 3 |
| 1.1. What is Flash N' Go™? | 3 |
| 1.2. Flash N' Go™ intended use | 3 |
| 2. Safety with Flash N' Go™ | 3 |
| 2.1. With Flash N' Go™ Safety Comes First | 3 |
| 2.2. Is Flash N' Go™ for you? | 4 |
| 2.3. When to avoid using Flash N' Go™? | 6 |
| 2.4. Precautions – How to use Flash N' Go™ Safely | 7 |
| 2.5. Reducing the risk of injury | 8 |
| 2.6. Possible Side Effects using Flash N' Go™ | 9 |
| 3. Long Term Hair Removal the Flash N' Go™ Way | 11 |
| 3.1. How does light remove hair? | 11 |
| 3.2. How does the hair growth cycle impact light-based hair removal? | 11 |
| 3.3. What is different about Home Pulsed Light™ technology in Flash N' Go™? | 12 |
| 3.4. Plan your Flash N' Go™ hair removal for best results | 12 |
| 4. First Steps with Flash N' Go™ | 14 |
| 4.1. Flash N' Go™ device description | 14 |
| 4.2. Package contents | 14 |
| 4.3. Flash N' Go™ Setup | 15 |
| 4.4. Treating with Flash N' Go™ for the first time | 15 |
| 4.5. What to Expect when treating with Flash N' Go™? | 17 |
| 4.6. After treating with Flash N' Go™ | 18 |
| 5. Maintenance of Flash N' Go™ | 19 |
| 5.1. Cleaning Flash N' Go™ | 19 |
| 5.2. Replacing the Lamp Cartridge | 19 |
| 6. Troubleshooting | 21 |
| 6.1. "My Flash N' Go™ does not start." | 21 |
| 6.2. "A light pulse is not emitted when I press the PULSE SWITCH." | 21 |
| 7. Customer Service | 22 |
| 8. Frequently Asked Questions | 23 |
| 9. Specifications | 26 |
| 10. Labeling | 27 |
| EC Declaration of Conformity | 28 |

1. Before You Start

Before using Flash N' Go™ for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures.

We recommend you re-familiarize yourself with this User Manual before each use of Flash N' Go™.



Flash N' Go™ is a powerful electrical device. As such, it should be used with special attention to safety.

Please read all warnings and safety precautions before use, and strictly follow them when using Flash N' Go™

1.1. What is Flash N' Go™?

Flash N' Go™ is a light-based device for long-term hair removal designed for home-use.

1.2. Flash N' Go™ intended use

Flash N' Go™ is intended for removal of unwanted hair.

Flash N' Go™ may be used to remove unwanted body hair. Ideal body areas for Flash N' Go™ use includes the underarms, bikini line, arms and legs.

2. Safety with Flash N' Go™

2.1. With Flash N' Go™ Safety Comes First

- **HPL™ technology in Flash N' Go™ - Superior safety with lower energy level**

Home Pulsed Light™ technology is able to achieve long-term hair removal results at a fraction of the energy level used in other light-based hair removal equipment. The low energy used in Flash N' Go™ reduces its potential to cause harm or complications, and contributes to your overall safety.

- **Flash N' Go™ protects your eyes**

The Flash N' Go™ *APPLICATOR* has a **built-in safety feature** for eye protection. It has been designed so that a light pulse can not be emitted when the *APPLICATOR*

is facing open air. The safety switch is activated only when the *APPLICATOR TIP* is in full contact with the tissue and pressed.

- **Flash N' Go™ protects your skin**

Flash N' Go™ comes with a *SKIN COLOR SENSOR* enabling use only on lighter skin complexions to ensure skin safety.

Furthermore, the 2X3cm² *OPTICAL LENS* through which pulses of light are delivered is recessed inside the *LAMP CARTRIDGE* at the *APPLICATOR TIP*. This enables Flash N' Go™ to protect your skin by avoiding direct contact between the *OPTICAL LENS* and the skin.

2.2. Is Flash N' Go™ for you?

Important Safety Information – Read Before Use!

Flash N' Go™ is not designed for everyone. Please read and consider the information in the following section before use. For further information and personalized advice you may also visit www.silk'n.com or your local Flash N' Go™ domain.

If you are unsure if Flash N' Go™ is safe for you to use, please ask your doctor or dermatologist!

- **DO NOT USE Flash N' Go™ on naturally dark skin complexion!**

Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™.

Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

Not Safe



DO NOT USE Flash N' Go™ on naturally dark skin!

A unique *SKIN COLOR SENSOR* is embedded in Flash N' Go™ to measure the treated skin complexion at the beginning of each session and occasionally

during the session. *SKIN COLOR SENSOR* ensures that pulse will be emitted on suitable skin tones.

- **DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure!**

Tanned skin particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk when using Flash N' Go™.

Using Flash N' Go™ on skin that has been recently exposed to sunlight can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).



DO NOT USE Flash N' Go™ on tanned skin or after recent sun exposure! Such use can cause serious burns or skin injury. Avoid exposure to the sun for 4 weeks before your Flash N' Go™ treatment!

The unique Flash N' Go™ *SKIN COLOR SENSOR* will also help you avoid treating sun exposed skin.

- **Flash N' Go™ is not recommended for use on the face or neck**

Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair.

Using Flash N' Go™ to remove facial hair may cause serious eye injury and may stimulate facial hair growth.



Flash N' Go™ is not recommended for use on the face or neck.

- **NOTE! Flash N' Go™ is not effective on light hair**

The Flash N' Go™ device is not effective on naturally white, grey, blond and red body hair. If your body hair is of these colors, Flash N' Go™ will not work on you.

2.3. When to avoid using Flash N' Go™?

Certain conditions may prevent the use of Flash N' Go™ temporarily. **DO NOT USE Flash N' Go™** if any of the following currently apply to you:

- If you are pregnant or nursing (lactating).
- If you were exposed to strong sunlight or an artificial tanning machine during the past 28 days.
- If you have a tattoo or permanent makeup on the area to be treated.
- If you have dark brown or black spots, such as large freckles, birth marks, moles or warts on the area to be treated.
- If you have eczema, psoriasis, lesions, open wounds or active infections, such as cold sore in the area to be treated. Wait for the effected area to heal before using Flash N' Go™.
- If you have a history of keloidal scar formation, a known sensitivity to light (photosensitivity) or are taking medication that makes the skin more sensitive to light, including non-steroidal anti-inflammatory agents, (e.g., aspirins, ibuprofens, acetaminophen), tetracyclines, phenothiazines, thiazide, diuretics, sulfonluraes, sulfonamides, DTIC, fluorouracil, vinblastine, griseofulvin, Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, Accutane® and/or topical retinoids.
- If you have abnormal skin conditions caused by diabetes, for example, or other systemic or metabolic diseases
- If you are currently or have recently been treated with Alpha-Hydroxi Acids (AHAs), Beta-Hydroxi Acids (BHAs), Retin-A®, topical retinoids or azelaic acid.
- If you have been treated with Accutane® (isotretinoin) within the past 6 months.
- If you have been on a steroid regimen within the past 3 months.
- If you have a history of herpes outbreaks in the area of treatment, unless you have consulted your physician and received preventative treatment before using Flash N' Go™.
- If you suffer from epilepsy.
- If you have an active implant, such as a pacemaker, incontinence device, insulin pump, etc.
- If you have a disease related to photosensitivity, such as porphyria, polymorphic light eruption, solar urticaria, lupus, etc.
- If you have a history of skin cancer or areas of potential skin malignancies.
- If you have received radiation therapy or chemotherapy treatments within the past 3 months.
- If you have any other condition which in your physician's opinion would make it unsafe for you to be treated.

If you are unsure if Flash N' Go™ is safe for you to use, please consult with your doctor or dermatologist.

2.4. Precautions – How to use Flash N' Go™ Safely

- **Choose your energy levels CAREFULLY!**

Energy level refers to the intensity of the light pulse that is projected on your skin during use, from the lowest level (-) to the highest level (+). *INDICATOR LIGHTS* on the *CONTROL PANEL* illustrate the energy level at which the machine is set. As the energy level increases, so do the results of Flash N' Go as well as the risk of side effects (see "Possible Side Effects" below).

Always begin your first use of Flash N' Go™ at the lowest energy setting (one light at "-")!

Only if you experience little or no discomfort during and after use of Flash N' Go™ at the lowest energy setting, raise the energy level by one *INDICATOR LIGHT* the next time you use Flash N' Go™, and so on for each subsequent hair removal session.

For detailed instructions on energy level setting see "Energy Level" box in section 4.4 - "Treating with Flash N' Go™ for the first time".

- **Avoid adverse effects!**

Do not treat the same area of skin more than once per hair removal session!

Try to avoid overlapping pulses!

If your skin blisters or burns, **STOP USE IMMEDIATELY!**

- **Avoid complications after use of Flash N' Go™!**

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last hair removal session.

- Always shave the area to be treated and make sure that the skin is clean and dry before using the Flash N' Go™.
- Cover birthmarks and tattoos before Flash N' Go™ application.
- Cover dark brown or black spots, such as large freckles, birth marks, moles or warts before Flash N' Go™ application.
- Never look directly at the light coming from the Flash N' Go™ *APPLICATOR* and *LAMP CARTRIDGE*.
- Do not use Flash N' Go™ on nipples and genitals (male or female).
- Do not use Flash N' Go™ on any body site where you might later want hair.
- Do not use Flash N' Go™ for any purpose other than hair removal.
- Never point the Flash N' Go™ *APPLICATOR* in an attempt to emit a light pulse into open space. Always make sure that the *APPLICATOR* is pointed at, and in full contact with the skin during application.
- Remove the Flash N' Go™ *APPLICATOR* from the skin if either the skin or the *APPLICATOR* is too hot.
- Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone to clean the skin before using Flash N' Go™.
- Use of Flash N' Go™ may cause temporary pigmentation changes (See "Possible Side Effects" below).

- Keep this device out of the reach of children. Do not use Flash N' Go™ on children or allow children to use it.

2.5. Reducing the risk of injury

As with any electrical device, certain precautions must be taken in order to ensure your safety when using Flash N' Go™.

- **Keep Flash N' Go™ away from water!**

Flash N' Go™ is an electrical device. As such it should always be kept away from water.

Do not place or store Flash N' Go™ where it can fall or be pushed into a tub, sink or any other vessel containing water. Do not place in, or drop into water or any other liquid.

This may cause severe electrocution.

Do not use Flash N' Go™ while bathing.

Do not use Flash N' Go™ if it becomes damp or wet.

Do not reach for Flash N' Go™ if it has fallen into water.

Unplug Flash N' Go™ immediately if it has fallen into water.



Keep Flash N' Go™ away from water!

- **Never open Flash N' Go™!**

Do not attempt to open or repair your Flash N' Go™ device. Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

- Trying to open Flash N' Go™ may also damage the device and will void your warranty.
Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.
- Use Flash N' Go™ only for its intended use and as described in its manual.
- Flash N' Go™ should never be left unattended when plugged into an outlet.
- Do not operate Flash N' Go™ if it has a damaged cord or plug and keep the power cord away from heated surfaces.
- Do not use Flash N' Go™ if you see or smell smoke when it is in use.

- Do not use Flash N' Go™ if it is not working properly or if it appears damaged.
- Do not use Flash N' Go™ if the fan vent in its *APPLICATOR* is cracked, coming off or missing.
- Do not use Flash N' Go™ if the *SKIN COLOR SENSOR* in its *APPLICATOR* is cracked, or broken.
- Do not use Flash N' Go™ if the outer shell is cracked or is coming apart.
- Do not use Flash N' Go™ with a damaged *DISPOSABLE LAMP CARTRIDGE*, or if its *OPTICAL LENS* is cracked, chipped or missing
- Always unplug Flash N' Go™ from the electrical outlet immediately after use.
- Unplug Flash N' Go™ before cleaning.
- Do not use Flash N' Go™ with any attachments or accessories not recommended by Home Skinovations Ltd.

2.6. Possible Side Effects using Flash N' Go™

When used according to the instructions, side effects and complications associated with use of Flash N' Go™ are uncommon. However every cosmetic procedure, including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that can occur with pulsed light hair removal systems designed for home use.

- **Minor Skin Discomfort**
Although home pulsed light hair removal is generally very well-tolerated, most users do feel some mild discomfort during use, usually described as being a mild stinging sensation on the treated skin areas. The stinging sensation usually lasts during the time of the application itself or for a few minutes thereafter. Anything beyond this minor discomfort is abnormal and means that either you should not continue to use Flash N' Go™ because you are unable to tolerate the hair removal application, or that the energy level setting is too high.
- **Skin Redness**
Your skin may become red right after using Flash N' Go™ or within 24 hours of using Flash N' Go™. Redness generally clears up within 24 hours. See your doctor if redness does not go away within 2 to 3 days.
- **Increased Sensitivity of the Skin**
The skin of the treated area is more sensitive so you may encounter dryness or flaking of the skin.
- **Skin Wounds and Burns**
Very rarely, burns or wounds to the skin can occur following the application. The burn or wound can require a few weeks to heal and, extremely rarely, may leave a noticeable permanent scar.
- **Scarring**
Although very rare, permanent scarring may occur. Usually when scarring occurs it is in the form of a flat and white lesion on the skin (hypotrophic). However, it can be large and red (hypertrophic) or large and extend beyond the margins of the injury itself (keloid). Subsequent aesthetic treatments may be required to improve the appearance of the scar.
- **Pigment Changes**

Flash N' Go™ targets the hair shaft, in particular the pigmented cells in the hair follicle and the hair follicle itself. Nevertheless there is risk of temporary hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. This risk of changes in skin pigmentation is higher for people with darker skin tones. Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.

- **Excessive Redness and Swelling**

In rare cases treated skin may become very red and swollen. This is more common in sensitive areas of the body. The redness and swelling should subside within 2 to 7 days and should be treated with frequent applications of ice. Gentle cleansing is OK, but one should avoid exposure to sun.

- **Infection**

Infection of the skin is exceedingly rare but is still a possible risk following a skin burn or wound caused by Flash N' Go™.

- **Bruising**


Very rarely, use of Flash N' Go™ may cause blue-purple bruising that can last 5 to 10 days. As the bruising fades, there may be a rust-brown discoloration of the skin (hyperpigmentation) that can be permanent.

3. Long Term Hair Removal the Flash N' Go™ Way

Flash N' Go™ is a personal light-based device for long-term hair removal. The process of laser and light-based hair removal is well known and established. It has been proven in clinical use around the world for over 15 years as a safe and effective way to achieve long-term hair reduction.

3.1. How does light remove hair?

Light-based hair removal is based on the theory of selective photothermolysis in which optical energy is used to disable hair growth. In order to achieve such thermal effect the hair shaft needs to selectively absorb light energy and transform it into heat. This selectivity is achieved when optical energy that is delivered to the tissue is mostly absorbed by hair shaft pigment, while the skin and the surrounding tissue stays cool. Melanin is the pigment in the hair shaft that is responsible for the absorption of the light, which generates the heat that eventually disables hair growth. Therefore the more melanin present in the hair (i.e. the darker the hair) the more light that can be absorbed and the more effective light can be at removing hair.




Before Flash N' Go™ Application
Go™ Application

During Flash N' Go™ Pulse After Flash N'

3.2. How does the hair growth cycle impact light-based hair removal?


Every hair in our body goes through the three phases of the hair growth cycle: Anagen, Catagen and Telogen. These phases have an important impact on how the process of light-based hair removal works.

Anagen is the hair growth phase while Catagen and Telogen, both are resting phases.



Anagen–Growth phase Catagen–Resting phase Telogen–Resting phase

The time it takes to complete a full hair growth cycle varies from person to person and the location of the hair on the body, but is typically 18-24 months. At any given moment the majority of the hair follicles in any skin area are in the resting phases. These resting hairs cannot be affected by Flash N' Go™. However, hairs in the growing Anagen phase will respond to Flash N' Go™ applications. It is important to understand that it may take a full hair growth cycle to realize complete hair removal results with Flash N' Go™.



3.3. *What is different about Home Pulsed Light™ technology in Flash N' Go™?*

Flash N' Go™ utilizes the patent pending Home Pulsed Light™ technology (HPL™) developed by the scientists at Home Skinovations exclusively for safe and effective removal of unwanted hair at home. The technology allows for Flash N' Go™ to be small, portable and easy to operate compared to other similar systems. And like other light and laser systems HPL™ uses light energy to target melanin pigment in the hair.

3.4. *Plan your Flash N' Go™ hair removal for best results*

A typical full hair growth cycle may take 18-24 months. During this time multiple Flash N' Go™ sessions may be required in order to achieve long term hair removal.

The efficiency of hair removal varies from person to person according to body area, hair color, and how Flash N' Go™ is used.

Typical Flash N' Go™ hair removal plan during a full hair growth cycle:

- The first **3-4 hair removal sessions** with Flash N' Go™ will be approximately **two weeks apart**.
- Hair removal **sessions 5-7** with Flash N' Go™ will be approximately **four weeks apart**.
- After that you will typically use Flash N' Go™ again **from time to time** if and when needed, until long-term results are achieved.

Typical maintenance with Flash N' Go™ after achieving long-term hair removal:

Due to hormonal or other physiological changes dormant hair follicles may become active. Maintenance hair removal sessions with Flash N' Go™ may be required from time to time.


4. First Steps with Flash N' Go™

4.1. Flash N' Go™ device description

Your Flash N' Go™ consists of a **BASE UNIT** and an **APPLICATOR**.

On the **BASE UNIT** you can find the Flash N' Go™ **CONTROL PANEL** including **POWER ON/OFF SWITCH** and **POWER INDICATOR LIGHT**, the **ENERGY LEVEL SETTING BUTTONS** and **ENERGY LEVEL INDICATOR LIGHTS**, the **AUDIO ON/OFF SWITCH**, **READY INDICATOR**, **SYSTEM WARNING INDICATOR LIGHT** and **CARTRIDGE 90% and 100% USAGE INDICATOR LIGHTS**.

On the Flash N' Go™ **APPLICATOR** you can find the **PULSE BUTTON**. The **DISPOSABLE LAMP CARTRIDGE** is located at the **APPLICATOR TIP**.



1. BASE UNIT
2. APPLICATOR
3. PULSE BUTTON
4. DISPOSABLE LAMP CARTRIDGE
5. SKIN COLOR SENSOR
6. CONTROL PANEL
7. POWER ON/OFF SWITCH
8. POWER INDICATOR LIGHT
9. ENERGY LEVEL SETTING BUTTONS
10. ENERGY LEVEL INDICATOR LIGHTS
11. READY INDICATOR LIGHT
12. AUDIO ON/OFF SWITCH
13. CARTRIDGE 90% USAGE INDICATOR LIGHT
14. CARTRIDGE 100% USAGE INDICATOR LIGHT
15. SYSTEM WARNING INDICATOR LIGHT

4.2. Package contents

Upon opening the Flash N' Go™ package, you will find the following parts:

- Flash N' Go™ **BASE UNIT** and **APPLICATOR**
- An AC cord
- A second **DISPOSABLE LAMP CARTRIDGE**
- This User's Manual and a Quick Start Guide leaflet
- An Instructional DVD

4.3. Flash N' Go™ Setup

1. Remove Flash N' Go™ *BASE UNIT*, *APPLICATOR* and other components from box.
 2. Insert the *APPLICATOR* into its cradle in the Flash N' Go™ *BASE UNIT*.
 3. Verify that the *DISPOSABLE LAMP CARTRIDGE* is inserted correctly into the *APPLICATOR*.
 4. Plug the power cord into the Flash N' Go™ *BASE UNIT* socket.
 5. Plug the other end of the power cord into an electrical outlet.
- Your Flash N' Go™ is now ready to start.

4.4. Treating with Flash N' Go™ for the first time

The skin should be **shaved, clean, dry** and **free** of any powders, antiperspirants or deodorants.

1. Press the *POWER ON/OFF SWITCH* to turn Flash N' Go™ on.
2. Shortly after pressing the *POWER ON/OFF SWITCH*, the *READY INDICATOR LIGHT* will turn on, accompanied by a beep, to indicate that your Flash N' Go™ is ready for its first pulse.
3. Your Flash N' Go™ is automatically set to the lowest energy level.

Energy Level

Energy level determines the intensity of the Flash N' Go™ light pulse delivered to your skin, from the lowest level (-) to the highest level (+). Corresponding *INDICATOR LIGHTS* on the *CONTROL PANEL* represent the increases in energy. As energy level increases, so does hair removal results as well as the risk of possible side effects and complications.

Always start your first Flash N' Go™ hair removal session at the lowest energy setting.

If you experience little or no discomfort during and after the hair removal session using the lowest energy level, raise the energy level by one *INDICATOR LIGHT* at the next session, and so on for each subsequent session.



***Do not raise the energy level if you experience abnormal discomfort during or after treating with Flash N' Go™ (See section 2.6 - "Possible Side Effects").
Do not raise the energy level during hair removal session even if you experience no discomfort.
Discomfort may also appear some time after the session.***

Whenever Flash N' Go™ is turned on its energy level will automatically be reset to the lowest energy level. Only one *ENERGY LEVEL INDICATOR LIGHT* will be on.

To set the energy level, press the *ENERGY LEVEL SETTING BUTTONS* using "-" or "+" to respectively decrease or increase the energy level. The number of *ENERGY LEVEL INDICATOR LIGHTS* will coincide with the change in energy level.


For your safety, when using Flash N' Go™ for the very first time, the system is automatically set to deliver the first 50 pulses at the lowest energy level, and the

next 200 pulses at up to level 3.

To cancel these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 2 consecutive "beep" sounds.

To **restore** these safety settings press and hold the "+" and "-" *ENERGY LEVEL SETTING BUTTONS* simultaneously until you hear 3 consecutive "beep" sounds.

4. Using firm pressure, apply the *APPLICATOR TIP* to the skin. Be sure that the skin area you are about to treat is spread evenly and smoothly, and that the *APPLICATOR TIP* is in full contact with the skin.



Firmly press the *APPLICATOR TIP* against the skin


5. Press the *PULSE BUTTON* to emit a pulse. The *SKIN COLOR SENSOR* will briefly measure your skin tone. If the measured skin tone is light enough, the first pulse will immediately be emitted. You will see a flash of light and hear a subtle pop-like sound. This is a normal light pulse of Flash N' Go™. If the measured skin tone is too dark for Flash N' Go application, the 5 green energy level LED will blink for 3 seconds. *SKIN TONE SENSOR* will measure again in your next attempt to emit a pulse.

Note: As this intense pulse may be disturbing to your eyes, it is recommended to close them shortly before pressing the *PULSE BUTTON*.

6. After emitting a pulse, remove the *APPLICATOR TIP* from the skin. Flash N' Go™ will recharge for the next pulse.
7. After approx. 3.5 seconds the *READY INDICATOR LIGHT* will turn on again, accompanied by a beep, to indicate that your Flash N' Go™ is ready for another pulse.
8. Continue to administer pulses until the desired skin area is fully treated.

Covering skin areas

Flash N' Go™ pulses should be administered in rows, starting at one end of each row and progressing sequentially towards the other end. This technique allows better control of skin coverage, and helps you avoid treating the same area more than once or overlapping skin areas



When applied to the skin, the Flash N' Go™ *APPLICATOR TIP* is designed to create temporary pressure marks on the treated area. These visible marks can be used for proper positioning of the next pulse.



Try to avoid overlapping pulses!

Do not treat the same area of skin more than once per hair removal session!



If your skin blisters or burns, STOP USE IMMEDIATELY!

4.5. What to Expect when treating with Flash N' Go™?

For many people, using Flash N' Go™ may be their first experience with a light-based device designed for home use. Flash N' Go™ is simple to use, and hair removal sessions go by quickly. During a Flash N' Go™ session it is **normal** to experience and feel:

- **A Flash of Light** – The bright light of Flash N' Go™ will not harm the eyes when applied to non-facial sites, and special eye-protection is not needed when using Flash N' Go™.
- **A Fan Noise** – The cooling fan in Flash N' Go™ makes noise similar to a hairdryer. This is normal.
- **A Pop Sound with Each Pulse** – When a pulse of HPL™ light is activated, it is normal to hear a subtle pop sound simultaneously with the flash of light.
- **Moderate Pressure of the *APPLICATOR*** – This is necessary and helpful for placement of adjacent pulses of light, and is part of the unique safety feature of Flash N' Go™.
- **A Sensation of Warmth and Tingling** – During each pulse of light it is normal to feel a mild sensation of warmth and tingling from the light energy. Remember it is important to always use low energy settings for initial hair removal sessions. You may feel some warmth for up to an hour after your Flash N' Go™ session.
- **Some Mild Red or Pink Color** – During and just after your Flash N' Go™ session it is not uncommon to see some very mild, pink-like color of the skin. This is usually most noticeable around the hairs themselves. However if you see full redness of the skin, blistering or burns stop use of Flash N' Go™ immediately.

4.6. After treating with Flash N' Go™

- When Flash N' Go™ session has been completed turn Flash N' Go™ off by pressing the *POWER ON/OFF SWITCH*. (Be sure to remember the last energy level setting you used, as it will not be restored when turning Flash N' Go™ on again.)
- Unplug the power cord from the electrical outlet.
- After each hair removal session it is recommended that you clean your Flash N' Go™ device, especially the *APPLICATOR TIP* (See section 5.1: "Cleaning Flash N' Go™").
- After cleaning, it is recommended to store your Flash N' Go™ device in its original box, and keep it away from water.



Skin care following hair removal session

Do not expose treated areas of skin to the sun. Be sure to carefully protect the treated skin with sunscreen, throughout the hair removal period and for at least 2 weeks following the last Flash N' Go™ session.

Side effects and complications

Some patients may experience pigmentation changes resulting from treating with Flash N' Go™. These effects, if they occur, are generally mild and transient.



In case you experience any complication (See section 2.6 Possible Side Effects using Flash N' Go™) ***please contact your physician immediately.***

5. Maintenance of Flash N' Go™

5.1. Cleaning Flash N' Go™

After each hair removal session, it is recommended to clean your Flash N' Go™ device, and especially the *APPLICATOR TIP*.

- Unplug Flash N' Go™ before cleaning.
- Use a dry, clean cloth and a specially formulated cleaner for electronic equipment to gently wipe Flash N' Go™ surface, and especially the *APPLICATOR TIP*.

Never immerse Flash N' Go™ or any of its parts in water!

5.2. Replacing the LAMP CARTRIDGE


The Flash N' Go™ *DISPOSABLE LAMP CARTRIDGE* can fire 750 light pulses that would typically cover 2 Legs, Arms, Under arms, and Bikini line. Pulse intensity is determined only according to the energy level setting of the device. There is no decrease of energy during the usable lifetime of the *LAMP CARTRIDGE*.

When a *DISPOSABLE LAMP CARTRIDGE* has reached 90% of its possible lifetime the *CARTRIDGE 90% USAGE INDICATOR LIGHT* on the *CONTROL PANEL* will be activated, indicating that the *DISPOSABLE LAMP CARTRIDGE* should soon be replaced.

Once all 750 pulses in a *DISPOSABLE LAMP CARTRIDGE* have been used the *CARTRIDGE 100% USAGE INDICATOR LIGHT* will be activated, indicating that further pulses can not be emitted.

In order to continue the hair removal session, it will now be necessary to replace the *DISPOSABLE LAMP CARTRIDGE*.

Note: A *DISPOSABLE LAMP CARTRIDGE* should also be replaced if large spots appear inside it or if its *LENS* is broken.



To replace a *DISPOSABLE LAMP CARTRIDGE*:

1. Press the *POWER ON/OFF SWITCH* to turn Flash N' Go™ off.
2. Unplug the power cord from the electrical outlet.
3. Grasp the used *DISPOSABLE LAMP CARTRIDGE* on both sides, pull it out of the

socket and discard with normal trash.

4. Unwrap a new *DISPOSABLE LAMP CARTRIDGE*.

5. Line up the holes in the *APPLICATOR* socket with the metal and plastic hinges in the new *DISPOSABLE LAMP CARTRIDGE*, and push it gently into place.

If correctly installed, the *DISPOSABLE LAMP CARTRIDGE* will spring back when pushed in and released.

6. Troubleshooting

6.1. "My Flash N' Go™ does not start."

- Make sure the power cord is properly connected to the Flash N' Go™ device.
- Make sure the power cord is plugged into an electrical outlet on the wall.

6.2. "A light pulse is not emitted when I press the PULSE SWITCH."

- Make sure that you have good contact with the skin and that the *APPLICATOR TIP* is evenly and firmly pressed to the skin. For your safety, the *PULSE BUTTON* will activate a pulse only if the *APPLICATOR TIP* is firmly pressed against the skin.
- Skin color sensor stopped the light pulse, your skin may be too dark for using the Flash N' Go. Try to apply Flash N' Go™ again.
- Check the *CARTRIDGE 100% USAGE INDICATOR*. If it is on, disconnect Flash N' Go™ and replace the *DISPOSABLE LAMP CARTRIDGE*.
- Make sure that the *READY INDICATOR* is on.
 - a. If within 10 seconds the *READY INDICATOR* remains off turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - b. If the problem persists, contact your local Flash N' Go™ Customer Service Center.
- Check the *SYSTEM WARNING INDICATOR LIGHT*.
 - c. If the light is on turn Flash N' Go™ off and back on by pressing the *POWER ON/OFF SWITCH* twice.
 - d. If the problem persists, contact your local Flash N' Go™ Customer Service Center.

If these problems persist, contact your local Flash N' Go™ Customer Service Center.



Do not attempt to open or repair your Flash N' Go™ device. Only authorized Flash N' Go™ repair centers are permitted to perform repairs.

Opening Flash N' Go™ may expose you to dangerous electrical components and to pulsed light energy, either of which may cause serious bodily damage and/or permanent eye injury.

Trying to open Flash N' Go™ may also damage the device and will void your warranty.

Please contact Flash N' Go™ Customer Service if you have a broken or damaged device in need of repair.

7. Customer Service

For more information about Flash N' Go™ please enter www.silkn.com.

If your Flash N' Go™ is broken, damaged, in need of repair, or for any other Flash N' Go™ user assistance, please contact Flash N' Go™ Customer Service:

1-877-DO-SILKN / 1-877-367-4556

contact@silkn.com

8. Frequently Asked Questions

1. ***Does Flash N' Go™ really work?***

Yes. In clinical trials held by physicians, Flash N' Go™ was proven to safely achieve long-term hair removal results.

2. ***Where on my body can I use Flash N' Go™?***

The Flash N' Go™ device has been designed for body hair removal anywhere below the neck. The most common areas treated with Flash N' Go™ are: legs, underarms, arms and bikini line. It is not recommended to use Flash N' Go™ on the face. Flash N' Go™ is not suitable for everyone. We recommend that you read all the Flash N' Go™ warnings and contraindications in this User's Manual.

3. ***How long does a Flash N' Go™ treatment session take?***

The time can vary depending on the area of the body treated. A full leg can take up to 30 minutes, or two underarms could take 10 minutes. Because Flash N' Go™ runs on regular electric power it can be used for as long as needed to complete a full hair removal session of the desired body part(s).

4. ***Is Flash N' Go™ safe?***

Flash N' Go™ has been designed with your safety in mind, and tested and approved by top dermatologists and plastic surgeons to meet their safety standards for a home-use device. But like any skin product or electronic device, one must use according to the operating instructions and user precautions.

5. ***Will Flash N' Go™ hurt?***

When used properly most users of Flash N' Go™ report feeling a slight sensation of heat at the time of the pulse of light. Users with thicker and darker hairs may feel slightly more discomfort, but this discomfort subsides once the hair removal session is completed. For your convenience Flash N' Go™ has five setting levels that can be used according to your sensitivity.

6. ***How often should I use Flash N' Go™?***

Hair removal sessions with Flash N' Go™ should be spaced every two weeks for the first three to four sessions. After that hair removal sessions should be done if hairs have grown back, until the desired results are achieved.

7. ***Is Flash N' Go™ effective on white, grey or blonde hairs?***

Flash N' Go™ works best on darker hair types, or hair that contains more melanin. Melanin is the pigment that gives hair and skin their color, and will absorb light energy. Black and dark brown respond the best, although brown and light brown hairs will also respond but typically require more hair removal sessions. Red may show some response. White, grey or blonde hairs usually don't respond to Flash N' Go™ though some users have noted results after multiple hair removal sessions.

8. ***Can I use Flash N' Go™ on brown or black skin***

Do not use Flash N' Go™ on naturally dark skin complexion! Flash N' Go™ removes unwanted hair by selectively addressing hair pigment. Varied quantities of pigment also exist in the surrounding tissue of skin. The quantity of pigment in a particular person's skin, which is manifested by their skin complexion, determines the degree of risk they are exposed to using Flash N' Go™. Treating dark skin with Flash N' Go™ can result in adverse effects such as burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

9. ***When will I see results from Flash N' Go™?***

As with any light-based or laser hair removal device, results are not immediate, and in fact you may not think anything happened at all. Hair may sometimes appear to be growing back after a hair removal session, but typically after two weeks many of these hairs will simply fall out. Additionally hairs grow in three different stages and only hairs in an active growth stage will be affected by Flash N' Go™. This is one of the main reasons that multiple sessions are required to achieve the desired result.

10. **Can a man use Flash N' Go™?**

Though designed exclusively for women, Flash N' Go™ may be suitable for use by men. However, hairs on men, typically those on the chest, will require more hair removal sessions than that of women to get the desired results. As with women, using Flash N' Go™ on men's facial hair such as beard and mustache, is not recommended.

11. **Why is my hair growing, even though I treated it a week ago?**

It is quite common for hair to appear as if it is still growing up to two weeks after a hair removal session with Flash N' Go™. This process is known as "ejection" and at around two-weeks you'll see that these hairs simply fall out or slide out with a slight tug. (We don't however recommend pulling on the hairs – just let them come out naturally.) It is also possible that some hairs, due to missed application or different stages of growth, were not affected by Flash N' Go™. These hairs will be treated in follow-up sessions, and hence the reason multiple hair removal sessions are needed to get the best result with Flash N' Go™.

12. **I've heard that some hairs grow back lighter and finer after light based hair removal?**

This phenomenon is well documented amongst aestheticians and doctors using light and laser devices for hair removal. It is possible that some hairs will grow back lighter and finer after hair removal with Flash N' Go™. Usually these hairs are a fraction of what was originally there, and continued use may have a desirable effect on them.

13. **Why can't I treat myself if I have an "active" suntan?**

Do not use Flash N' Go™ on tanned skin or after sun exposure! Tanned skin and particularly following sun exposure, contains large quantities of the pigment Melanin. This applies to all skin types and complexions, including those which don't seem to tan quickly. The presence of large quantities of Melanin exposes the skin to higher risk of adverse effects when using Flash N' Go™ including burns, blisters, and skin color changes (hyper- or hypo-pigmentation).

14. **Is long-term use of Flash N' Go™ dangerous for my skin?**

The use of light and laser energy in aesthetic medicine has been well documented for over 15 years in professional peer-reviewed journals, and by well respected institutions like the Mayo Clinic. These journals and institutions have not reported any side-effects or damage from long-term use of light and laser device.

15. **Can I use Flash N' Go™ to remove my chin hair or elsewhere on my face?**

It is not recommended to use Flash N' Go™ on the face or neck. Unlike body hair, most of women's potential facial hair remains inactive and hidden, while only a fraction of it grows and surfaces. Exposing facial hair to the light pulses of Flash N' Go™ may remove apparent hair but simultaneously stimulate unwanted growth of hidden hair. In addition, using Flash N' Go™ to remove facial hair may cause serious eye injury.

16. **How long should I wait to treat with Flash N' Go™ after unprotected exposure to the sun?**

One should wait 4 weeks before using Flash N' Go™ after unprotected exposure to the sun. However, if there is ever any uncertainty about sun exposure please contact your physician or Home Skinovations customer support.

17. **Should I do anything before using Flash N' Go™?**

Before any Flash N' Go™ session it is important to avoid sun exposure on the treated area for at least four weeks. A high level UV Sun Screen (SPF 50+) will help, as will clothing covering the treated area. The area to be treated should also be cleaned with mild soap and water, and the hairs shaved down to skin level.

18. **How should I care for the treatment area after using Flash N' Go™?**

The area treated with Flash N' Go™ can be cleaned and maintained with standard skin care products. Special care must be taken to avoid unprotected sun exposure. Strong sunscreens (50+ SPF) and covering clothing are suitable for protection from the sun.

19. **Should I pull the hairs out after treatment?**

No, let the hairs gradually fall out on their own. This may take up to 2 weeks.

9. Specifications

| | |
|---------------------------------|--|
| Spot size | 2cm x 3cm [6cm ²] |
| Speed | 1 pulse every 3.5 second: 1.7 cm ² /sec |
| Technology | New Home Pulsed Light™ |
| Max Energy Level | Max 5J/cm ² |
| Wavelength | 475-1200nm |
| Charging time / Power source | Continuous operation |
| Electrical requirements | 100-240VAC, 2A |
| Time needed to treat lower legs | 30 minutes |
| Operation and Safety | SKIN COLOR SENSOR seamlessly ensures use only on appropriate skin types. Safety sensor Tip enables self-trigger for maximum safety & control. Easy to use on all body parts. |
| package size | Height 9 inch, Width 9 inch, Depth 5.2 inch |
| System weight | 4 pounds |
| | |

10. Labeling



Class II equipment

Degree of protection against electric shock: type BF applied part

WEEE - Waste Electrical and Electronic Equipment

CE Mark

Follow operating instructions

This device comply with part 15 of the (FCC) Federal Communications Commission.

CSA Mark for USA and Canada

Degree of protection against ingress water: ordinary

This device is not suitable for use in the presence of flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.



HomeSkinovations

EC Declaration of Conformity

We Home Skinovations Ltd. declare under our sole responsibility that the home Skin treatment products

Brand Name: **Home Skinovations**
Model Name: **Flash N' Go, Disposable Lamp Cartridge**
Models Number: **AS100001E, AS100674B, AS100756A**

are fully in conformity with the essential requirements of Council Directive 89/336/EEC and 73/23/EEC amended by Directives 93/68/EEC.

This declaration is based on the full compliance of the products with the following European standards:

| | |
|--|--|
| EMC: | EN 61000-3-2:2000; (A2:2005), EN 61000-3-3:1995; (A1:2001) EN 55014-1:2000; (A2:2002), EN 55014-2:1997; (A1:2001) |
| Electrical safety: | EN 60335-1:2002; (A1:2004; A12:2006), EN 60335-2-23:2003, EN 60335-2-27:2003 |
| Human Exposure to electromagnetic field radiated | EN 50364:2001 |

| | |
|---|--|
| Manufacturer: | Authorized Representative: |
| Home Skinovations Ltd. Apolo Building, Shaar Yokneam, POB 533, Yokneam 20692, ISRAEL Tel: +972-4-9097470 Fax: +972-4-9097471 | Home Skinovations GMBH Dr. Kurt Huber Str. 6 D-82031 Grünwald, GERMANY Tel. +49-89-64919530 Fax +49-89-64919531 |

Signed:

Signed:

Gabi Lavi
General Manager
Yokneam, February 4, 2008

Dr. Amir Waldman
VP Regulatory & Clinical Affairs
Yokneam, February 4, 2008


For Information: On the basis of this declaration, these products and packaging will bear the following mark:



287

Flash N' Go™ is a **physician recommended** and innovative light-based device for hair removal in the **privacy of your home**. **Safe, easy to use** and **cost effective**, Flash N' Go™ achieves **marked clinical results** while offering personal convenience.

Flash N' Go™ utilizes the highly sophisticated and new **Home Pulsed Light™ (HPL™)** technology developed by **Home Skinovations Ltd.**



Home Skinovations Ltd.
Apolo Building
Soltam Industrial Zone.
Yokneam, 20692, Israel

Flash N' Go - Risk Analysis Report

DO1008941

289

| Approved By | Name | Signature | Date |
|-------------|------|-----------|------|
| (b) (4) | | | |
| | | | |

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Home Skinovations


Revisions

| Rev. | Description | Approval | Signature | Date |
|------|-------------|----------|-----------|------|
| | | 46 | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

290

Contents

| | | |
|------|----------------------------|----|
| 1 | Introduction | 5 |
| 1.1 | Purpose | 5 |
| 1.2 | Scope | 5 |
| 1.3 | Overview | 5 |
| 1.4 | Applicable Standards | 5 |
| 1.5 | Product Short Description | 5 |
| 1.6 | Core team | 5 |
| 2 | Risk Method | 6 |
| 3 | Risk Analysis | 8 |
| 3.1 | Chemical Hazards | 8 |
| 3.2 | Mechanical Hazards | 9 |
| 3.3 | Fire Hazards | 10 |
| 3.4 | Environmental Hazards | 12 |
| 3.5 | Biological Hazard | 13 |
| 3.6 | Electrical Hazards | 14 |
| 3.7 | Electromagnetic Hazards | 15 |
| 3.8 | Software Hazards | 16 |
| 3.9 | System Malfunction Hazards | 18 |
| 3.10 | Clinical Hazards | 20 |
| 3.11 | Other Hazards | 22 |



Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

| | | |
|------|--|----|
| 3.12 | Specific Hazards of EN60601-1, Sub-Clause 52.2; 52.4.1 prohibits | 23 |
| 4 | Information supplied by the manufacturer | 24 |
| 4.1 | Contraindications | 24 |
| 4.2 | Precautions: (treat with caution if patient has any of following risk factors) | 25 |
| 4.3 | Possible Side Effects | 25 |
| 5 | Risk Estimation of Hazards Associated to the Flash N' Go | 26 |
| 5.1 | CONCLUSION | 26 |

272

1 Introduction

1.1 Purpose

This document is used to evaluate the risk associated with the *Flash N' Go* system.

1.2 Scope

This Hazard Analysis report applies to the Flash N' Go system.

1.3 Overview

This document describes the safety hazards associated with the Flash N' Go. It also describes the possible consequences of these hazards, their probability of occurrence, the Risk Level, how they will be mitigated and how and where that mitigation will be implemented.

1.4 Applicable Standards

ISO14971: Dec 2000
IEC 60601-1:1998
EN 60601-1-4

1.5 Product Short Description

The Flash N' Go is intended to be used for home use hair removal applications. It's based on Optical radiation.

1.6 Core team

Dr. Amir Waldman
Gabi Lavi
Emil
Suhair Francis-Najjar

VP QA & Regulation
General Manager
Software Engineer
Software Consultants

293

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel



2 Risk Method

Levels of Risk:

| Rating | Risk Level Definitions |
|------------|--|
| Negligible | Nuisance only / Cosmetic in nature / User Inconvenience / Will not cause personal injury or harm |
| Minor | Loss of efficacy but patient safety not compromised / Will not cause personal injury or harm |
| Moderate | Non-functional or poor functional performance of device requiring additional medical intervention to compensate / Injury to patient or user requiring no or minor medical intervention to compensate |
| Major | Incapacitating / Serious injury or death |

Methods of Control:

| Level | Methods of control |
|-------|---------------------------|
| 1 | Inherent safe design |
| 2 | Manufacturing control |
| 3 | Adequate USER information |

294

Home Skinovations Ltd.
 Apolo Building.
 Soltam Industrial Zone.
 Yokneam, 20692, Israel

HomeSkinovations

Risk levels:

| | Negligible(1) | Minor(2) | Moderate(3) | Major (4) |
|-------------------------------------|---------------|----------|-------------|-----------|
| Adequate USER information (3) | IV | III | II | I |
| Manufacturing control (2) | IV | III | III | II |
| Inherent safe design (1) | IV | IV | IV | III |

RISK IV: Acceptable Risk
 RISK III: Medium risk
 RISK II: High Risk.
 RISK I: Very High Risk, Intolerable risk

Mitigation Strategy (methods of Control):

| | Mitigation Strategy | Description |
|-----|--|--|
| (1) | Inherent Safety by Design | This strategy includes system functionality and features that have been especially incorporated to alleviate safety hazards. |
| (2) | Manufacturing Control, ATP, Preventive Maintenance | This includes protective measures in the device itself/ in the manufacturing process, and testing |
| (3) | Information needed for | This includes the provision of adequate user information, training, and labeling |

295

Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

HomeSkinovations

3 Risk Analysis

3.1 Chemical Hazards

| Cause Of the Hazard | | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---------------------|--|---------------|-------------------|------------|--|-------------------|----------------|-----|
| 3.1.1 | Hazard due to the presence of dangerous gases. | Moderate | (1)(3) | IV | No gases are used, or generated during normal operation | Minor | IV | |
| 3.1.2 | Hazard due to the presence of dangerous liquids. | Moderate | (1)(3) | IV | The only liquid is sealed in approved electrical capacitors | Minor | IV | |
| 3.1.3 | Hazards of Implosion and Explosion (Electrolytic Capacitors) | Moderate | (1) | II | The Electrolytic Capacitors are the only components thought to be at risk of exploding. The Flash N' Go enclosure is believed to be stronger than needed to contain any ejected parts | Negligible | IV | |
| 3.1.4 | Hazard due to the use of solid materials | Moderate | (1) | II | All solid materials such as metals, plastics, terminal strips, PCB etc. are common materials are proven safe by long use or approved by UL (U94-0) No materials known to generate dangerous vapours or gases when heated or burned. | Negligible | IV | |

296

3.2 Mechanical Hazards

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---|---------------|-------------------|------------|--|-------------------|----------------|-----|
| 3.2.1 Sharp edge | Minor | (1) | III | Mechanical hazards have been avoided and minimized by construction and good engineering practice and testing according to EN60601-1. | Negligible | IV | |
| 3.2.2 Accidental mechanical damage | Minor | | III | | Minor | III | |
| 3.2.3 Moving parts (The only used moving part are Fans) | Minor | (1)(3) | | Fan is protected with grid. User not allowed opening covers. | Negligible | IV | |
| 3.2.4 Vibration | Minor | (1) | III | Vibration absorbed by mounts. | Negligible | IV | |
| 3.2.5 Mechanical shock | Minor | (1) | III | The unit complies with EN60601-1 Standard Clause 24 requirements | Negligible | IV | |

297

3.3 Fire Hazards

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|-----------------------------------|-----------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.3.1 Electrical Fire Hazards | <i>Moderate</i> | (1)(2) | II | Electrical fire hazards are avoided and minimized by construction and type testing according EN60601-1, and good engineering practice such as: Proper selection fuses, thermal fuses, type and wires size (UL listed). Operating all electrical components within their ratings. Protection of wires from sharp edges. All electrical components (IC, Capacitors, Transformers, resistor, etc) were choosing so it operates below its rated (Temperature, Current, Voltage, etc). To verify this Home Skinovations engineers has selected the components that are hot, and tested them with thermocouple during testing of design verification. | <i>Minor</i> | III | |
| 3.3.2 Ignition of thermo-plastics | <i>Moderate</i> | (1)(2) | II | All thermo-plastic materials used in unit have the adequate flame protection according to UL94 | <i>Minor</i> | III | |

298

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|--|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.3.3 Ignition of Transformer Windings | Moderate | (1) | II | Use of a thermal switch in the transformer | Minor | III | |
| 3.3.4 Other Fire Hazards, Working beam (inside the applicator) | Moderate | (1) | II | If something flammable intrudes into the air paths of the working beam, it could be ignited. The probability for it to happen is very small because of two reasons: 1) The air paths is very narrow, 2) The air paths is covered by plastic parts and it unlikely that something will intrudes into the air paths. Never the less if something can be ignited it's the parts are self-extinguishing. | Minor | III | |
| 3.3.5 Fires when the device is not in use | Moderate | (1)(3) | IV | The user manual instruct the user to shut down the system when it not in used, and to unplug the line cord, in this position it is believed to be impossible for fire to start in the Flash N' Go, if the Flash N' Go is plugged in but the switch is off, it's unlikely that fire can be start. | Minor | IV | |

299

3.4 Environmental Hazards

| | Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|-------|--|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.4.1 | Storage or operation outside prescribed environmental conditions | Minor | (3) | IV | System to be used in home under normal room environmental condition | Negligible | IV | |
| 3.4.2 | Contamination due to waste procedure and/or disposal | Minor | (3) | IV | Cleaning procedure are defined in the User Manual prior to using the system | Negligible | IV | |
| 3.4.3 | Inadequate supply of power Restriction of cooling | Minor | (1)(2) | III | Self test and status polling detects power instability | Negligible | IV | |
| 3.4.4 | Incompatibility with other devices (by radiation and conducting lines) | Minor | (1) | III | Tested for compatibility with RFI interference, Tested were done per EN 60601-1-2 | Negligible | IV | |

3.5 Biological Hazard

| | Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|-------|---|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.5.1 | Bio-burden/ contamination (applicator parts) | Moderate | (1) | II | The parts in contact with the patient comply with bio-compatibility requirements of EN 30993-1 Standard | Negligible | IV | |
| 3.5.2 | Bio incompatibility (and Inability to maintain hygienic safety) | Moderate | (1)(3) | IV | | Negligible | IV | |
| 3.5.3 | toxicity | Moderate | (1) | II | The parts in contact with the user comply with bio-compatibility requirements of EN 30993-1 Standard | Negligible | IV | |
| 3.5.4 | allergenicity | Moderate | (1) | II | No metal parts are in contact with the user | Negligible | IV | |
| 3.5.5 | (Cross) Infection | Moderate | (1)(3) | IV | | Moderated | IV | |

3.6 Electrical Hazards

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.6.1 Electrical Shock Hazards in normal use. | Major | (1)(2) | I | Shock hazards are avoided and minimized by construction and production testing in accordance with EN60601-1, good engineering and good manufacturing practices. | Minor | III | |
| 3.6.2 Unstable electrical supply | Minor | (1) | III | Self-test is done after power electronics is initiated. Power on does not turned on automatically, system is turned on in Stand By mode Pooling of electrical parameters every 100[ms] | Negligible | IV | |
| 3.6.3 Heat due electrical short | Minor | (1) | III | See "Electrical Fire Hazards" | Negligible | IV | |
| 3.6.4 Mechanical force causing electrical short | Minor | (1) | III | Electrical short by mechanical force are avoided and minimized by construction and production testing in accordance with EN60601-1, good engineering and good manufacturing practices | Negligible | IV | |
| 3.6.5 Electrostatic phenomena (Electrostatic charge build-up) | Minor | (1) | III | Flash N' Go is complies with electrostatic discharge requirements of EN 60601-1-2 Standard | Negligible | IV | |

302

3.7 Electromagnetic Hazards

| Cause Of the Hazard | | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---------------------|---|-----------------|-------------------|------------|--|-------------------|----------------|-----|
| 3.7.1 | Emission s of electromagnetic interference | <i>Moderate</i> | (1) | II | The Flash N' Go is designed, built and tested to conform to: 1) EN 60601-1-2 (EN55011 Conducted) 2) EN 60601-1-2 (EN55011 Radiated) 3) CFR 47 Part 15 Class B | <i>Minor</i> | III | |
| 3.7.2 | Susceptibility to electromagnetic interference | <i>Moderate</i> | (1) | II | The Flash N' Go is designed, built and tested to conform to: 4) EN 61000-4-2 (ESD) 5) EN 61000-4-4 (Radiated Susceptibility) 6) EN 61000-4-4 (Electrical fast Transients) 7) EN 61000-4-5 (Surge) 8) EN 61000-4-5 (Conducted RF) 9) EN 61000-4-8 (Power frequency Magnetic Field) 10) EN 61000-4-11 (Voltage dips and interrupts) | <i>Minor</i> | III | |
| 3.7.3 | Human exposure to electromagnetic fields radiated | <i>Moderate</i> | (1) | II | The Flash N' Go is designed, built and tested to conform to: EN 50364:2001 | <i>Minor</i> | III | |

3.8 Software Hazards

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|--|---------------|-------------------|------------|---|-------------------|----------------|---------------------|
| 3.8.1 Firmware Code corruption (total code corruption). | Minor | (1) | III | MCU load firmware software; but watchdog reset the MCU. | Minor | III | SW V&V 9.4.1 |
| 3.8.2 Firmware Code corruption (errors in a few bytes) | Minor | (1) | III | MCU load firmware but software checks CRC and error led turn on. | Minor | III | SW V&V 9.2.1, 9.2.2 |
| 3.8.3 Omission – the system produce no flow at a time when one would be expected | Minor | (1) | III | Hardware watchdog expires, the watchdog de-energizes the system. | Minor | III | SW V&V 9.4.4 |
| 3.8.4 Commission – the system produces a flow when a perfect system would not have | Minor | (1) | III | Watchdog uses logical program flow monitoring, the watchdog de-energizes the system. | Minor | III | SW V&V 9.4.4 |
| 3.8.5 The software fire pulses before attachment to the skin | Moderate | (1) | IV | Hardware control pulse activation only when three conditions are met: 1. HP push button. 2. Cartridge switch is pressed. 3. A software pulse is enabled. | Minor | IV | SW V&V 9.4.2 |
| 3.8.6 Buttons are not functioning | Negligible | (1) | IV | Design reviews ; and production test | Negligible | IV | |
| 3.8.7 Error in downloading the software | Minor | (1) | III | Verification test in production | Minor | III | SW V&V 9.1.5 |

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---|---------------|-------------------|------------|---|-------------------|----------------|-------------------------------|
| 3.8.8 Software failure in Monitoring the HW status | Minor | (1) | III | Continuous monitoring of the correct behaviour of the hardware, by high rate scanning the signals of all the sensors; detection of errors as soon as they occur and shutting down the system. | Minor | III | SW V&V 9.3.1, 9.3.2, 9.3.3 |
| 3.8.9 User Confusion on the status of the system | Minor | (1) | III | Error indication shall have explanation in the operation manual. | Minor | III | Operation manual |
| 3.8.10 Software does not indicate if the Tip was replaced during treatment. | Minor | (1),(3) | IV | Software routinely displays indication if Tip is disconnected Tip removal is forbidden during system power on. | Minor | IV | SW V&V 9.4.3 |
| 3.8.11 Late – a flow is produced after it should | Minor | (1) | III | Hardware watchdog expires, the watchdog de-energizes the system. | Minor | III | Watchdog |
| 3.8.12 Value (d) – Reading data with error from RFID Tag | Minor | (1) | III | CRC protocol and checksum to check communication error | Minor | III | SW V&V 9.1.6 |
| 3.8.13 Value (u) – a flow contains a value that is incorrect, but in a way which the recipient cannot detect. | Minor | (1) | III | Code review for error handling strategy in software. | Minor | III | Code review |

25

3.9 System Malfunction Hazards

| | Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|-------|--|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.9.1 | Power delivery is not as displayed: a. Damage of light source. b. Line voltage problems, c. Software fault. | Moderate | (1)(2) | III | | Moderate | III | |
| 3.9.2 | Optical power is higher by 20% of nominal. | Moderate | (1)(2) | III | Capacitor voltage are monitored by the system, and if fault are detected the system stops charger power supply | Minor | III | |
| 3.9.3 | Acoustic pressure | Moderate | (1) | II | Closed console design, rubber buttons for vibrating modules. | Negligible | IV | |
| 3.9.4 | Inadequate performance characteristics for the intended use | Moderate | (1)(3) | II | Applicator lifetime is limited. | Moderate | IV | |
| 3.9.5 | System is not working as desired (Specifications do not meet the application.). | Moderate | (3) | IV | Specification is validated with the clinical data as an on going process. | Minor | IV | |
| 3.9.6 | Loss of mechanical integrity (poor mechanical design or lack of QC) | Minor | (1) | II | The machine was designed by competent mechanical engineers and manufactured by high quality shops. QC checks each module and the whole machine in the critical points of the manufacturing process. | Negligible | IV | |

306

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|------------------------------|-----------------|-------------------|------------|--|-------------------|----------------|-----|
| 3.9.7 Inadequate packaging | <i>Minor</i> | (1) | III | Packaging procedures were written and controlled by the QC personnel. All packaging are specially treated and disinfected against termites, fungus and fire according to the Israeli standard IS262 | <i>Negligible</i> | IV | |
| 3.9.8 Liquid spillage | <i>Minor</i> | (1) | III | The upper side of the system is sealed to prevent spillage of liquids. The enclosure provides a protection of IP20 according to IEC529, degrees of protection provided by enclosures | <i>Negligible</i> | IV | |
| 3.9.9 Skin color malfunction | <i>Moderate</i> | (1) | IV | Lading including user manual and instructional DVD are clearly stating that the device should be used on light skin completion, (skin color chart is included) | <i>Minor</i> | IV | |

3.10 Clinical Hazards

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|-----------------------|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.10.1 Tissue Hazards | Moderate | (1)(3) | II | Flash N' Go is intended to be used non-invasively, as hair removal. The hazard of cutting major blood vessels, or major nerves, or perforating bowel which are present for cutting instruments, are not presents. | Minor | IV | |
| | | | | (b) (4) | | | |

308

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|--|---------------|-------------------|------------|---|-------------------|----------------|-----|
| | | | | c. The Flash N' Go emits a visual warning that it is in Ready mode. (b) (4) | | | |
| 3.10.2 Eye Hazards (Aiming directly into the eye) | Major | (1)(3) | II | | Minor | III | |
| 3.10.3 Excessive or Insufficient Output Power from pulsed light applicator | Moderate | (1)(3) | II | Each applicator is assembled by imbedded power meter. the light output power is adjusted to be $\pm 10\%$ of the correct value. | Minor | IV | |
| 3.10.4 Treatment is unintentionally interrupted (Unstable electrical supply) | Moderate | (1)(3) | II | Self-test is done after power electronics is initiated. Power does not turn on automatically, system is turned on in Stand By mode (b) (4) | Minor | IV | |

3.11 Other Hazards

| Cause Of the Hazard | | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---------------------|-------------------------|---------------|-------------------|------------|---|-------------------|----------------|-----|
| 3.11.1 | Human-Machine Interface | Moderate | (3) | IV | Considerable thought has been devoted to making the interface natural and simple. The read-out device is a set of LEDs, with good brightness. In addition, an audible emission tone alerts the user of an emission in process. All Human-machine requirements of EN60601-1 have been met. | Negligible | IV | |

3.12 Specific Hazards of EN60601-1, Sub-Clause 52.2; 52.4.1 prohibits

| Cause Of the Hazard | Level of Risk | Method of Control | Risk Level | Risk Reduction Method | New Level of Risk | New Risk Level | V&V |
|---|---------------|-------------------|------------|--|-------------------|----------------|-----|
| 3.12.1 Emission of flames, molten metal, and poisonous or ignitable gas | Moderate | (1)(2) | III | The Flash N' Go will not do this, These matters have been addressed in Fire Hazards , and they were the subjects of various abnormal tests during the EN60601-1 Safety Tests. | Negligible | IV | |
| 3.12.2 Excessive deformation of covers | Moderate | (1)(2) | III | The strength of the Flash N' Go covers was evaluated by the ITL engineers for EN60601-1 | Negligible | IV | |
| 3.12.3 Excessive temperatures | Moderate | (1)(2) | III | The temperature was measured during Normal and Abnormal tests made by Home Skinovations and ITL inspecting engineer for the EN60601-1 safety tests. | Negligible | IV | |
| 3.12.4 System is used not according inclusion criteria | Moderate | (3) | II | Intended use is described in user manual.. | Negligible | | |
| 3.12.5 System is used in contraindicated cases | Major | (3) | I | Contraindications are described in user manual. | Moderate | IV | |
| 3.12.6 Wrong treatment parameters are used | Moderate | (3) | II | Energy produced by the system is limited and can not exceed the preset significantly. User manual describe the treatment procedure to avoid significant damage | Negligible | IV | |

3/1

4 Information supplied by the manufacturer

| Hazard Event | Risk Reduction Method | Risk Level |
|---|---|-------------------|
| Manufacturer not identified | Manufacturer Identified in the User Manual | <i>Negligible</i> |
| Lack of safety information in the instructions manual | Safety instructions are defined in the User Manual | <i>Negligible</i> |
| Use of symbols non harmonized | Use of symbol according to EN980 | <i>Negligible</i> |
| Intended use not obvious for user | Intended use is defined in the User Manual | <i>Negligible</i> |
| Lack of check-list following installation | User Manual | <i>Negligible</i> |
| Lack of preventive maintenance instructions to keep the device safe | Preventive maintenance instruction are defined in the User Manual | <i>Negligible</i> |

4.1 Contraindications

- Pacemaker or internal defibrillator or metallic implants in the treatment area.
- Current or History of cancer, especially skin cancer, or pre-malignant moles.
- Pregnancy and nursing.
- Any active infection in the treatment area, such as sores, psoriasis, eczema and rash as well as very dry and fragile skin.
- Diseases which may be stimulated by light at the wavelengths used, such as AIDS and HIV, or use of immunosuppressive medications.
- Poorly controlled endocrine disorders such as diabetes.
- History of bleeding coagulopathies, or use of anticoagulants.
- Excessively tanned skin from sun, sun-beds or tanning creams.

4.2 Precautions: (treat with caution if patient has any of following risk factors)

- Use of medication and herbs known to induce photosensitivity to light exposure at the wavelengths used, such as Isotretinoin (Accutane), tetracyclines, or St. John's Wort.
- Medical Conditions that impair healing.
- In case of preventive Aspirin, treat with care or stop a few days before and after treatment.

4.3 Possible Side Effects

- Pain
- Skin redness (Erythema)
- Damage to natural skin texture (crust, blister, burn)
- Change of pigmentation (hyper- and hypo-pigmentation)
- Scarring
- Excessive Swelling (Edema)
- Fragile skin
- Bruising

5 Risk Estimation of Hazards Associated to the Flash N' Go

Through hazard analysis, we determined that electric shock and eye/tissue burns and system used in contraindicated cases are the three significant hazards present with the Flash N' Go. Neither of the hazards identified has been determined to be unacceptably high.

The risk of electric shock hazard is considered to be LOW. It is mostly limited to Service personnel because the only hazardous live parts are contained within the Flash N' Go console; which would require removal of the side cover while the system is turned on at the Key Switch.

The risks of eye/tissue burn hazards are considered to be MODERATE because the user could be harmed by such a hazard if they do not take adequate measures to protect against said hazards. Adequate measures are defined in the User Manual Safety section.

5.1 CONCLUSION

No known hazards have been identified. The hazard prevention features does not limit the product efficiency and defined specification

Flash N' Go - Software System Test Plan and Results

DO1008951

| Role | Name | Signature | Date |
|---------|------|-----------|------|
| (b) (4) | | | |
| | | | |
| | | | |



Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

Revisions

| Rev. | Page | Description | Approval | Date |
|---------|------|-------------|----------|------|
| (b) (4) | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Table of Contents

| | |
|---|-----------|
| 1. INTRODUCTION | 4 |
| 1.1 PURPOSE | 4 |
| 1.2 SCOPE | 4 |
| 1.3 DEFINITIONS..... | 4 |
| 1.4 REFERENCES | 4 |
| 2. TEST PLAN TYPE & ATTRIBUTE REFINEMENT | 5 |
| 2.1 TEST PLAN TYPE | 5 |
| 2.2 LEVEL | 5 |
| 2.3 OBJECTIVE | 5 |
| 2.4 APPROACH | 5 |
| 2.5 READINESS CRITERIA | 6 |
| 2.6 INPUT WORK PRODUCTS | 6 |
| 2.7 TEST ENVIRONMENT | 6 |
| 2.8 ROLES & RESPONSIBILITIES..... | 6 |
| 2.9 FORMAL REVIEWING | 8 |
| 3. TEST ITEMS..... | 8 |
| 4. FEATURES TO BE TESTED | 8 |
| 5. FEATURES NOT TO BE TESTED | 8 |
| 6. TESTER SKILLS | 9 |
| 7. SAFETY ISSUES | 9 |
| 8. TEST TERMINATION CRITERIA..... | 9 |
| 9. TEST PROCEDURES | 10 |
| 9.1 VERIFY INTERFACES..... | 10 |
| 9.2 VERIFY FUNCTIONALITY | 12 |
| 9.3 VERIFY SYSTEM RECOVERY..... | 19 |
| 9.4 VERIFY SAFETY..... | 24 |
| 10. TEST PLAN RESULTS SUMMARY..... | 26 |

1. Introduction

1.1 Purpose

A system test plan for the Flash N' Go should support the following:

To detail the activities required to prepare for and conduct the system test.

To communicate to all responsible parties the tasks that they are to perform and schedule to be followed in performing the tasks.

To define the source of information used to prepare the plan.

To define the test tools and environment needed to conduct the system test

1.2 Scope

This test plan covers a full systems test of the Flash N' Go Device. This includes operate and user procedures, as well as programs. In addition to comprehensively testing firmware functionality, hardware interfaces, performance, load test, download procedure, integrity, recovery and usability.

Ideal system testing presumes that all components have been previously, successfully, integrated.

1.3 Definitions

| Term | Definition |
|------|------------------|
| STP | System Test Plan |
| | |
| | |

1.4 References

| Document | Revision | Location |
|---|----------|----------|
| Software requirement specification | | |
| Software verification and validation plan | | |
| Software test strategy | | |
| | | |

2. Test Plan Type & Attribute Refinement

2.1 Test Plan Type

This test plan is functional, stress, performance, safety, testing used to verify the Flash N' Go SW release.

2.2 Level

System Level

2.3 Objective

The objectives of the system test plan are:

1. Functional testing: testing based on a comparison of a component's actual behavior with its specified behavior without regard to the component's implementation.
2. Performance testing: tests that explore the consistency of a system with its performance objectives such as response time or throughput
3. Stress testing - testing in which a system is subjected to unrealistically harsh inputs or load with inadequate resources with the intention of breaking it
4. Recovery testing - testing that explores how well the system meets its objectives in recovering from programming errors, hardware failures, and data errors.
5. Safety testing – aims at showing that an unsafe state cannot be reached although that state could be invalid (traceable to hazard analysis document)
6. Regression testing – it assumed that several iterations of the system test will be done in order to test program modifications made during the system period. A regression test will be performed for each new version of the system to detect unexpected impact resulting from program modifications. The regression test will be done by running all of the tests on a new version that were run on the previous version and then comparing the results.

2.4 Approach

Testing methods used to achieve the objective of this test plan:

Black box technique: Data-Flow testing.

Black box technique: Domain testing & boundary value analysis.

The developer will assist in developing the test cases; this will help ensure that the tests represent the production use of the system.

2.5 Readiness Criteria

The Unit and Integration testing shall have been passed and completed successfully.

The software is stable (means a reliable software with no major issues).

The SW Packages under test have been collected under a system test build in configuration management.

All test Hardware and environments must be in place, and free for system test use.

ALL Human resources must be assigned and in place

2.6 Input Work Products

The input work product used for the system test plan is the:

Software requirements specification documents – refer to DHF.

Software Detailed Design documents – refer to DHF.

Software System Verification and Validation Plan – refer to DHF

2.7 Test Environment

The tests require a system with the software fully installed and the options (hardware and software) available.

2.7.1 Hardware

Flash N' Go device consists of a Base unit and an Applicator.

2.7.2 Software

Software version: 1.0.0

2.7.3 Security

NA.

2.7.4 Tools

2.7.4.1 Error reporting Tools

An Excel file named: bugsreport.xls.

2.7.4.2 Publications

NA.

2.8 Roles & Responsibilities

2.8.1 Test Plan generation

Functional test plan is generated by the system engineer.

Test Procedures for safety related components shall be defined by a system engineer and a team regulatory& quality member.

2.8.2 Test Execution

Test Procedures shall be executed by a team member not directly involved in the design and implementation of the component.

2.8.2.1 Test Phases and Cycles

There will be two main stages of testing for the Software during system testing:

System Testing

Acceptance Testing

2.8.2.1.1 System Testing Cycles

The main thrust of the approach is to intensively test the SW in two releases, thus raising approximately 80% of errors in this period. With the majority of these errors fixed, standard and/or frequently used actions will be tested to prove individual elements and total system processing. Regression testing of outstanding errors will be performed on an ongoing basis.

When all errors (which potentially impact overall processing) are fixed, an additional set of test cases are processed to ensure the system works in an integrated manner. It is intended that Release Flash N' GoSW 1.2.0 be the final proving of the system as a single application. There should be no severity 1 or 2 errors outstanding prior to the start of Release.

| Revision number | Creation date | System Test Cycle |
|------------------|---------------|-------------------|
| Flash N' Go1.0.0 | Jul 27; 2008 | First Cycle |
| Flash N' Go1.0.0 | Jul 27; 2008 | Second Cycle |

2.8.2.1.2 Acceptance-Testing Cycles

The main thrust of the approach is to test the usability of the device after approval of the system test.

2.8.2.2 Software Delivery

During System Test the release of new versions of the software will be co-ordinates between the Development SW engineer and the Development System engineer. However, unless it concerns a fix to a very serious error, new versions should only be released when agreed targets have been reached (i.e. the next version contains fixes to X or more numbers of bugs).

Notes:

It is intended that 80% of the functionality will have been tested in full prior to the Phase 3 Release.

All the functionality must be present in the Phase 2 Release. No previously undelivered functionality will be accepted for testing after Phase 2.

2.9 Formal Reviewing

There will be several formal review points before and during system test, including the review of this document. This is a vital element in achieving a quality product.

- Design Documentation - Requirements Specification & Functional Specification
- System Test Plans & Test conditions.
- System Test results
- Post system test review.

3. Test Items

All items that makeup the software Flash N' Go device will be tested during the system test. The versions to be tested will be downloaded to the device

4. Features to be tested

- Resets and Initialization
- Set energy level
- Pulse handler
- Monitor the Capacitor s voltage.
- Hardware checks.

5. Features not to be tested

NA

6. Tester Skills

The following staff is needed to carry out this testing project:

System engineer: he provides the overall management of the testing and technical testing expertise.

SW engineer: he transmits the system to be tested and responds to the system test Incident reports he does any program debugging that required. It also supplies the version control.

All the staff must be trained to perform the test.

7. Safety Issues

Safety issues have been analyzed. This test plan takes into account this analyze and all the tests associated to a safety issue are noted "safety related" in the title refer to risk analysis management file.

8. Test Termination Criteria

This test plan terminates normally if the execution of each test procedures results in a passed condition. The test plan may be terminated abnormally under the following conditions:

- Upon reaching a "Fail" condition on for any test case within a test procedure, the tester may terminate the test procedure.
- Upon reaching a failed condition in one test procedure, the tester may continue testing if the preconditions of subsequent test procedures can be satisfied. Otherwise, the tester may terminate execution of the test plan.

If the preconditions of a test procedure cannot be satisfied, the tester may proceed to execute subsequent test procedure as long as the preconditions of those procedures can be satisfied. Otherwise, the tester may terminate execution of the test plan."

9. Test Procedures

9.1 Verify Interfaces

9.1.1 Verify User Interfaces

| ID | Action | Expected result | Test result | Comments and Recommendations |
|------------------------|---------|-----------------|-------------|------------------------------|
| <i>General actions</i> | | | | |
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |

9.1.2 Verify Self test

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---|-----------------|-------------|------------------------------|
| 1 | Disconnect Power Supply Switch (disconnect R21) | ---- | ---- | |
| 2 | (b) (4) | | | |

9.1.3 Verify Polling

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---|-----------------|-------------|------------------------------|
| 1 | Turn off the system and install new valid cartridge | ---- | ---- | |
| 2 | Turn on the system and set to energy level (1), wait for Ready Led to turn ON | ---- | ---- | |
| 3 | (b) (4) | | | |

(b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.1.5 Verify firmware download process

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |

(b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.2 Verify Functionality

9.2.1 Verify system initialization

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.2.2 Verify system initialized with corrupted firmware

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.2.3 Verify system set energy level (general actions)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|------------------------|---------|-----------------|-------------|------------------------------|
| <i>General actions</i> | | | | |
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | (b) (4) | | | |
| 5 | | | | |
| 6 | | | | |

9.2.4 Verify system set energy level (Treatment)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|-----------------------|-----------------|-------------|------------------------------|
| | <i>Treatment Mode</i> | | | |
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 7 | | | | |
| 8 | | | | |
| 9 | | | | |
| 10 | | | | |

| <i>ID</i> | <i>Action</i> | <i>Expected result</i> | <i>Test result</i> | <i>Comments and Recommendations</i> |
|-----------|---------------|------------------------|--------------------|-------------------------------------|
| 11 | (b) (4) | | | |
| 12 | | | | |
| 13 | | | | |
| 14 | | | | |
| 15 | | | | |

9.2.5 Verify system calibration energy level

| <i>ID</i> | <i>Action</i> | <i>Expected result</i> | <i>Test result</i> | <i>Comments and Recommendations</i> |
|-----------|---------------|------------------------|--------------------|-------------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 7 | | | | |
| 8 | | | | |
| 9 | | | | |



Home Skinovations Ltd.
Apolo Building.
Soltam Industrial Zone.
Yokneam, 20692, Israel

| | | |
|----|---------|--|
| 10 | (b) (4) | |
| 11 | | |
| 12 | | |

9.2.6 Verify new system level limitation

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 7 | | | | |
| 8 | | | | |
| 9 | | | | |
| 10 | | | | |
| 11 | | | | |
| 12 | | | | |
| 13 | | | | |
| 14 | | | | |
| 15 | | | | |

| | | |
|----|---------|--|
| 10 | (b) (4) | |
| 11 | | |
| 12 | | |

9.2.7 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 3 | | | | |
| 4 | | | | (b) (4) |

9.3 Verify System Recovery

9.3.1 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.3.2 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.3.3 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |

9.3.4 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |

| | | | | |
|---|---------|-------|---------|--|
| | | apart | | |
| 4 | (b) (4) | | | |
| 5 | | | | |
| 6 | | | (b) (4) | |

9.3.5 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |

9.3.6 Verify system ignores wrong energy level calibration value

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 7 | | | | |
| 8 | | | | |

i

9.3.7 Verify system response to miss-matched color

| <i>ID</i> | <i>Action</i> | <i>Expected result</i> | <i>Test result</i> | <i>Comments and Recommendations</i> |
|-----------|---------------|------------------------|--------------------|-------------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 5 | | | | |
| 4 | | | | (b) (4) |

9.4 Verify Safety

9.4.1 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |

9.4.2 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |

9.4.3 (b) (4)

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |
| 2 | | | | |
| 3 | | | | |

9.4.4

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | ----- | ----- | |
| 2 | (b) (4) | (b) (4) | (b) (4) | |

10. Test Plan Results Summary

| ID | Action | Expected result | Test result | Comments and Recommendations |
|----|---------|-----------------|-------------|------------------------------|
| 1 | (b) (4) | | | |

Tester Name: _____

Date: _____

Signature: _____

Comments: _____

Approved by: _____

Date: _____

Signature: _____

Comments: _____

Appendix 7 – Clinical Results

Background

(b) (4)

(b) (4) In this clinical study we demonstrated that the Flash N' Go is safe for all users. (b) (4)

(b) (4)

Name and Intended Use of Device

Device Name: Flash N' Go

Common: Pulsed Light Device

Intended Use: The Flash N' Go is intended for removal of unwanted hair.

Description of Experiment

Tests were carried out on (b) users with a variety of skin types, (b) (4) and on a variety of body locations. (b) (4)

(b) (4) No adverse effects were noted on all of the (b) users. (b) (4)

(b) (4)

Results

(b) (4)

341

(b) (4)



Conclusion

The Flash N' Go device (b) (4) is safe for use on all skin colors. (b) (4)


(b) (4)



**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS
PREMARKET NOTIFICATION TRUTHFUL AND
ACCURACY STATEMENT**

(As required by 21 CFR 807.87(j))

I certify that in my capacity as VP of Regulatory Affairs of Home Skinovations Ltd., I believe, to the best of my knowledge that all data and information submitted in this premarket notification are truthful and accurate and that no material fact has been omitted.



Dr. Amir Waldman

August 10, 2008

510(k) Number (if known)_____.

Device Name Flash N' Go

Indications For Use:

The Flash N' Go is an over the counter device intended for the removal of unwanted hair.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use_____
(Per 21 CFR 801.109)

OR

Over The Counter Use_____

(Optional Format 1-2-96)

Performance & clinical data:

The device complies with the following U.S. Food and Drug Administration performance standards: 21CFR § 1040.10 & 1040.11. Clinical safety data was collected in a prospective safety study.


Substantial Equivalence:

The Flash N' Go system is substantial equivalent to its predicate devices. The data in this 510k submission demonstrate that the Flash N' Go system has identical output as the Silk'n, and shares the similar intended use as other predicate devices. Therefore is substantial equivalent to its predicate devices. Details are provided in chapter 3: Substantial equivalent of this submission.

Based upon an analysis of the overall performance characteristic for the device, Home Skinovations Ltd. believes that no significant differences exist. Therefore the Flash N' go should raise no new issues of safety or effectiveness.

August 10, 2008

Date



Dr. Amir Waldman,
VP Regulatory Affairs
Home Skinovations Ltd.

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Home SkinovationsLtd.

Flash N' Go

This summary of safety and effectiveness information is being submitted in accordance with the requirements of the SMDA 1990 and 21 CFR 807.92.

Submitter's information

Name: Home Skinovations Ltd.
Address: Apolo building, POB 533, Yokneam 20692, Israel
Contact: Dr. Amir Waldman VP Regulatory Affairs

Device information

Trade/Proprietary name: Flash N' Go
Common/Usual name: Light based hair removal device
Classification name: Laser surgical instrument for use in general and plastic surgery and in dermatology (21CFR §878.4810)
Product code: GEX

Predicate devices

- ABC hair removal system, (K060839)
Palomar Medical Technologies Inc.
- Spectra hair removal laser system, (K053527).
SpectraGenics, Inc.
- Silk'n, (K072906).
Home Skinovations Ltd.

Intended use:

The Flash N' Go is intended for removal of unwanted hair.

Device Description:

The Flash N' Go hair removal system is composed of a base unit and hand held applicator. Details are provided in the chapter 1: Device description, of this submission.



COVER SHEET MEMORANDUM

From: Reviewer Name

Subject: 510(k) Number

To: The Record

Kareem S. Burney
K103184

Please list CTS decision code

SE

- ☐ Refused to accept (Note: this is considered the first review cycle, See Screening Checklist http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_5631/Screening%20Checklist%207%202%2007.doc)
- ☐ Hold (Additional Information or Telephone Hold).
- ☒ Final Decision (SE, SE with Limitations, NSE, Withdrawn, etc.).

| Please complete the following for a final clearance decision (i.e., SE, SE with Limitations, etc.): | | YES | NO |
|---|--------------------------------------|-----|----|
| Indications for Use Page | Attach IFU | ✓ | |
| 510(k) Summary /510(k) Statement | Attach Summary | ✓ | |
| Truthful and Accurate Statement. | Must be present for a Final Decision | ✓ | |
| Is the device Class III? | | | ✓ |
| If yes, does firm include Class III Summary? | Must be present for a Final Decision | | ✓ |
| Does firm reference standards? (If yes, please attach form from http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf). | | | ✓ |
| Is this a combination product? (Please specify category _____, see http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/COMBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC) | | | ✓ |
| Is this a reprocessed single use device? (Guidance for Industry and FDA Staff – MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices, http://www.fda.gov/cdrh/ode/guidance/1216.html) | | | ✓ |
| Is this device intended for pediatric use only? | | | ✓ |
| Is this a prescription device? (If both prescription & OTC, check both boxes.) | | | ✓ |
| Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? | | | ✓ |
| Is clinical data necessary to support the review of this 510(k)? Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If not, then applicant must be contacted to obtain completed form.) | | | ✓ |
| Does this device include an Animal Tissue Source? | | | |
| All Pediatric Patients age ≤ 21 | | | |
| Neonate/Newborn (Birth to 28 days) | | | |
| Infant (29 days - < 2 years old) | | | |
| Child (2 years - < 12 years old) | | | |
| Adolescent (12 years - < 18 years old) | | | |
| Transitional Adolescent A (18 - < 21 years old) Special considerations are being given to this group, different from adults age ≥ 21 (different device design or testing, different protocol procedures, etc.) | | | |

Transitional Adolescent B (18 -<= 21; No special considerations compared to adults => 21 years old)

Nanotechnology

Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance, <http://www.fda.gov/cdrh/comp/guidance/169.html>) Contact OC.

Regulation Number

878.4810

Class*

II

Product Code

OHT

(*If unclassified, see 510(k) Staff)

Additional Product Codes:

Review:



(Branch Chief)


GSD3

(Branch Code)

11/10/00

(Date)

Final Review:


(Division Director)

(Date)

ACK DIA
11/10/00

**DEPARTMENT OF HEALTH AND HUMAN SERVICES****MEMORANDUM**

Food and Drug Administration
Office of Device Evaluation
9200 Corporate Boulevard
Rockville, MD 20850

**Premarket Notification [510(k)] Review
Traditional/Abbreviated**

K103184

Date: November 9, 2010
To: The Record
From: Kareem S. Burney, M.S.

Office: ODE
Division: DSORD

510(k) Holder: Home Skinovations LTD
Device Name: Silk'N Flash N Go
Contact: Amir Waldman
Phone: 972 547 404616 (Israel)
Fax:
Email:

Grass
10/11/10

I. Purpose and Submission Summary

The 510(k) holder would like to introduce the Silk'N Flash N Go into interstate commerce.

III. Administrative Requirements

| | Yes | No | N/A |
|---|-----|----|-----|
| Indications for Use page (Indicate if: Prescription or OTC) | X | | |
| Truthful and Accuracy Statement | | X | |
| 510(k) Summary or 510(k) Statement | X | | |
| Standards Form | | | X |

III. Device Description

| | Yes | No | N/A |
|---|-----|----|-----|
| Is the device life-supporting or life sustaining? | | X | |
| Is the device an implant (implanted longer than 30 days)? | | X | |
| Does the device design use software? | X | | |
| Is the device sterile? | | X | |
| Is the device reusable (not reprocessed single use)? | X | | |
| Are "cleaning" instructions included for the end user? | | | |

- o The Flash N go hair removal system is pulsed light system composed of a base unit and hand help applicator.

II. Indications for Use

Flash N Go is also intended for permanent reduction in hair regrowth defined as a long-term, stable reduction in hair counts following a treatment regime.

This device is for over the counter use.

III. Predicate Device Comparison

Flash N' Go (K082298)

TRIA Hair Removal System (K090820)

EpiLight and PhotoDem HR K991935

The sponsor has provided a substantial equivalence table. The Flash N Go for this device has a similar wavelength, energy as the Flash N Go (K082298) The wavelength is 475-1200nm and a energy of 5 J/cm². The indication for K082298 is just for the removal of unwanted hair. **Because these devices are similar the sponsor should also provide the spot size and pulse repetition rate. This information is not in the submission.**

The TRIA has a wavelength of 810nm and energy of 20 J/cm². The indication is intended for the adjunctive use with shaving fo rhair ernoval sustained with periodic treatments and for permanent reduction in hair regrowth defined as a long term stable reduction in hair counts following a treatment regime.

The EpiLight/PhotoDerm has a wavelength range of 645-1200nm and energy of 45 J/cm. The indication is for the removal of unwanted hair. To stable long-term or permanent hair reduction in skin types I – V through selective targeting of melanin in hair follicles. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs regrowing after a treatment regime.

IV. Labeling

The sponsor states that the only change to the labeling compared to previously cleared Flash N Go is the new indication for use.

The sponsor has made some changes to the instruction for use such as providing a device description, information about the package contents an disposable lamp cartridge. This information was in the instruction for use but the sponsor moved it up to the front of the instruction for use.

After reviewing the instruction for use for both the current submission and K082298's instruction for use. The instruction for use is the same as K082298. The sponsor will not have to provide usability study, label comprehension on this instruction for use. The sponsor will have to provide the box labeling.

However there are some inconsitencies with the labeling. On page 19 the sponsor states that the cleaning section is 5.1 when it is actually 11.1.

The sponsor put in the labeling that the clincial data is based on 1,3, and 18 month follow up data.

V. Sterilization/Shelf Life/Reuse

The sponsor states that the sterilization information is the same as the K082298. Here is the information from the K082298 submission:

Cleaning instructions have been included in the instruction for use. The cleaning instructions are as follows: " Unplug Flash N Go before cleaning. Use a dry clean cloth and specially formulated cleaner for electronic equipment to gently wipe Flash N Go surface and especially the Applicator tip. Never immerse Flash N G or any of its parts in water.

VI. Biocompatibility

The sponsor states that the biocompatibility information is the same as the K082298. Here is the information from the K082298 submission:

The aperture comes in contact with the patient. The materials that are in direct contact with the skin are identical to the Silk'N K072096. Sponsor has stated that parts in contact with the patient but needs to state whether this device complies the requirement of ISO 10993 or EN 30993-1.

VII. Software

The sponsor has stated that the software is the same as the predicate device K082298. The software is a moderate level of concern. Please refer to K082298 for information about the software.

| | | |
|---------------------------------------|-----|----|
| Version: | | |
| Level of Concern: Moderate | | |
| | Yes | No |
| Software description: | | |
| Device Hazard Analysis: | | |
| Software Requirements Specifications: | | |
| Architecture Design Chart: | | |
| Design Specifications: | | |
| Traceability Analysis/Matrix: | | |
| Development: | | |
| Verification & Validation Testing: | | |
| Revision level history: | | |
| Unresolved anomalies: | | |

VIII. Electromagnetic Compatibility and Electrical, Mechanical and Thermal Safety

The EC Declaration of Conformity showing that is in full compliance with the following standards have been provided: EN 61000-3-2:2000; EN 61000-3-3:1995; EN 55014-1:2000; EN 55014-2: 1997; EN 60335-1:2002; EN 60335-2-23:2003; EN 60335-2-27:2003; and EN 50364:2001

IX. Performance Testing – Bench

X. Performance Testing – Animal

XI. Performance Testing – Clinical

The sponsor stated that they carried out clinical trials to verify the safety and effectiveness of the Silk'N Flash and Go.

(b) (4)



The treatment test spot is as follows:

- 1) In the area to be treated hairs should be trimmed to 3/32 inch length, (1-2mm) or shaved 3 days prior to the treatment.
- 2) Clean the skin from trimmed hairs.
- 3) Hold the applicator to ensure that the light output window is clean.
- 4) Switch on the system and set the energy level.
- 5) Place the Applicator on the treatment area and apply slight pressure. Press the trigger switch to emit light pulse.
- 6) Move the applicator to the next spot to ensure full coverage of the area to be treated.
- 7) Repeat the test spot procedure for 3 pulses.
- 8) If side effects are noted after 15 minutes, reduce the energy level by one level to lower setting.
- 9) Transient Erythema and follicular redness are normal reactions as well as smell of burned hair.
- 10) If no side effects are noted after 15 minutes, complete the treatment to cover entire area to be treated.

Post Treatment Care:

Apply moisturizer after the treatment.

Use high factor sunscreen and protect the treated area from sunlight for at least 2 days after the treatment. Tanning after a treatment may cause hyper-pigmentation.

Repetitive Treatments:

(b) (4)

Complications:

If an unanticipated adverse event occurs at any time during or after the use of the system, the Investigator must report it to the sponsor. A serious adverse event or complication directly attributable to the use of the System will be reported to the Sponsor as soon as possible within twenty-four (24) hours or sooner.

Data Analysis:

Hair loss at each follow up session will be calculated as the ratio of the number of hairs after the treatment compared with the hair count prior to treatment. Hair count will be performed by the investigator.

This is the protocol from K072906. The original submission for the Flash N' Go:

1. In the area to be treated, hairs should be trimmed to 3/32 inch length, (1-2 mm), or shaved 3 days prior to the treatment.
2. Clean the skin from the trimmed hairs.
3. Hold the applicator and ensure that the light output window is clean.
4. Switch on the system and set the energy level.
5. Place the Applicator on the treatment area and apply slight pressure. Press the trigger switch to emit light pulse.
6. Move the applicator to the next spot to ensure full coverage of the area to be treated.
7. Repeat the test spot procedure for 3 pulses.
8. If side effects are noted after 15 minutes, reduce the energy level by one level to lower setting.
9. Transient Erythema and follicular redness are normal reactions as well as smell of burned hairs.
10. If no side effects are noted after 15 minutes, complete the treatment to cover entire area to be treated.

Post Treatment Care:

Apply moisturizer after the treatment.

Use high factor sunscreen and protect the treated area from sunlight for at least 2 days after the treatment. Tanning after a treatment may cause hyper-pigmentation.

(b) (4)

Complications:

If an unanticipated adverse event occurs at any time during or after the use of the system, the Investigator must report it to the sponsor. A serious adverse event or complication directly attributable to the use of the System will be reported to the Sponsor as soon as possible within twenty-four (24) hours or sooner.

Data Analysis:

Hair loss at each follow up session will be calculated as the ratio of the number of hairs after the treatment compared with the hair count prior to treatment. Hair count will be performed by the investigator.

(b) (4)



The sponsor stated that the only adverse event was edema, peri follicular, edema and light erythema which were resolved in 60 minutes. No other side effects were noted by the sponsor.

(b) (4)



XII. Substantial Equivalence Discussion

| | Yes | No |
|--|-----|-------------------------------------|
| 1. Same Indication Statement? | | If YES = Go To 3 |
| 2. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness? | | If YES = Stop NSE |
| 3. Same Technological Characteristics? | | If YES = Go To 5 |
| 4. Could The New Characteristics Affect Safety Or Effectiveness? | | If YES = Go To 6 |
| 5. Descriptive Characteristics Precise Enough? | | If NO = Go To 8 If YES = Stop SE |
| 6. New Types Of Safety Or Effectiveness Questions? | | If YES = Stop NSE |
| 7. Accepted Scientific Methods Exist? | | If NO = Stop NSE |
| 8. Performance Data Available? | | If NO = Request Data |
| 9. Data Demonstrate Equivalence? | | Final Decision: |

Note: See

http://erom.fda.gov/eRoomReg/Files/CDRH3/CDRHPremarketNotification510kProgram/0_4148/FLOWCHART%20DECISION%20TREE%20.DOC for Flowchart to assist in decision-making process. Please complete the following table and answer the corresponding questions. "Yes" responses to questions 2, 4, 6, and 9, and every "no" response requires an explanation.

1. Explain how the new indication differs from the predicate device's indication:
2. Explain why there is or is not a new effect or safety or effectiveness issue:
3. Describe the new technological characteristics:
4. Explain how new characteristics could or could not affect safety or effectiveness:
5. Explain how descriptive characteristics are not precise enough:
6. Explain new types of safety or effectiveness question(s) raised or why the question(s) are not new:
7. Explain why existing scientific methods can not be used:
8. Explain what performance data is needed:
9. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

XIII. Deficiencies

XIV. Contact History

XV. Recommendation

I find no safety and effective issues with the changes presented in 510(k) K103184. Therefore, I find the proposed device to be substantially equivalent to the predicate.

XII. Substantial Equivalence Discussion

| | Yes | No | |
|--|-----|----|-------------------------------------|
| 1. Same Indication Statement? | ✓ | | If YES = Go To 3 |
| 2. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness? | | | If YES = Stop NSE |
| 3. Same Technological Characteristics? | ✓ | | If YES = Go To 5 |
| 4. Could The New Characteristics Affect Safety Or Effectiveness? | | | If YES = Go To 6 |
| 5. Descriptive Characteristics Precise Enough? | ✓ | ✓ | If NO = Go To 8 If YES = Stop SE |
| 6. New Types Of Safety Or Effectiveness Questions? | | | If YES = Stop NSE |
| 7. Accepted Scientific Methods Exist? | | | If NO = Stop NSE |
| 8. Performance Data Available? | ✓ | | If NO = Request Data |
| 9. Data Demonstrate Equivalence? | ✓ | | Final Decision: SE |

Note: See

http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_4148/FLOWCART%20DECISION%20TREE%20.DOC for Flowchart to assist in decision-making process. Please complete the following table and answer the corresponding questions. "Yes" responses to questions 2, 4, 6, and 9, and every "no" response requires an explanation.

1. Explain how the new indication differs from the predicate device's indication:
2. Explain why there is or is not a new effect or safety or effectiveness issue:
3. Describe the new technological characteristics:
4. Explain how new characteristics could or could not affect safety or effectiveness:
5. Explain how descriptive characteristics are not precise enough:
6. Explain new types of safety or effectiveness question(s) raised or why the question(s) are not new:
7. Explain why existing scientific methods can not be used:
8. Explain what performance data is needed:
9. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

XIII. Deficiencies**XIV. Contact History**

XV. Recommendation

I find no safety and effective issues with the changes presented in 510(k) K103184. Therefore, I find the proposed device to be substantially equivalent to the predicate.

Regulation Number: 878.4810

Regulation Name: Laser Instrument, Surgical, Powered

Regulatory Class: II

Product Code: OHT

Kareem S. Bury
Reviewer

Neil A. Ozden
Branch Chief

11/9/10
Date

11/10/10
Date

Burney, Kareem

From: Amir Waldman [amirw@silkn.com]
Sent: Tuesday, November 09, 2010 8:40 AM
To: Burney, Kareem
Subject: RE: FDA review of Flash N Go (K103184)

The number of treatment in the current study was 6.
 On the K082298/K072906 the number of treatment was 3.
 Amir

From: Burney, Kareem [mailto:Kareem.Burney@fda.hhs.gov]
Sent: Tuesday, November 09, 2010 3:04 PM
To: 'Amir Waldman'
Subject: RE: FDA review of Flash N Go (K103184)

Amir,
 I was looking at your clinical data and I have a question. What was the total number of treatments done in this study? Can you compare this to the total number of treatment done in the study for K082298/K072906? Thanks.

*Kareem S. Burney, M.S.
 Biomedical Engineer
 Division of Surgical, Orthopedic and Restorative Devices
 U.S. Food and Drug Administration
 WO66-1439
 10903 New Hampshire Ave
 Silver Spring, MD 20993
 Phone # (301)796-6388
 E-mail: kareem.burney@fda.hhs.gov*

From: Amir Waldman [mailto:amirw@silkn.com]
Sent: Monday, November 08, 2010 3:10 AM
To: Burney, Kareem
Subject: RE: FDA review of Flash N Go (K103184)

Chapter 2 on page 4 Flash N GO intended use changed to :

2. Flash N' Go™ Intended Use

Flash N' Go™ is intended for removal of unwanted hair. Flash N' Go™ is also intended to effect long term, or permanent hair reduction. Permanent hair reduction is defined as a long-term stable reduction in the number of hairs re-growing after a treatment regime.

Flash N' Go™ may be used to remove unwanted body hair. Ideal body areas for Flash N' Go™ use includes the underarms, bikini line, arms and legs.

(b) (4)

(b) (4)

Amir

From: Burney, Kareem [mailto:Kareem.Burney@fda.hhs.gov]
Sent: Friday, November 05, 2010 4:49 PM
To: 'Amir Waldman'
Subject: RE: FDA review of Flash N Go (K103184)

(b) (4)

*Kareem S. Burney, M.S.
Biomedical Engineer
Division of Surgical, Orthopedic and Restorative Devices
U.S. Food and Drug Administration
WO66-1439
10903 New Hampshire Ave
Silver Spring, MD 20993
Phone # (301)796-6388
E-mail:kareem.burney@fda.hhs.gov*

From: Amir Waldman [mailto:amirw@silkn.com]
Sent: Friday, November 05, 2010 1:39 AM
To: Burney, Kareem
Subject: RE: FDA review of Flash N Go (K103184)

Mr. Burney,

Attached is chapter 2 clinical information in a pdf format, is this OK or you prefer that I'll send printed color copy?

Sorry for not putting the detailed terms but here they are:

(b) (4)

11/9/2010

17

Amir

From: Burney, Kareem [mailto:Kareem.Burney@fda.hhs.gov]
Sent: Thursday, November 04, 2010 4:52 PM
To: 'Amir Waldman'
Subject: RE: FDA review of Flash N Go (K103184)

(b) (4)



*Kareem S. Burney, M.S.
Biomedical Engineer
Division of Surgical, Orthopedic and Restorative Devices
U.S. Food and Drug Administration
WO66-1439
10903 New Hampshire Ave
Silver Spring, MD 20993
Phone # (301)796-6388
E-mail:kareem.burney@fda.hhs.gov*

From: Amir Waldman [mailto:amirw@silkn.com]
Sent: Wednesday, November 03, 2010 2:08 PM
To: Burney, Kareem
Subject: RE: FDA review of Flash N Go (K103184)

Yes,

(b) (4)



Amir Waldman
Home Skinovations Ltd.

From: Burney, Kareem [mailto:Kareem.Burney@fda.hhs.gov]
Sent: Wednesday, November 03, 2010 7:59 PM
To: 'Amir Waldman'
Subject: RE: FDA review of Flash N Go (K103184)

(b) (4)



*Kareem S. Burney, M.S.
Biomedical Engineer
Division of Surgical, Orthopedic and Restorative Devices
U.S. Food and Drug Administration*

11/9/2010

18

WO66-1439
10903 New Hampshire Ave
Silver Spring, MD 20993
Phone # (301)796-6388
E-mail:kareem.burney@fda.hhs.gov

From: Amir Waldman [mailto:amirw@silkn.com]
Sent: Wednesday, November 03, 2010 1:25 PM
To: Burney, Kareem
Subject: RE: FDA review of Flash N Go (K103184)

(b) (4)

Amir

From: Burney, Kareem [mailto:Kareem.Burney@fda.hhs.gov]
Sent: Wednesday, November 03, 2010 6:52 PM
To: 'amirw@silkn.com'
Subject: FDA review of Flash N Go (K103184)

Amir,
I have done look at some of this submission and I have the following questions. Are the patient contact materials and sterilization information the same as K082298? What we mean by sterilization is whether there are sterilization information or the device can be cleaned like K082298.

Kareem S. Burney, M.S.
Biomedical Engineer
Division of Surgical, Orthopedic and Restorative Devices
U.S. Food and Drug Administration
WO66-1439
10903 New Hampshire Ave
Silver Spring, MD 20993
Phone # (301)796-6388
E-mail:kareem.burney@fda.hhs.gov

Standard Operating Procedures for 510(k) Summaries


- 1) If a 510(k) submitter decides to submit a 510(k) summary, as per 21 CFR 807.87(h), they need to follow the content requirements as per 21 CFR 807.92.
- 2) If during the review of the 510(k), the submitter decides to switch to a 510(k) statement, as per 21 CFR 807.87(h), they may do so while the 510(k) is under review. They must follow 21 CFR 807.93 for a 510(k) statement.
- 3) The 510(k) summary is written by the 510(k) submitter, but FDA will agree with the content prior to clearing any 510(k).
- 4) The 510(k) summary needs to agree with the final classification decision of FDA, e.g., the predicate classification needs to match the FDA classification decision. The submitter may need to revise the summary while the 510(k) is under review. In other words, the 510(k) summary will need to reflect the predicate(s) and decision made by FDA. The submitter will then need to revise the summary while the 510(k) is under review.
- 5) The IFU provided in the 510(k) summary needs to match the IFU statement that is determined to be substantially equivalent
- 6) The 510(k) summary should reflect all the testing done by the 510(k) submitter to demonstrate substantial equivalence. This may include testing that FDA did/would not require to demonstrate substantial equivalence, but would include all testing that FDA does/would require to demonstrate substantial equivalence.
- 7) If the 510(k) summary is deficient, the deficiency(ies) may be put in an AI letter or handled through interactive review.
- 8) The 510(k) may not be found to be substantially equivalent until the 510(k) summary meets the regulatory requirements of 21 CFR 807.92.
- 9) If, after interactions, a 510(k) submitter does not revise the 510(k) summary as requested and we disagree with their rationale, the 510(k) may be found not substantially equivalent for lack of required data/information.
- 10) Neither the 510(k) summary, nor 510(k) statement, are needed if the decision is other than SE.
- 11) The 510(k) summary will go on FDA's website approximately the 5th of the month following any SE decision.

510(k) SUMMARY REQUIREMENTS CHECKLIST 21 CFR 807.92

| | | Y | N | N/A |
|---|--|---|---|-----|
| All 510(k) summaries shall contain the following information: | | | | |
| 1 | The submitter's name, address, telephone number, a contact person, and the | ✓ | | |

| | | | | |
|---|--|---|--|--|
| | date the summary was prepared | | | |
| 2 | The name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name | ✓ | | |
| 3 | An identification of the legally marketed device(s) to which the submitter claims equivalence. | ✓ | | |
| 4 | A description of the device that is the subject of the 510(k), including an explanation of how the device functions, the scientific concepts that form the basis for the device, and the significant physical and performance characteristics of the device (e.g., device design, material used, and physical properties) | ✓ | | |
| 5 | A statement of the indications for use of the device that is the subject of the 510(k), including a general description of the diseases or conditions that the device will diagnose, treat, prevent, cure, or mitigate, including a description, where appropriate, of the patient population for which the device is indicated. Or, if the indication statements are different from those of the legally marketed device(s) identified in paragraph (3) of this section, an explanation as to why the differences are not critical to the intended therapeutic, diagnostic, prosthetic, surgical or other use of the device, and why the differences do not affect the safety and effectiveness of the device when used as indicated. | ✓ | | |
| 6 | If the device has the same technological characteristics (i.e., design, material, chemical composition, energy source, etc.) as the predicate device(s) identified in paragraph(3) of this section, a summary of the technological characteristics of the new device in comparison to those of the predicate device(s). Or, if the device has different technological characteristics from the predicate device(s), a summary of how the technological characteristics of the device compare to a legally marketed device(s) identified in paragraph (3) of this section. | ✓ | | |
| 510(k) summaries for those 510(k)s in which a determination of substantial equivalence is also based on an assessment of performance data shall contain the following information | | | | |
| 7 | A brief discussion of the nonclinical tests submitted, referenced, or relied on in the 510(k) for a determination of substantial equivalence | | | |
| 8 | A summary discussion of the clinical tests submitted, referenced, or relied on in the 510(k) for a determination of substantial equivalence. This discussion shall include, where applicable, a description of the subjects upon whom the device was tested, a discussion of the safety or effectiveness data obtained from the testing, with specific reference to adverse effects and complications, and any other information from the clinical testing relevant to a determination of substantial equivalence. (There can not be any patient identifier information in the summary.) | | | |
| 9 | The conclusions drawn from the nonclinical and clinical tests that demonstrate that the device is as safe, as effective, and performs at least as safely and effectively as the legally marketed device identified in paragraph(3) of this section. | | | |

510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



❖ 510
ma
❖❖ Thi
❖❖❖ Dat



医课汇
公众号
专业医疗器械资讯平台
WECHAT OF
HLONGMED



hlongmed.com
医疗器械咨询服务
MEDICAL DEVICE
CONSULTING
SERVICES



医课培训平台
医疗器械任职培训
WEB TRAINING
CENTER



医械宝
医疗器械知识平台
KNOWLEDG
ECENTEROF
MEDICAL DEVICE



MDCPP.COM
医械云专业平台
KNOWLEDG
ECENTEROF MEDICAL
DEVICE

if the relationship between
metimes required.