



Basic UDI-DI & UDI-PI attributes

Basic UDI-DI set of data in UDI database

Principle: Each UDI-DI inherits the attributes of its linked Basic UDI-DI and devices DI

Basic UDI-DI

Applicable legislation (IVDR) (*)

- **2. Basic UDI-DI value (*)**
- **2b Basic UDI-DI Issuing entity (*)**;
- **6. Manufacturer SRN (*)**
- **5. Name and address of manufacturer**
- **7. Name and address and SRN of AR**
- **9. Risk class (*)**
 - **A.2.14 Intended for self-testing (Y/N) (*)**
 - **A.2.14 Intended for near-patient-testing (Y/N) (*)**
 - **Companion diagnostics (Y/N) (*)**
 - **Instrument(Y/N) (*)**
 - **Reagent(Y/N) (*)**
 - **Professional testing (Y/N) (*)**
- **11. A. Name and/or, if applicable, device model that identifies the device(s) with this BASIC UDI-DI in the technical documentation and/or certificate or declaration of conformity**
(Name and/or model shall be provided)

UDI-DIs

- **0. UDI-DI value (*)**
- **0b. UDI-DI Issuing entity (*)**
- **Secondary DI (value and issuing entity)**
- **11.B. Reference, Article or Catalogue number (*)**
- **Device with Direct marking (Y/N) (*)**
- **Direct marking UDI-DI value (*)**
- **Direct marking UDI-DI issuing entity (*)**
- **1. Quantity of device(s) (*)**
- **3. Type of UDI-PI (*)**
- **4. Unit of use UDI-DI (*)**
- **13. Storage/handling conditions**
- **10-14. Name(s)/Trade name(s) (including languages)**
- **12. Additional product description**
- **19. URL for additional information**
- **15. Labelled as single use (YN) (*)**
- **16. Maximum number of reuse (*)**
- **17. Device labelled sterile (Y/N) (*)**
- **18. Need for sterilisation (Y/N) (*)**
- **20. Critical warnings or contra-indications**
- **8. Medical device nomenclature (CND) code (1)**
- **21. Status**
- **27 (A.2.10). In the case of devices designed and manufactured by another legal or natural person as referred in Article 10(14), the name, address and contact details of that natural/legal person**

UDI-DIs (container package DI)

- **0. UDI-DI value (*)**
- **0b. Issuing entity (*)**
- **1. Quantity per package (*)**
- **21. Status**

(1) Nomenclature decision:

https://ec.europa.eu/doc_sroom/documents/34264

(*) may not be changed

- **Mandatory**
- **Mandatory if applicable**
- **Optional**

Other Device Data attributes

Basic UDI-DI

- **A.2.2 Certificate IDs (with NB, type .. Link);**
- **A.2.11 SSP;**
- **A.2.9 Performance study IDs (..link);**
- **A.2.5 Presence of Human tissues/Cells (Y/N) (*);**
- **A.2.6 Presence of Animal tissues/Cells (Y/N) (*);**
- **A.2.7 Presence of Substances/cells of microbial origin (Y/N) (*);**
- **Kit (Y/N) (*);**

UDI-DIs

- **A.2.13 New Device (Y/N) (*);**
- **A.2.3 Member State of the Placing on the EU Market of the Device (*);**
- **A.2.4 Member State(s) where the Device is made available in the Country;**



Provided by
Art 26(2) pr
confirmed by

Version April



医课汇
公众号
专业医疗器械资讯平台
WECHAT OF
HLONGMED



hlongmed.com
医疗器械咨询服务
MEDICAL DEVICE
CONSULTING
SERVICES



医课汇培训平台
医疗器械任职培训
WEB TRAINING
CENTER



医械宝
医疗器械知识平台
KNOWLEDG
ECENTEROF
MEDICAL DEVICE



MDCPP.COM
医械云专业平台
KNOWLEDG
ECENTEROF MEDICAL
DEVICE

(*) may not be changed



Mandatory



Mandatory if applicable



Optional